



Jun 21,
2019

Barcode plasmid library cloning

[Benjamin Emert](#)¹

¹University of Pennsylvania

In devel.

[dx.doi.org/10.17504/protocols.io.4hgg3w](https://doi.org/10.17504/protocols.io.4hgg3w)

RajLab



Benjamin Emert
University of Pennsylvania



ABSTRACT

Protocol for cloning Time Machine barcode plasmid library.

MATERIALS

NAME	CATALOG #	VENDOR
Gibson Assembly Master Mix - 10 rxns	E2611S	New England Biolabs
Electroporation Cuvette 1mm	1652089	BioRad Sciences
Monarch® PCR & DNA Cleanup Kit (5 µg)	T1030	New England Biolabs
LB plates with 100 µg/ml ampicillin		
LB Broth with 100ug/mL ampicillin	244620	Gibco, ThermoFisher
Lucigen Endura electrocompetent cells	60242	Lucigen
ssDNA barcode template	View	IDT
Linearized vector backbone	View	
EndoFree Plasmid maxiprep kit	12362	Qiagen

Gibson assembly

1

Perform Gibson assembly using 100ng linearized vector and 5:1 molar excess of ssDNA template per 20 microliter reaction.

1.1 Pre-heat thermal cycler to 50°C with lid set to 65°C.

50 °C Pre-heat thermal cycler

1.2 Combine 100ng linearized vector with ~21.7 fmol ssDNA template (1.08 microliters of 100nM oligo) to a final volume of 10 microliters in nuclease free water. In parallel, perform a negative control reaction with plasmid and nuclease free water without the ssDNA insert. Keep on ice.

100 ng Linearized vector

1.08 µl ssDNA oligo (160nt)

 **10 µl Final volume in NF water**

1.3 Thaw Gibson assembly master mix on ice then add 10 microliters to plasmid plus insert. Pipette mix then spin down if needed.

 **10 µl Gibson Assembly master mix**

1.4 Incubate on thermal cycler at 50°C for 1 hour.

 **50 °C**  **01:00:00 incubation**

2 Purify Gibson reaction using PCR clean-up column according to manufacturer's protocol eluting in 20 microliters of

2.1 Add 40 microliters of DNA binding buffer to 20 microliters of Gibson assembly.

 **40 µl DNA binding buffer**

 **20 µl Gibson assembly**


2.2 Re-load eluant and re-spin.

2.3 Add 200 microliters DNA wash buffer. Make sure ethanol has been added to wash buffer before use.

 **200 µl Wash buffer**

2.4 Centrifuge at 12,000 x g for 1 minute.

 **12000 x g**

 **00:01:00**

2.5 Repeat wash and centrifugation (steps 2.4 and 2.5)

2.6 Transfer column to a new eppendorf or collection tube then centrifuge at 12,000 x g for 2 minutes to dry column

 **12000 x g**  **00:02:00**



Sometime residual liquid can collect on a lip inside the column. It may be helpful to rotate the column 180° after a 1 minute spin, then repeat another 1 minute spin to get rid of this liquid.

2.7 Add 10 microliter of nuclease free water pre-heated to 50°C to column and incubate for ≥ 2 minutes

 **10 µl nuclease free water**  **50 °C Pre-heat water**  **00:02:00 Incubate**

2.8 Centrifuge at 12,000 x g for 1 minute to elute

 **12000 x g**  **00:01:00**

2.9 Optionally nanodrop the recovered DNA. Expected concentration is 10-15ng/uL.

Electroporation

3 Electroporation of competent cells

3.1 Pre-warm Lucigen recovery media to 37°C

 **37 °C Pre-warm recovery media**

3.2 If using Biorad Gene Pulser, pre-chill cuvettes to 4°C.

 **4 °C Pre-chill electroporation cuvettes**

3.3 Thaw Lucigen Endura electrocompetent cells on ice. Each DUO vials contains 50uL of cells, enough for 2 electroporations. Use one reaction's worth of cells to electroporate the negative control Gibson assembly.

3.4 Pipet vial of electrocompetent cells into 2 25 microliter aliquots.

3.5 Add 2 microliters of cleaned Gibson assembly to 25 microliters of electrocompetent cells.

 **2 µl Gibson assembly** per  **25 µl electrocompetent cells**

3.6 Transfer electrocompetent cells to pre-chilled cuvette. Keep on ice.

3.7 Using the GenePulserXCell, electroporate with the following manual pulse setting: PL - 25msec, C - 10uF, PC- 600Ohms, V-1800



Pulse length - 25msec, Capacitance - 10 microFaraday, Resistance - 600 ohms, Voltage - 1800 Volts.



Expected time constant around 5-6 msec.

3.8 Immediately add 975 microliters of pre-warmed recovery media to cuvette then transfer to clean eppendorf or vial.

3.9 Incubate on shaker at 250 rpm at 37°C for 1 hour.
🔥 37 °C ⌚ 01:00:00 On shaker at 250rpm

3.10 Plate dilutions of electroporated cells for plating and calculating electroporation efficiency.
Prepare first dilution by adding 📄 10 µl electroporated cells to 📄 990 µl LB Broth
Prepare second dilution by adding 📄 100 µl First dilution to 📄 900 µl LB broth
Plate 📄 100 µl first dilution on pre-warmed LB/ampicillin agar plate for a 1:1000 dilution plate.
Plate 📄 100 µl second dilution on pre-warmed LB/Ampicillin agar plate for a 1:10,000 dilution plate.



After incubating the plates overnight, expect to see >1,000 colonies on the 1:1000 dilution plate and >200 colonies on the 1:10,000 plate.

3.11 Add the remainder of the electroporated cells to 150mL LB broth containing 100ug/mL ampicillin. Incubate 12-14 hours on shaker at 32°C.
📄 150 ml LB broth
🔥 32 °C
⌚ 12:00:00 - ⌚ 14:00:00 shaker

3.12 Pellet the overnight culture by centrifugation in 50mL falcon tubes then store pellets at -20°C or proceed immediately with maxiprep.
⚙️ 6000 x g ⌚ 00:15:00 🔥 4 °C Pellet e. coli cultures.
🔥 -20 °C Store e. coli pellets



If centrifuge capable of 6,000 x g is unavailable, can pellet at $\geq 3,000$ x g.

4 PCR screen colonies to confirm barcode insertion.

5 Maxiprep barcode plasmid library.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited