

# Chlorella Virus Plaque Assay Version 2

Dr. Steven Wilhelm, Samantha Coy

## Abstract

Adapted from: Van Etten, J. (n.d.). Titering of *Chlorella* Viruses. Retrieved from <http://ncv.unl.edu/vanettenlab/>  
Contact Dr. Steven Wilhelm (wilhelm@utk.edu) or Samantha Coy (srose16@vols.utk.edu) for additional information regarding this protocol.

**Citation:** Dr. Steven Wilhelm, Samantha Coy Chlorella Virus Plaque Assay. **protocols.io**

[dx.doi.org/10.17504/protocols.io.hgqb3vw](https://doi.org/10.17504/protocols.io.hgqb3vw)

**Published:** 29 Mar 2017

## Protocol

### MBBM Soft Agar

#### Step 1.

Melt MBBM soft agar and dispense in 3.0 mL aliquots and hold at 45°C - 50°C in a water bath

### Concentrate Chlorella

#### Step 2.

Concentrate Chlorella to  $4.0 \times 10^8$  cells/mL at 5,000 rpm for 5 min at 4°C and resuspend with fresh MBBM

\*Concentrate so that 0.3 mL can be used per plate

### Dilute Chlorella Virus

#### Step 3.

Dilute virus with 50 mM Tris-HCl, pH 7.8 in 1/10 serial dilutions

\*Fresh lysate contains approximately  $10^6$ - $10^7$  PFU/mL (check the titer prior to plaque assay)

\*Dilute sample to have 20-200 plaque forming units (PFU) per plate

### Titering

#### Step 4.

To each 3.0 mL soft agar aliquot, add 0.1 mL Chlorella virus and 0.3 mL Chlorella

\*Mix by rolling between palms and pour onto MBBM plate

\*Tilt the plate gently to allow the entire surface to be covered

\*Allow plate to solidify

## Incubation

### Step 5.

Incubate at 25°C in continuous light for 3-4 days or until plaques are visible to count

\*Use a sharpie to count plaques once there is a visible contrast