

Extraction method A (FMS and CR)

Faezah Mohd Salleh, Jazmin Ramos-Madrigal, Fernando Penalzoza, Shanlin Liu, Mikkel-Holger S Sinding, Riddhi P Patel, Renata Martins, Dorina Lenz, Jorns Fickel, Christian Roos, Mohd Shahir Shamsir, Mohammad Shahfiz Azman, Burton K Lim, Stephen J Rossiter, Andreas Wilting, M Thomas P Gilbert

Abstract

This protocol allows for adequate DNA extraction from fresh tissue samples.

Citation: Faezah Mohd Salleh, Jazmin Ramos-Madrigal, Fernando Penalzoza, Shanlin Liu, Mikkel-Holger S Sinding, Riddhi P Patel, Renata Martins, Dorina Lenz, Jorns Fickel, Christian Roos, Mohd Shahir Shamsir, Mohammad Shahfiz Azman, Burton K Lim, Stephen J Rossiter, Andreas Wilting, M Thomas P Gilbert Extraction method A (FMS and CR). **protocols.io** dx.doi.org/10.17504/protocols.io.im7cc9n

Published: 05 Jul 2017

Before start

Separate PCR-free facility

Materials

Buffer AW1 [19081](#) by [Qiagen](#)

Buffer AW2 [19072](#) by [Qiagen](#)

Buffer AE [19077](#) by [Qiagen](#)

Buffer ATL [19076](#) by [Qiagen](#)

✓ Ethanol by Contributed by users

Protocol

Extraction

Step 1.

Cut tissue into small pieces.

Extraction

Step 2.

Digestion of ca. 25 mg sample with Buffer ATL (Qiagen, Hilden, DE) and Proteinase K overnight at 56 °C.

 **AMOUNT**

25 mg Additional info: Sample

🔌 **NOTES**

GigaScience Database 26 Jun 2017

Incubation was performed on a thermomixer with 800 rpm

Extraction

Step 3.

After pelleting remaining material, remove supernatant and added to Buffer AL and ethanol.

Extraction

Step 4.

Transfer the solution to DNeasy Mini spin column (Qiagen, Hilden, DE) and centrifuge at 8000 rpm.

Extraction

Step 5.

DNA purification following manufacturer's protocol with Buffer AW1 and Buffer AW2 (Qiagen, Hilden, DE).

Extraction

Step 6.

Prior to elution, incubate DNA in the membrane with elution Buffer AE for 5 min at room temperature.

Extraction

Step 7.

Elution proceeds by centrifugation at 8000 rpm for 1 min.

Extraction

Step 8.

Measure DNA concentration on Nanodrop (Thermo Fischer Scientific, Darmstadt, DE).