



Oct 23, 2019

## PCR V.2

Claudia Troncone Clemente<sup>1</sup><sup>1</sup>Universidad Complutense de Madrid

1

Works for me

dx.doi.org/10.17504/protocols.io.8nyhvw

AEGIS - Madrid iGEM 2019



Laura Armero ⚡

## ABSTRACT

Our aim with this protocol is to amplify DNA. This protocol has been optimized has a general amplification

As the quantity of DNA is exponentially increased during the performance of the selection, further modification in the numbers of cycle will be needed to be implemented.

## GUIDELINES

We sure to have all the surfaces and materials clean before the start.

All the procedures must be done in an sterile environment to avoid contamination.

## MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Speedy Supreme Green Master Mix	MB39102	NZYtech
Agarose (LM-ultrapure grade)	MB123	NYZtech

## MATERIALS TEXT

- Aptamer library (order to IDT)

5' - G TTG CTC GTA TTT AGG GAA TG N<sub>40</sub> ACA CCA GTC TTC ATC CGC TTT<sub>6</sub> - 3'

- Forward primer (order to IDT):

G TAG GCG AAA<sub>6</sub> - Cy3 - 5'

- Reserve primer (order to IDT):

5' - BiodTG TTG CTC GTA TTT AGG GAA TG

- Thermocycler
- TAE buffer

- 1 Prepare the PCR reaction mixture following the specifications below:

Component	Positive control (V; ul)	Negative control (V; ul)
<i>Template</i>	5	0
<i>Fwd primer</i>	1.25	1.25
<i>Rev primer</i>	1.25	1.25
<i>dH<sub>2</sub>O</i>	15	20
<i>DMSO</i>	2.5	2.5
<i>Master Mix</i>	25	25

- 2 Perform the amplification in a general thermocycler in the following conditions. Adjust the annealing temperature according to the primers used, and the hotstart to the specifications of your polymerase:

Hot start	95 °	5 min	
Amplification cycles			
Denaturing	95 °	30s	
Annealing	52 °	30s	X 15 cycles
Extension	72 °	30s	
Final extension	72 °	3 min	
Hold	4°		

- 3 Prepare a 3% agarose gel. Load the samples and perform the electrophoresis at 90V for 50 min.
- 4 Remove the gel and observe the bands under UV light.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited