

# DNA Gels

David Dunigan and Irina Agarkova

## Abstract

**Citation:** David Dunigan and Irina Agarkova DNA Gels. [protocols.io](https://doi.org/10.17504/protocols.io.erubd6w)

[dx.doi.org/10.17504/protocols.io.erubd6w](https://doi.org/10.17504/protocols.io.erubd6w)

**Published:** 13 Jun 2016

## Guidelines

### FIGURE OUT WHAT YOU NEED:

A. Type of gel: regular agarose or low melting point (LMP) agarose

B. Concentration of gel

1. Mini gels (50 mL gels):

0.7% = 0.35 gm per 50 mL

0.8% = 0.40 gm per 50 mL

1.2% = 0.60 gm per 50 mL

2. Small gels (100 mL gels):

0.7% = 0.70 gm per 100 mL

0.8% = 0.80 gm per 100 mL

1.2% = 1.20 gm per 100 mL

3. Large gels (250 mL gels):

0.7% = 1.75 gm per 250 mL

0.8% = 2.00 gm per 250 mL

1.2% = 3.00 gm per 250 mL

The concentration of the gel will depend on the size fragments you are working with.

Low concentration gels are used for large DNA fragment and high concentration gels are used for small DNA fragments.

C. Type of buffer

1. TPE (Tris/Phosphate/EDTA), 10X stock

2. TBE (Tris/Borate/EDTA), 10X stock

3. TAE (Tris/Acetate/EDTA), 10X stock

TPE buffer is a general use buffer. TBE and TAE buffers are generally used when electroeluting fragments from gels.

D. Size of gel

1. Mini - 50 ml gel + 250 ml running buffer

2. Small - 100 ml gel + 900 ml running buffer

3. Large - 250 ml gel + 2000 ml running buffer (slope end gel boxes) or 2500 ml running buffer

(square end gel boxes)

E. Gel combs

1. Mini - 8 lane combs (to 25 µl each well)
2. Small - 14 lane or 20 lane combs (to 25 µl each well)
3. Large - 20 lane or 30 lane combs (to 20 µl each well)

**GEL BUFFERS AND ELECTROPHORESIS MATERIALS:**

1) TPE buffer, 2 liters, 10X

216.0 gm of Tris base  
30.2 mL of 85% Phosphoric acid  
80.0 mL of 500 mM EDTA, pH 8.0

2) TBE buffer, 2 liters, 10X

216.0 gm of Tris base  
110.0 gm of Boric acid  
100.0 mL of 500 mM EDTA, pH 8.0

Filter the 10X buffer through Whatman #1 filter paper before placing the buffer in the bottle.

3) TAE buffer, 2 liters, 10X

193.76 gm of Tris base  
108.86 gm of Sodium acetate  
26.9 gm of disodium EDTA

Adjust the pH to 8.0 with glacial acetic acid.

**Protocol**

Gel Preparation and Electrophoresis

**Step 1.**

Place the agarose in an appropriately sized flask (usually a flask that holds at least twice the volume of the gel you are preparing) with a stir bar, buffer, and d-H<sub>2</sub>O.

Gel Preparation and Electrophoresis

**Step 2.**

Melt the agarose in a microwave oven, autoclave, or on a hot plate/stirrer.

Gel Preparation and Electrophoresis

**Step 3.**

When the agarose has been melted, place the flask on a stir plate and gently stir the contents until the agarose cools to 50-60°C.

🔗 **NOTES**

**Irina Agarkova** 25 Mar 2016

The flask should be warm, but not hot, when holding the flask.

Gel Preparation and Electrophoresis

**Step 4.**

Pour the gel at room temperature for regular agarose or in the cold room for low melting point agarose.

#### Gel Preparation and Electrophoresis

##### **Step 5.**

Tape the ends of the gel tray.

#### Gel Preparation and Electrophoresis

##### **Step 6.**

Place the comb in the tray.

#### Gel Preparation and Electrophoresis

##### **Step 7.**

Level the gel tray so that the gel will solidify evenly.

#### Gel Preparation and Electrophoresis

##### **Step 8.**

Let the gel solidify (30 min).



**DURATION**

00:30:00

#### Gel Preparation and Electrophoresis

##### **Step 9.**

Remove the tape.

#### Gel Preparation and Electrophoresis

##### **Step 10.**

Pour the running buffer into the gel apparatus.

#### Gel Preparation and Electrophoresis

##### **Step 11.**

Place the gel in the apparatus.

#### Gel Preparation and Electrophoresis

##### **Step 12.**

Remove the comb from the gel.



**NOTES**

**Irina Agarkova** 25 Mar 2016

The comb is easier removed from the gel when the gel is submerged in buffer; it is possible to tear the wells trying to remove the comb before the gel is submerged in buffer.

#### Gel Preparation and Electrophoresis

##### **Step 13.**

Load the gel.

#### Gel Preparation and Electrophoresis

##### **Step 14.**

Connect the leads from the gel apparatus to the power supply.

#### Gel Preparation and Electrophoresis

##### **Step 15.**

Electrophorese the samples toward the positive (red) electrode.



**NOTES**

**Irina Agarkova** 25 Mar 2016

Mini gels are generally run for 100-150 volt-hours, small gels are generally run for 300-350 volt-

hours, and the large gels are generally run for 900-1000 volt-hours. These are general conditions and may be changed for particular experiments.

#### Gel Preparation and Electrophoresis

##### **Step 16.**

Stain the gels with 0.5 µg/mL of ethidium bromide for 30 min.

🕒 **DURATION**

00:30:00

#### Gel Preparation and Electrophoresis

##### **Step 17.**

Stain the regular gels at room temperature and the low melting point gels at 4°C.

#### Gel Preparation and Electrophoresis

##### **Step 18.**

Destain the gels with d-H<sub>2</sub>O and view.