

Lox WT gene PCR: amplification from genomic DNA

Miguel Torres, Tara deBoer

Abstract

This protocol is a standard polymerase chain reaction, optimized to amplify a () bp gene flanked by two loxP sites. (?)

Refer to page 32, April 5, 2017 entry in Miguel's laboratory notebook for the original document.

Citation: Miguel Torres, Tara deBoer Lox WT gene PCR: amplification from genomic DNA. **protocols.io**
dx.doi.org/10.17504/protocols.io.hnwb5fe

Published: 14 Apr 2017

Materials

Phusion high-fidelity PCR kit [F553S](#) by [Thermo Scientific](#)

Protocol

PCR reaction

Step 1.

Into a 2.5 mL Eppendorf® tube add:

50 uL H₂O

16 uL GC buffer

1.6 uL dNTPs solution

4 uL 10 uM Forward primer

4 uL 20 uM Reverse primer

1.5 uL DNA to amplify (lox)

0.8 uL Phusion polymerase

Total=80 uL

Add 20 uL into four 0.2 mL tubes

--Ready for thermal cycler (vid. step 2).

Thermal cycler protocol

Step 2.

95°C, 5 min

95°C, 30 seconds

60°C 30 seconds

72°C, 30 seconds

30 X

72°C, 5 min

4°C infinite hold

Result example

Step 3.