

Making AA+ medium for Synechococcus sp. PCC 7002

Anne Vogel, Martin Hohmann-Marriott, Rahmi Lale

Abstract

This protocol describes how to make growth medium for Synechococcus sp. PCC 7002, including the stock solutions needed.

Citation: Anne Vogel, Martin Hohmann-Marriott, Rahmi Lale Making AA+ medium for Synechococcus sp. PCC 7002. protocols.io

dx.doi.org/10.17504/protocols.io.kg3ctyn

Published: 16 Jul 2018

Protocol

Make stock solutions

Step 1.

Make 1000x FeCl₂·6H₂O

- Add 0.405g FeCl₂·6H₂O to 100mL MilliQ
- Filtersterilize. **Do not autoclave!**

Step 2.

Make 100x Tris 8.2

- 9.99g Tris in 100mL MilliQ
- Adjust pH to 8.2 with HCl
- Autoclave

Step 3.

Make 1000 x Trace mineral solution

- Use the following compounds:

Compound	Amount (g/L)
H ₃ BO ₃	2.86
MnCl ₂ x 4 H ₂ O	1.81
ZnSO ₄ x 7 H ₂ O	0.222
Na ₂ MoO ₄ x 2 H ₂ O	0.390
CuSO ₄ x 5 H ₂ O	0.079
Co(NO ₃) ₂ x 6 H ₂ O	0.0494

- Autoclave or filtersterilize

Step 4.

Make 100x AA+ solution

- Use the following compounds, autoclave

Compound	500mL	250mL	100mL
KH ₂ PO ₄	1.25g	0.625g	0.25g
Na ₂ -EDTA x 2 H ₂ O	1.5g	0.750g	0.3g

Compound	500mL	250mL	100mL
KCl	30g	15g	6g
NaNO ₃	50g	25g	10g
CaCl ₂	6.65g	3.325g	1.33g

Make the AA+ medium

Step 5.

- Add compounds as described in table below to desired amount of miliQ.
- Add glycerol only when non-phototrophic growth is needed

Compound	1L	500mL	250mL
NaCl	18g	9g	4.5g
MgSO ₄ ·7H ₂ O	5g	2.5g	1.25g
Glycerol	1.1mL	550μL	275μL

- If you want to make plates, add 15g/L agar agar
- Autoclave

Step 6.

- Add compounds (table below) in the flow cabinet to the autoclaved and *hand-warm* solution
- Only add the B₁₂ if you are immediately using the medium. If some medium with B₁₂ is left after you used it, store it in a dark, cold place for a short time

Compound	1L	500mL	250mL
100x Tris pH 8.2	10mL	5mL	2.5mL
1000x FeCl ₂ ·6H ₂ O	1mL	500μL	250μL
1000x Trace mineral solution	1mL	500μL	250μL
100x AA+ solution	10mL	5mL	2.5mL
vit B12 (1mg/mL)	4μL	2μL	1μL