



Apr 15, 2019

Working

DNA isolation from Formalin-Fixed, Paraffin-Embedded (FFPE) material

Hendrik F van Essen¹¹Amsterdam UMC

dx.doi.org/10.17504/protocols.io.zjwf4pe

Hendrik F HFvE. van Essen
Amsterdam UMC

ABSTRACT

This protocol describes the isolation of DNA from paraffin embedded tissue. The resulting DNA can be used for CGH array hybridizations.

After removal of the paraffin, the tissue is lysed and the DNA is released from the cells. Then the DNA is isolated by using a silica-based column that binds the DNA.

STEPS MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Buffer ATL	19076	Qiagen
Proteinase K	E00491	Thermo Fisher Scientific
Buffer AL, Lysis buffer	19076	Qiagen
EtOH		
Buffer AW1	19081	Qiagen
Buffer AW2	19072	Qiagen
Buffer AE	19077	Qiagen

H&E staining from paraffin blocks

- 1 Cut 3 µm (1x), 10 µm (2-3x, or more if needed depending on tumour area)
- 2 Mount all the section on slides with BSA 0.1%, and dry them o/n at 37°C (to a maximum of one week)
- 3 Place slides in xylene ⌚ 00:07:00 3 times
- 4 Hydrate by passing sequentially through 100%, 96%, 70% EtOH, and distilled-water
- 5 Stain slides with haematoxylin ⌚ 00:02:00
(if solution is not fresh staining time can be extended, max 4 min.)
- 6 Wash with running tap water ⌚ 00:05:00

- 7 Stain slides with eosin 🕒 00:02:00
- 8 Wash quickly in water
- 9 Dehydrate by passing sequentially through 70%, 96%, 100% EtOH, xylene
- 10 Remove slides from xylene and use DePex to mount a covering glass.
- 11 Let the 3 µm H&E slide be judged by the pathologist for dissection.
Mark the tumor area and tumor percentage

DNA isolation

- 12 Place slides in xylene 3 x 7 min. 🕒 00:07:00 3 times
- 13 Hydrate by passing sequentially through 100%, 96%, 70% EtOH, and water
- 14 Stain slides with haematoxylin 1-2 min. 🕒 00:02:00
(if solution is not fresh staining time can be extended, max 3 min.)
- 15 Wash with running tap water 🕒 00:05:00
- 16 Place slides in distilled water
- 17 Dissect the tissue from the slides (when they are still slightly wet) with a scalpel or a needle
- 18 Place the material in a clean safelock eppendorf cup
- 19 Spin down in a centrifuge at full speed and get rid of the water layer 🕒 00:05:00
- 20 Shake and invert samples (do not vortex).
- 21 Add 160 µl ATL buffer



Buffer ATL
by [Qiagen](#)
Catalog #: 19076

22 Add 40 µl prot K (20 mg/mL)



Proteinase K
by [Thermo Fisher Scientific](#)
Catalog #: E00491

23 vortex sample ⌚ 00:00:15
Make sure all tissue is in the liquid

24 Incubate o/n at 56°C in heat-block or waterbath and vortex regularly! ⌚ 08:00:00 🌡 56 °C
Check if all tissue is digested

25 Vortex samples and spin down

26 Incubate at 98°C 10 min. 🌡 98 °C ⌚ 00:10:00

27 Spin down at 16600 rcf. 1 min. ⌚ 00:01:00

28 Add 200 µl AL buffer and mixvery well by vortexing!




Buffer AL, Lysis buffer
by [Qiagen](#)
Catalog #: 19076

29 Add 200 µl Ethanol 100% and mix very well by vortexing!



EtOH


- 30 Incubate at RT ⌚ 00:05:00
- 31 Spin down the eppendorf cups at full speed 1 min.
- 32 Transfer up to 600 µl lysate to the QIAamp MinElute Column.
- 33 Spin down at 8.000 rpm ⌚ 00:01:00
- 34 Place the QIA Column in a clean 2 ml collection tube, discard the flow-through make sure the column tip is clear of droplets
- 35 Add 500 µl AW1 buffer to the column
- 36 Spin down for at 8.000 rpm and place the column in a new tube ⌚ 00:01:00
- 37 Add 500 µl AW2 buffer (QIAamp micro-kit) to the column
- 38 Spin down at 8.000 rpm and place the column in a new tube ⌚ 00:01:00
- 39 Spin down at full speed to dry the membrane ⌚ 00:03:00
- 40 Place the QI column in a properly marked eppendorf cup, discard collection tube
- 41 Add 20 to 50 µl (depending on sample size) of AE buffer to the column and incubate. ⌚ 00:05:00 🌡 20 °C room temperature



Buffer AW1

by Qiagen


Catalog #: 19081



Buffer AW2

by Qiagen

Catalog #: 19072



Buffer AE

by Qiagen

Catalog #: 19077

42 Spin down at full speed. ⌚ 00:03:00

43 Throw away the column, close the eppendorf cup and store the DNA at 🧊 4 °C



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited