

Immunogold labeling of bacterial cells for transmission electron microscopy (TEM)

Jessica Sacher

Abstract

Citation: Jessica Sacher Immunogold labeling of bacterial cells for transmission electron microscopy (TEM).
protocols.io

dx.doi.org/10.17504/protocols.io.mv4c68w

Published: 30 Jan 2018

Materials

- ✓ 1X PBS (Phosphate-buffered saline) by Contributed by users
- ✓ 70% Ethanol by Contributed by users
- Paraformaldehyde Aqueous Solution -16% [15700](#) by [Electron Microscopy Sciences](#)
- Anti-Rabbit IgG (whole molecule)-Gold antibody produced in goat [Sigma 0162787](#) by [Sigma](#)
- ✓ 5% BSA in PBST by Contributed by users
- ✓ Formvar-coated 300-mesh electron microscopy grids by Contributed by users

Protocol

Step 1.

Cut a strip of Parafilm (20-30 cm long), place flat on bench.

Step 2.

Prepare one 1-mL tube each of 70% ethanol and mqH2O for washing tweezers between transfers.

Step 3.

Harvest cells by flooding an overnight plate with 5 mL nutrient broth or buffer, using a spreader to dislodge cells, pipetting into a conical flask and mixing gently by inversion (avoid pipetting up and down, vortexing and centrifuging to preserve flagella).

Step 4.

Set cells to an OD600 of 1-3 in NZCYM broth (1-mL volume), place on ice until use.

Step 5.

Pipette 100 uL cells onto parafilm, use tweezers to float Formvar-coated copper (300 mesh) grid atop droplet, Formvar-side-down (should do this in duplicate, ie. 2 grids per drop), incubate 45 min at RT (between transferring grids, particularly between transfer of grids containing different samples, swirl tweezers 2-3 seconds in EtOH, then in H2O, then wipe dry with a fresh Kimwipe)

Step 6.

Wash grids 1x by transferring to a 100-uL drop of PBS, incubating 3-5 min.

Step 7.

Incubate 20 min in 2.5% paraformaldehyde/PBS to kill cells and fix them to grids.

Step 8.

Wash 3x in PBS (100 uL drops, 3-5 min each).

Step 9.

Block in 5% BSA/PBST (100 uL, 35 min).

Step 10.

*If doing Far Western blot, ie. testing protein binding to cells, incubate in desired protein diluted 1/25 in blocking solution, 45 min, 100 uL drop

Step 11.

Wash 3x in blocking solution (100 uL drops, 3-5 min each).

Step 12.

Incubate in primary antibody (anti-Gp047, rabbit) at 1/50 (diluted in blocking solution), 45 min, 100 uL.

Step 13.

Wash 3x in blocking solution (100 uL drops, 3-5 min each).

Step 14.

Incubate in secondary antibody (goat anti-rabbit IgG-gold) at 1/50 (diluted in blocking solution), 45 min, 100 uL.

Step 15.

Wash 3x in blocking solution (100 uL drops, 3-5 min each).

Step 16.

Wash 3x in PBS (100 uL drops, 3-5 min each).

Step 17.

Wash 3x in water (100 uL drops, 3-5 min each) to remove salts, which can form crystals under EM.

Step 18.

Dry completely on Whatman filter paper in a Petri dish (at least 20 min, preferably overnight). OK to store grids in this way at RT for several days prior to TEM, but ideally visualize within 1 day.

Step 19.

Analyze grids via TEM.