

Nov 06, 2019

Transient transformation of Ostreococcus species (OTTH595, RCC809 and RCC802) and Bathycoccus V.3

francois-yves bouget1, Francois-Yves Bouget1

¹Laboratoire d'Oceanographie Microbienne



Protist Research to Optimize Tools in Genetics (PROT-G)



ABSTRACT


This protocol describes the preparation of cells and introduction of DNA into the cells by electroporation. For selection of stable transformants or measure of transient gene expression see related protocols.

MATERIALS TEXT

Ostreococcus lucimarinus (RCC802)

pHAPT:Luc vector (Djouani Tahri et al., PLOS ONE 2011 https://doi.org/10.1371/journal.pone.0028471)

pHAPT:Luc map



1

PHAPT:Luc sequence

O. tauri pHAPT:Luc6240 bp ds-DNAcircular DEFINITIONsynthetic circular DNA

FEATURESLocation/Qualifiers

source1..6240

/organism="synthetic DNA construct"

/mol_type="other DNA"

CDScomplement(377..1237)

/codon_start=1

/gene="bla"

/product="beta-lactamase"

/label=AmpR

/note="confers resistance to ampicillin, carbenicillin, and

related antibiotics"

/ translation = "MSIQHFRVALIPFFAAFCLPVFAHPETLVKVKDAEDQLGARVGYI"

ELDLNSGKILESFRPEERFPMMSTFKVLLCGAVLSRIDAGQEQLGRRIHYSQNDLVEYS

PVTEKHLTDGMTVRELCSAAITMSDNTAANLLLTTIGGPKELTAFLHNMGDHVTRLDRW

EPELNEAIPNDERDTTMPVAMATTLRKLLTGELLTLASRQQLIDWMEADKVAGPLLRSA

LPAGWFIADKSGAGERGSRGIIAALGPDGKPSRIVVIYTTGSOATMDERNROIAEIGAS

LIKHW"

promotercomplement(1238..1342)

/gene="bla"

/label=AmpR promoter

primer_bind1816..1832

/label=M13 fwd

/note="common sequencing primer, one of multiple similar

variants"

promoter1883..2002

/label=Ostreococcus tauri High affinity phosphate transporte promoter

CDS2005..3657

/codon_start=1

/gene="luc+"

/product="firefly luciferase"

/label=luciferase

/note="enhanced luc+ version of the luciferase gene"

/translation = "MEDAKNIKKGPAPFYPLEDGTAGEQLHKAMKRYALVPGTIAFTDA"

 ${\sf HIEVDITYAEYFEMSVRLAEAMKRYGLNTNHRIVVCSENSLQFFMPVLGALFIGVAVAP}$

ANDIYNERELLNSMGISQPTVVFVSKKGLQKILNVQKKLPIIQKIIIMDSKTDYQGFQS

 ${\tt MYTFVTSHLPPGFNEYDFVPESFDRDKTIALIMNSSGSTGLPKGVALPHRTACVRFSHAMMS} \\$

 ${\tt RDPIFGNQIIPDTAILSVVPFHHGFGMFTTLGYLICGFRVVLMYRFEEELFLRSLQDYK}$

IQSALLVPTLFSFFAKSTLIDKYDLSNLHEIASGGAPLSKEVGEAVAKRFHLPGIRQGY

GLTETTSAILITPEGDDKPGAVGKVVPFFEAKVVDLDTGKTLGVNQRGELCVRGPMIMS

 ${\tt GYVNNPEATNALIDKDGWLHSGDIAYWDEDEHFFIVDRLKSLIKYKGYQVAPAELESIL}$

 ${\tt LQHPNIFDAGVAGLPDDDAGELPAAVVVLEHGKTMTEKEIVDYVASQVTTAKKLRGGVV}$

FVDEVPKGLTGKLDARKIREILIKAKKGGKIAV"

terminator3820..4017

/label=TEF terminator

/note="Ashbya gossypii TEF terminator"

CDScomplement(4026..4835)

/codon_start=1

/gene="aph(3')-la"

/product="aminoglycoside phosphotransferase"

/label=KanR

/note="confers resistance to kanamycin"

/translation="MGKEKTHVSRPRLNSNMDADLYGYKWARDNVGQSGATIYRLYGKP

DAPELFLKHGKGSVANDVTDEMVRLNWLTEFMPLPTIKHFIRTPDDAWLLTTAIPGKTA FOVI FEYPDSGENIVDALAVELRRI HSIPVCNCPENSDRVERI AOAOSRMNNGI VDASD FDDERNGWPVEQVWKEMHKLLPFSPDSVVTHGDFSLDNLIFDEGKLIGCIDVGRVGIAD RYQDLAILWNCLGEFSPSLQKRLFQKYGIDNPDMNKLQFHLMLDEFF" primer_bindcomplement(5456..5472) /label=M13 rev /note="common sequencing primer, one of multiple similar variants' protein_bind5480..5496 /label=lac operator /bound_moiety="lac repressor encoded by lac!" /note="The lac repressor binds to the lac operator to inhibit transcription in E. coli. This inhibition can be relieved by adding lactose or isopropyl-beta-D-thiogalactopyranoside (IPTG)." promotercomplement(5504..5534) /label=lac promoter /note="promoter for the E. coli lac operon" protein_bind5549..5570 /label=CAP binding site /bound_moiety="E. coli catabolite activator protein" /note="CAP binding activates transcription in the presence of cAMP." rep_origincomplement(join(5858..6240,1..206)) /direction=LEFT /label=ori /note="high-copy-number ColE1/pMB1/pBR322/pUC origin of replication" **ORIGIN** 1 atgtaggegg tgctacagag ttcttgaagt ggtggcctaa ctacggctac actagaagga 61 cagtatttgg tatctgcgct ctgctgaagc cagttacctt cggaaaaaga gttggtagct 121 cttgatccgg caaacaaacc accgctggta gcggtggttt ttttgtttgc aagcagcaga 181 ttacgcgcag aaaaaaagga tctcaagaag atcctttgat cttttctacg gggtctgacg 241 ctcagtggaa cgaaaactca cgttaaggga ttttggtcat gagattatca aaaaggatct 301 tcacctagat ccttttaaat taaaaatgaa gttttaaatc aatctaaagt atatatgagt 361 aaacttggtc tgacagttac caatgcttaa tcagtgaggc acctatctca gcgatctgtc 421 tatttcgttc atccatagtt gcctgactcc ccgtcgtgta gataactacg atacgggagg 481 gettaceate tggeceeagt getgeaatga tacegegaga eccaegetea eeggeteeag 541 atttatcagc aataaaccag ccagccggaa gggccgagcg cagaagtggt cctgcaactt 601 tatccgcctc catccagtct attaattgtt gccgggaagc tagagtaagt agttcgccag 661 ttaatagttt gegeaaegtt gttgeeattg etaeaggeat egtggtgtea egetegtegt 721 ttggtatggc ttcattcagc tccggttccc aacgatcaag gcgagttaca tgatccccca 781 tgttgtgcaa aaaagcggtt agctccttcg gtcctccgat cgttgtcaga agtaagttgg 841 ccgcagtgtt atcactcatg gttatggcag cactgcataa ttctcttact gtcatgccat 901 ccgtaagatg cttttctgtg actggtgagt actcaaccaa gtcattctga gaatagtgta 961 tgcggcgacc gagttgctct tgcccggcgt caatacggga taataccgcg ccacatagca 1021 gaactttaaa agtgctcatc attggaaaac gttcttcggg gcgaaaactc tcaaggatct 1081 taccgctgtt gagatccagt tcgatgtaac ccactcgtgc acccaactga tcttcagcat 1141 cttttacttt caccagcgtt tctgggtgag caaaaacagg aaggcaaaat gccgcaaaaa 1201 agggaataag ggcgacacgg aaatgttgaa tactcatact cttccttttt caatattatt 1261 gaagcattta tcagggttat tgtctcatga gcggatacat atttgaatgt atttagaaaa 1321 ataaacaaat aggggttccg cgcacatttc cccgaaaagt gccacctgac gtctaagaaa 1381 ccattattat catgacatta acctataaaa ataggcgtat cacgaggccc tttcgtctcg

1441 cgcgtttcgg tgatgacggt gaaaacctct gacacatgca gctcccggag acggtcacag 1501 cttgtctgta agcggatgcc gggagcagac aagcccgtca gggcgcgtca gcgggtgttg 1561 gcgggtgtcg gggctggctt aactatgcgg catcagagca gattgtactg agagtgcacc 1621 atatgcggtg tgaaataccg cacagatgcg taaggagaaa ataccgcatc aggcgccatt 1681 cgccattcag gctgcgcaac tgttgggaag ggcgatcggt gcgggcctct tcgctattac

to protocols.io 3 11/06/2019

1741 gccagctggc gaaaggggga tgtgctgcaa ggcgattaag ttgggtaacg ccagggtttt 1801 cccagtcacg acgttgtaaa acgacggcca gtgaattcga gctcggtacc cttcatctcc 1861 ggttctgctg ctagcaagat ctaaatactt ttccaggaaa agcaggaaaa aggaatattc 1921 tttaatattc tcggaatatt tcacgttctt gcggcctggc gacgctgacg gcgttacttc 1981 acgtcgcgct ttacccgtct cgccatggaa gacgccaaaa acataaagaa aggcccggcg 2041 ccattctatc cgctggaaga tggaaccgct ggagagcaac tgcataaggc tatgaagaga 2101 tacgccctgg ttcctggaac aattgctttt acagatgcac atatcgaggt ggacatcact 2161 tacgctgagt acttcgaaat gtccgttcgg ttggcagaag ctatgaaacg atatgggctg 2221 aatacaaatc acagaatcgt cgtatgcagt gaaaactctc ttcaattctt tatgccggtg 2281 ttgggcgcgt tatttatcgg agttgcagtt gcgcccgcga acgacattta taatgaacgt 2341 gaattgctca acagtatggg catttcgcag cctaccgtgg tgttcgtttc caaaaagggg 2401 ttgcaaaaaa ttttgaacgt gcaaaaaaag ctcccaatca tccaaaaaat tattatcatg 2461 gattetaaaa eggattacea gggattteag tegatgtaca egttegteac ateteateta 2521 cctcccggtt ttaatgaata cgattttgtg ccagagtcct tcgataggga caagacaatt 2581 geactgatea tgaacteete tggatetaet ggtetgeeta aaggtgtege tetgeeteat 2641 agaactgcct gcgtgagatt ctcgcatgcc agagatccta tttttggcaa tcaaatcatt 2701 ccggatactg cgattttaag tgttgttcca ttccatcacg gttttggaat gtttactaca 2761 ctcggatatt tgatatgtgg atttcgagtc gtcttaatgt atagatttga agaagagctg 2821 tttctgagga gccttcagga ttacaagatt caaagtgcgc tgctggtgcc aaccctattc 2881 tccttcttcg ccaaaagcac tctgattgac aaatacgatt tatctaattt acacgaaatt 2941 gcttctggtg gcgctcccct ctctaaggaa gtcggggaag cggttgccaa gaggttccat 3001 ctgccaggta tcaggcaagg atatgggctc actgagacta catcagctat tctgattaca 3061 cccgaggggg atgataaacc gggcgcggtc ggtaaagttg ttccattttt tgaagcgaag 3121 gttgtggatc tggataccgg gaaaacgctg ggcgttaatc aaagaggcga actgtgtgtg 3181 agaggtccta tgattatgtc cggttatgta aacaatccgg aagcgaccaa cgccttgatt 3241 gacaaggatg gatggctaca ttctggagac atagcttact gggacgaaga cgaacacttc 3301 ttcatcgttg accgcctgaa gtctctgatt aagtacaaag gctatcaggt ggctcccgct 3361 gaattggaat coatcttgct coaacacccc aacatcttcg acgcaggtgt cgcaggtctt 3421 cccgacgatg acgccggtga acttcccgcc gccgttgttg ttttggagca cggaaagacg 3481 atgacqqaaa aaqaqatcqt qqattacqtc qccaqtcaaq taacaaccqc qaaaaaqttq 3541 cgcggaggag ttgtgtttgt ggacgaagta ccgaaaggtc ttaccggaaa actcgacgca 3601 agaaaaatca gagagatcct cataaaggcc aagaaggcg gaaagatcgc cgtgtaattc 3661 tagctaaaaa gtggaacgat cattcactat atatatatca atttatatat acgtatgtgt 3721 aattgaagaa agatacgttt tttcctctat tgagaggatc atcgatgaat tcgagctcgt 3781 tttcgacact ggatggcggc gttagtatcg aatcgacagc agtatagcga ccagcattca 3841 catacgattg acgcatgata ttactttctg cgcacttaac ttcgcatctg ggcagatgat 3901 gtcgaggcga aaaaaaatat aaatcacgct aacatttgat taaaatagaa caactacaat 3961 ataaaaaaac tatacaaatg acaagttett gaaaacaaga atetttttat tgtcagtact 4021 gattattaga aaaactcatc gagcatcaaa tgaaactgca atttattcat atcaggatta 4081 tcaataccat atttttgaaa aagccgtttc tgtaatgaag gagaaaactc accgaggcag 4141 ttccatagga tggcaagatc ctggtatcgg tctgcgattc cgactcgtcc aacatcaata 4201 caacctatta atttcccctc gtcaaaaata aggttatcaa gtgagaaatc accatgagtg 4261 acgactgaat ccggtgagaa tggcaaaagc ttatgcattt ctttccagac ttgttcaaca 4321 ggccagccat tacgctcgtc atcaaaatca ctcgcatcaa ccaaaccgtt attcattcgt 4381 gattgcgcct gagcgagacg aaatacgcga tcgctgttaa aaggacaatt acaaacagga 4441 atcgaatgca accggcgcag gaacactgcc agcgcatcaa caatattttc acctgaatca 4501 ggatattett etaatacetg gaatgetgtt ttgeegggga tegeagtggt gagtaaceat 4561 gcatcatcag gagtacggat aaaatgcttg atggtcggaa gaggcataaa ttccgtcagc 4621 cagtttagtc tgaccatctc atctgtaaca tcattggcaa cgctaccttt gccatgtttc 4681 agaaacaact ctggcgcatc gggcttccca tacaatcgat agattgtcgc acctgattgc 4741 ccgacattat cgcgagccca tttataccca tataaatcag catccatgtt ggaatttaat 4801 cgcggcctcg aaacgtgagt cttttcctta cccatgcatg gttgttttcg aacgtgtcga 4861 ttagtgacga cctgcgcgcg cggtcggggc cacggacgcg ggcgccggtg agttcgcggt 4921 gagttgaaaa tactccgtgc gagtcacgac ccatgaaatt ataggggtaa aaagttcatt 4981 tcacacaatc tttgatgcgt tgattgtttc taaatggcac acaatgagtc atgaaccaca 5041 atgcgtggat ttaaaaacaa tcaagtgatg catgcacgtg aaaaatgggc atttcaagta 5101 aatagagagt aatgactcgt gtatgacccc ctgctttcaa actcaccgac tcaccggcgc 5161 gegeggtgaa aegetegatt egategeaeg gegtegteae aactegaeag teaatatgge

5221 ccgtacaaag caaactgccc gtaagtccac cgggggcaag gccccgcgca agcaactcgc 5281 caccaaggeg geeegeaagt eegeeegge gaeeggegge gteaagaage egeaeaggta 5341 ccgcccgggt accgttgcgc tcgtgggatc tttcacctct tagtcttctg gcctcgaggt 5401 cgacggggga tcctctagag tcgacctgca ggcatgcaag cttggcgtaa tcatggtcat 5461 agctgtttcc tgtgtgaaat tgttatccgc tcacaattcc acacaacata cgagccggaa 5521 gcataaagtg taaagcctgg ggtgcctaat gagtgagcta actcacatta attgcgttgc 5581 geteactgee egettteeag tegggaaace tgtegtgeea getgeattaa tgaateggee 5641 aacgcgcggg gagaggcggt ttgcgtattg ggcgctcttc cgcttcctcg ctcactgact 5701 cgctgcgctc ggtcgttcgg ctgcggcgag cggtatcagc tcactcaaag gcggtaatac 5761 ggttatccac agaatcaggg gataacgcag gaaagaacat gtgagcaaaa ggccagcaaa 5821 aggccaggaa ccgtaaaaag gccgcgttgc tggcgttttt ccataggctc cgccccctg 5881 acgagcatca caaaaatcga cgctcaagtc agaggtggcg aaacccgaca ggactataaa 5941 gataccagge gttteeceet ggaageteec tegtgegete teetgtteeg accetgeege 6001 ttaccggata cetgtccgcc tttctccctt cgggaagcgt ggcgctttct caatgctcac 6061 getgtaggta teteagtteg gtgtaggteg ttegeteeaa getgggetgt gtgeacgaac 6121 cccccgttca gcccgaccgc tgcgccttat ccggtaacta tcgtcttgag tccaacccgg 6181 taagacacga cttatcgcca ctggcagcag ccactggtaa caggattagc agagcgaggt

Validation of O. lucimarinus transgenic lines by PCR (Figure 4A)

Forward primer: 1GAGCGCAACGGTACCCGGGCGGTACCTGTGCG Reverse primer: CTGGCGACGCTGACGGCGTTACTTCACGT

Bathycoccus prasinos (RCC4222)

pH4:KanMx pHAPT:Luc transgene (PCR product)

Construction of the transgene by PCR:

The transgene was generated by fusion PCR as described in Shevchuk et al. (2004) using DNA templates described below and oligonucleotides mapped on the sequence below.

Promoters of Histone H4 and high affinity phosphate (HAPT) transporters were amplified using (Fas1/Ras1) and (Fas4/Ras4) oligonucucleotides respectively. KanMx and Luciferase sequences were amplified from pHAPT: luc plasmid template using oligonucleotides (Fas2/Ras2) and (Fas3/Ras3). Final PCR was done on amplified fragments using (Fas5/Ras5) oligonucleotides.

Oligonucleotide sequences:

Fas1:CCGGCTTCGTGATGCCTTGGATGTTCTC

Ras1: TCGAAACGTGAGTCTTTTCCTTACCcatTGTGTTTGATTTATAATGAGGTTTTCTT

Fas2: AAGAAAACCTCATTATAAATCAAACACAatgGGTAAGGAAAAGACTCACGTTTCGA

Ras2: CCAAGAAGGGCGGAAAGATCGCCGTGTAAGAAATACCGTCTATCATCGATGAATTCGA


Fas3: CGAATTCATCGATGATAGACGGTATTTCTTACACGGCGATCTTTCCGCCCTTCTTGG

Ras4: AAGTGTACACGACCAAACGCGCGTCGAC

Fas5: CCACCTTTACCTCTGCCGGACATTGTGA

Ras5: GCGCGGTAATATCTACGAGGTAGCACGAG

Map of the transgene:



Sequence of the transgene:

4020 bp ds-DNAlinear

REFERENCE1 (bases 1 to 4020)
AUTHORSFY Bouget and JC Lozano

FEATURESLocation/Qualifiers

source1..4020

/organism="unspecified"

/mol_type="genomic DNA"

primer_bind1..30

/label=Fas1

primer_bind78..105

/label=Fas5

promoter101..805

/label=Bprasinos H4 p

primer_bindcomplement(778..833)

/label=Ras1

primer_bind778..833

/label=Fas2

CDS806..1615

/codon_start=1

/gene="aph(3')-la"

/product="aminoglycoside phosphotransferase"

/label=KanR

/note="confers resistance to kanamycin"

/translation="MGKEKTHVSRPRLNSNMDADLYGYKWARDNVGQSGATIYRLYGKP

DAPELFLKHGKGSVANDVTDEMVRLNWLTEFMPLPTIKHFIRTPDDAWLLTTAIPGKTA

 ${\tt FQVLEEYPDSGENIVDALAVFLRRLHSIPVCNCPFNSDRVFRLAQAQSRMNNGLVDASD}$

FDDERNGWPVEQVWKEMHKLLPFSPDSVVTHGDFSLDNLIFDEGKLIGCIDVGRVGIAD

RYQDLAILWNCLGEFSPSLQKRLFQKYGIDNPDMNKLQFHLMLDEFF"

terminator1624..1821

/label=TEF terminator

/note="Ashbya gossypii TEF terminator"

primer_bindcomplement(1867..1924)

/label=Ras2

primer_bind1868..1924

/label=Fas 3

CDScomplement(1896..3548)

/codon_start=1

/gene="luc+"

/product="firefly luciferase"

/label=luciferase

/note="enhanced luc+ version of the luciferase gene"

/translation="MEDAKNIKKGPAPFYPLEDGTAGEQLHKAMKRYALVPGTIAFTDA HIEVDITYAEYFEMSVRLAEAMKRYGLNTNHRIVVCSENSLQFFMPVLGALFIGVAVAP ANDIYNERELLNSMGISQPTVVFVSKKGLQKILNVQKKLPIIQKIIIMDSKTDYQGFQS MYTFVTSHLPPGFNEYDFVPESFDRDKTIALIMNSSGSTGLPKGVALPHRTACVRFSHA RDPIFGNQIIPDTAILSVVPFHHGFGMFTTLGYLICGFRVVLMYRFEEELFLRSLQDYK IQSALLVPTLFSFFAKSTLIDKYDLSNLHEIASGGAPLSKEVGEAVAKRFHLPGIRQGY GLTETTSAILITPEGDDKPGAVGKVVPFFEAKVVDLDTGKTLGVNQRGELCVRGPMIMS GYVNNPEATNALIDKDGWLHSGDIAYWDEDEHFFIVDRLKSLIKYKGYQVAPAELESIL LQHPNIFDAGVAGLPDDDAGELPAAVVVLEHGKTMTEKEIVDYVASQVTTAKKLRGGVV FVDEVPKGLTGKLDARKIREILIKAKKGGKIAV"

I VDEVERGETGREDARRIKEILIRARRGGRIAV

primer_bindcomplement(3520..3577)

/label=Ras3

primer_bind3520..3577

/label=Fas4

PROMOTER3549..3920

/label=B. prasinos High affinity phosphate transporter

primer_bindcomplement(3906..3934)

/label=Ras5

primer_bindcomplement(3978..4005)

/label=Ras4

ORIGIN

1 coggettegt gatgeettgg atgttgtete teaacacett geggtgtett ttggegeeae

61 ctttacctaa acctttacca cctttacctc tgccggacat tgtgaatatg tgttatattt

121 gtttttaaat gacgatttga acggacgttt tcgttttgtt tttcgtgttt ttttgtgttt

181 tttgcctgtt ttgtaggggt tacttctcag ttttgctata aagagtgcat ttgatgcatc

241 aaaacgagcg tttgagaagg tcgcgagtga aaaagctacg aagacgaaat tcaatcaacg

301 acaagaaaat acaacaaaaa gaacaacctc gatccttatc acttggaacg cgtatcacgt

361 tttaaaagaa ataattatac aatataactc ctcgtaacca tctatcatcc ctatatctat

421 cgcatcagaa ctttttagct tttcatttgc tctccgagtt cctcgagcgt cgaatacgtg

481 cgactagttg tgcatgtaca accetaaatg atagtgtgca aagtgcacct ttttttcgca

541 aatgctcgaa aatgcaagac cttaaatatt cggcgacatt aagcgcgctg aaaggtacgt

601 tactqttcac qqatttqatc qcqcqtqaac ttqttcqata qatatttaac atatctqtqt

661 ctttgagaag ttttgattca tatgaatcag tgatacatct gatgtcgtac aaaacgttaa

721 atcttgattt ctaaccttct tttgttttcg aatggttcaa atgaatcgtt gacacgaaag

781 aaaacctcat tataaatcaa acacaatggg taaggaaaag actcacgttt cgaggccgcg

841 attaaattcc aacatggatg ctgatttata tgggtataaa tgggctcgcg ataatgtcgg

901 gcaatcaggt gcgacaatct atcgattgta tgggaagccc gatgcgccag agttgtttct

961 gaaacatggc aaaggtagcg ttgccaatga tgttacagat gagatggtca gactaaactg

1021 getgaeggaa tttatgeete tteegaeeat eaageatttt ateegtaete etgatgatge

1081 atggttactc accactgcga tccccggcaa aacagcattc caggtattag aagaatatcc

1141 tgattcaggt gaaaatattg ttgatgcgct ggcagtgttc ctgcgccggt tgcattcgat

1201 tcctgtttgt aattgtcctt ttaacagcga tcgcgtattt cgtctcgctc aggcgcaatc

1261 acgaatgaat aacggtttgg ttgatgcgag tgattttgat gacgagcgta atggctggcc

1321 tgttgaacaa gtctggaaag aaatgcataa gcttttgcca ttctcaccgg attcagtcgt

1381 cactcatggt gatttctcac ttgataacct tatttttgac gaggggaaat taataggttg

1441 tattgatgtt ggacgagtcg gaatcgcaga ccgataccag gatcttgcca tcctatggaa

1501 ctgcctcggt gagttttctc cttcattaca gaaacggctt tttcaaaaat atggtattga

1561 taatcctgat atgaataaat tgcagtttca tttgatgctc gatgagtttt tctaataatc

1621 agtactgaca ataaaaagat tettgtttte aagaacttgt catttgtata gtttttttat 1681 attgtagttg ttetatttta ateaaatgtt agegtgattt atatttttt tegeetegae

106 Fattytagtty tictatitta atcadatytt agcytgatti atattitti tegeciegae

1741 atcatctgcc cagatgcgaa gttaagtgcg cagaaagtaa tatcatgcgt caatcgtatg

1801 tgaatgctgg tegetatact getgtegatt egatactaae geegeeatee agtgtegaaa

1861 acgagetega atteategat gatagaeggt atttettaca eggegatett teegeeette

1921 ttggccttta tgaggatctc tctgattttt cttgcgtcga gttttccggt aagacctttc

1981 ggtacttcgt ccacaaacac aactcctccg cgcaactttt tcgcggttgt tacttgactg

2041 gcgacgtaat ccacgatctc tttttccgtc atcgtctttc cgtgctccaa aacaacaacg

2101 geggegggaa gttcaccggc gtcatcgtcg ggaagacctg cgacacctgc gtcgaagatg 2161 ttggggtgtt ggagcaagat ggattccaat tcagcgggag ccacctgata gcctttgtac 2221 ttaatcagag acttcaggcg gtcaacgatg aagaagtgtt cgtcttcgtc ccagtaagct 2281 atgtctccag aatgtagcca tccatccttg tcaatcaagg cgttggtcgc ttccggattg 2341 tttacataac eggacataat cataggacet etcacacaca gttegeetet ttgattaacg 2401 cccagcgttt tcccggtatc cagatccaca accttcgctt caaaaaatgg aacaacttta 2461 ccgaccgcgc ccggtttatc atcccctcg ggtgtaatca gaatagctga tgtagtctca 2521 gtgagcccat atccttgcct gatacctggc agatggaacc tcttggcaac cgcttccccg 2581 acttccttag agaggggagc gccaccagaa gcaatttcgt gtaaattaga taaatcgtat 2641 ttgtcaatca gagtgctttt ggcgaagaag gagaataggg ttggcaccag cagcgcactt 2701 tgaatcttgt aatcctgaag gctcctcaga aacagctctt cttcaaatct atacattaag 2761 acgactcgaa atccacatat caaatatccg agtgtagtaa acattccaaa accgtgatgg 2821 aatggaacaa cacttaaaat cgcagtatcc ggaatgattt gattgccaaa aataggatct 2881 ctggcatgcg agaatctcac gcaggcagtt ctatgaggca gagcgacacc tttaggcaga 2941 ccagtagate cagaggagtt catgateagt gcaattgtet tgteectate gaaggaetet 3001 ggcacaaaat cgtattcatt aaaaccggga ggtagatgag atgtgacgaa cgtgtacatc 3061 gactgaaatc cctggtaatc cgttttagaa tccatgataa taattttttg gatgattggg 3121 agctttttt gcacgttcaa aattttttgc aacccctttt tggaaacgaa caccacggta 3181 ggctgcgaaa tgcccatact gttgagcaat tcacgttcat tataaatgtc gttcgcgggc 3241 gcaactgcaa ctccgataaa taacgcgccc aacaccggca taaagaattg aagagagttt 3301 teactgeata egacgattet gtgatttgta tteageceat ategttteat agettetgee 3361 aaccgaacgg acatttcgaa gtactcagcg taagtgatgt ccacctcgat atgtgcatct 3421 gtaaaagcaa ttgttccagg aaccagggcg tatctcttca tagccttatg cagttgctct 3481 ccagcggttc catcttccag cggatagaat ggcgccgggc ctttctttat gtttttggcg 3541 tcttccattt tgtatgtgtg tgtatgtata tatgctttgg gaatatatgt tcacagaatg 3601 acgactttga aagcgcgttt gaatttttaa acgaaaatct ccgtgtggct gatatttttt 3661 getttttget tttttteaac cacceggatt tttgettttt tteaaaacaa eecacegace 3721 gtaaatgttg tgtgttcttg tttctgttgg ggctgctttc ttttagagga gggaggatgc 3781 attcagagtt aatatattat atgtgctccg agatctgtgg tatacgagga gttggtgttg 3841 gctttttaat acacaaaata cgcctaaacg cgaggaggcg tcgttgaaac gtaaaggtac 3901 tattactcgt gctacctcgt agatattacc gcgcgtaatt agaagtcgtg ggagttgttg 3961 tegttgtegt tgtattegte gaegegegtt tggtegtgta eaetttaete gegeegegae //

Validation of B. prasinos transgenic lines by PCR (Figure 4b)

Oligonucleotide sequences

1(Fas5) CCACCTTTACCTCTGCCGGACATTGTGA

2TCGAAACGTGAGTCTTTTCCTTACC

3 TCGCCTCGACATCATCTGCCCAGATGC

4 (Ras5)GCGCGGTAATATCTACGAGGTAGCACGAG

Reference:

Nikolai A. Shevchuk, Anton V. Bryksin, Yevgeniya A. Nusinovich, Felipe C. Cabello, Margaret Sutherland, Stephan Ladisch / Nucleic Acids Research/, Volume 32, Issue 2, 16 January 2004, Page e19, https://doi.org/10.1093/nar/gnh014

Cell preparation

- 1 1) Starting from a culture of Ostreococcus tauri, RCC809 or Bathycoccus in stationary phase, innoculate cultures at 1 million cells/ml as determined by flow cytometry (Accuri C6 BD) in 200 ml plastic flasks in Artificial Seawater supplemented with Keller medium supplement (trace metals, vitamins, nitrate and Phosphate as described in Djouani Tahri et al., PloS ONE 2011). For each transformation (including control), you should plan on using 50 ml de culture in exponential phase.
 - 2) Grow cells for 4 to 5 days depending on the light conditions, until they reach densities of 30 to 40.10⁶ cells/ml.
 - 3) Count cells by flow cytometry. Check by SyBR Green II straining that bacterial contamination is below 2%.
 - 4) Transfer lcultures to 50 ml Falcon tubes.
 - 5) Centrifuge at 8000g for 10 min at 4°C.
 - 6) Remove the supenatant, resuspend the cell pellet in 1 ml de sorbitol 1M (pH 7.5) in H₂O MQ, at 4°C.
 - 7) Transférer the cell suspension to 1.5 ml eppendorf.
 - 8) Centrifuge at 8000g for 10 min at 4°C.
 - 9) Remove 900 µl of supernatant
 - 10) Resuspend cells by gently pipeting.

Electroporation of the transgene

- 1) Add 5μg of transgene DNA to cell suspension. Keep on ice for 5 minutes. The transgene consist of the high affinity phosphate promoter fused to the firefly luciferase (see Djouani Tahri et al., PloS one 2011).
 - 2) Transfer cells to a 2 mm electroporation cuvette (Biorad).
 - 3) Apply an electric field
 - For Ostreococcus tauri (OTTH595) : capacitance: $25\mu F$, resistance 600Ω , voltage 1.35 KV.
 - For Ostreococcus sp RCC809 : capacitance: $25\mu F$, resistance $600 \, \Omega$, voltage 1.4 KV.
 - For Bathycoccus (RC4222) : capacitance: $25\mu F$, resistance 600Ω , voltage 1.5 KV.
 - For Ostreococcus lucimarinus RCC802 : capacitance: 25μF, resistance 600 Ω, voltage 1.2KV.
 - 4) Add 1ml of fresh culture Medium to resuspend the cells.
 - 5) Add 40 ml of culture medium and transfer to a culture flask.
 - 6) Incubate at 20°C overnight in a light incubator.

At this stage, transient transgene expression can me measured or stable transformants can be selected (see relevant protocols).

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited