



Oct 18, 2019

## Making LB agar plates

Priota Islam<sup>1</sup>

<sup>1</sup>Imperial College London

1

Works for me

dx.doi.org/10.17504/protocols.io.8eqhtdw



Priota Islam  
Imperial College London



### ABSTRACT

Luria broth (LB) is a nutrient-rich media commonly used to culture bacteria in the lab. The addition of agar to LB results in the formation of a gel that bacteria can grow on, as they are unable to digest the agar but can gather nutrition from the LB within. When required the addition of an antibiotic to this gel allows for the selection of only those bacteria with resistance to that antibiotic - usually conferred by a plasmid carrying the antibiotic resistance gene.

See the link for more details: <https://www.addgene.org/protocols/pouring-lb-agar-plates/>

- 1 Purchase LB Agar from the media kitchen and autoclave it to melt

#### Content of LB Agar:

37 g pre-mixed powder consisting of:

- 5 g peptone
- 10 g peptone from casein
- 10 g sodium chloride
- 10 g agar-agar
- 1 L Sterile H<sub>2</sub>O

- 2 If the bacteria have a specific antibiotic resistance then add 1000x of that antibiotic to the melted LB agar before pouring, e.g. if we use ampicillin we add 100mg/ml to the LB in a way that the final concentration becomes 100µg/ml. So, for 200ml LB agar we add 200ul of 100mg/ul of stock ampicillin

Antibiotic	Recommended Stock Concentration	Recommended Working Concentration
Ampicillin	100 mg/mL	100 µg/mL
Bleocin	5 mg/mL	5 µg/mL
Carbenicillin	100 mg/mL	100 µg/mL
Chloramphenicol	25 mg/mL(dissolve in EtOH)	25 µg/mL
Coumermycin	25 mg/mL(dissolve in DMSO)	25 µg/mL
Gentamycin	10 mg/mL	10 µg/mL
Kanamycin	50 mg/mL	50 µg/mL
Spectinomycin	50 mg/mL	50 µg/mL
Tetracycline	10 mg/mL	10 µg/mL

Antibiotic Concentrations

- 3 Maintaining sterile techniques, pour the molten agar on large petri dishes (90mm) roughly about 30-45ml (don't need to be exact on the volume) and leave to dry overnight

- 4 Store the plates at 4C to be used later, ideally it should be used within 2 weeks but can be used for longer if not contaminated or dehydrated but never more than 3months



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited