

# Yellow fever virus real-time RT-PCR

Judy Northill, Alyssa Pyke, Ian Mackay

## Abstract

A real-time RT-PCR targeting the 5' untranslated region of Yellow fever virus. This protocol was designed and developed at this laboratory.

**Citation:** Judy Northill, Alyssa Pyke, Ian Mackay Yellow fever virus real-time RT-PCR. **protocols.io**

dx.doi.org/10.17504/protocols.io.rszd6f6

**Published:** 19 Jul 2018

## Guidelines

- If using a different brand or model of real-time thermocycler, check the concentration of ROX is adequate.
- Method assumes the user is familiar with the thermocycler and software used to run the protocol.

## Materials

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit 11732088 by Life Technologies

## Protocol

### Oligonucleotide sequences

#### Step 1.

Name	5'-3'
YF-TaqFor	TGTGCTAATTGAGGTGCATTGG
YFV wildrev	TCTCTGCTAATCGCTCAACGAA
YFV-Prob	6FAM-AATCGAGTTGCTAGGCAATAAACACATTTGGA-BHQ1

### Reaction set-up

#### Step 2.

- Assay has been used on both a Rotor-Gene 6000 / Rotor-Gene Q 5-plex using 100-place rotor discs.
- Total reaction volume is 20µL.
- Prepare sufficient for number of reaction plus a 'dead volume' usually 2 extra. Adjust as necessary if using a robotic dispenser.

Reagent	Vol (μL) x1	Final reaction concentration
Nuclease free water	4.41	
YF-TaqFor (200pmol/μl)	0.03	300nM
YFV wildrev (200pmol/μl)	0.09	900nM
YFV-Prob (100pmol/μl)	0.03	150nM
2 X Reaction mix <sup>1</sup>	10	1X
SuperScript® III/Platinum® Taq Mix <sup>1</sup>	0.4	
ROX Reference Dye (25μM) <sup>1</sup>	0.04	
<b>TOTAL</b>	<b>15</b>	

<sup>1</sup> Superscript III Platinum One-step qRT-PCR kit, Cat no. 11732088

Dispense 15μL to each reaction well.

Add 5μL of template, extracted RNA, controls or NTC (nuclease-free water).

Total reaction volume is 20μL

## Amplification

### Step 3.

The assay has been optimised and validated for the Rotor-Gene 6000 and Rotor-Gene Q thermocyclers.

RT-PCR

50°C	5min	
95°C	2min	
95°C	3s	40X
60°C	30s*	


\*Florescence acquisition step

## Result Analysis

### Step 4.

The definition used for a satisfactory positive result from a real-time fluorogenic PCR should include each of the following:

1. A **sigmoidal curve** – the trace travels horizontally, curves upward, continues in an exponential rise and followed by a curve towards a horizontal plateau phase
2. A **suitable level of fluorescence** intensity as measured in comparison to a positive control (y-axis)
3. A defined threshold ( $C_T$ ) value which the fluorescent curve has clearly exceeded (Fig.1 arrow) and which sits early in the log-linear phase and is <40 cycles
4. A flat or non-sigmoidal curve or a curve that crosses the threshold with a  $C_T$  value >40 cycles is considered a negative result
5. NTCs should not produce a curve



**Figure 1.** Examples of satisfactory sigmoidal amplification curve shape when considering an assay's fluorescent signal output. The crossing point or threshold cycle ( $C_T$ ) is indicated (yellow arrow); it is the value at which fluorescence levels surpass a predefined (usually set during validation, or arbitrary) threshold level as shown in this normalized linear scale depiction. LP-log-linear phase of signal generated during the exponential part of the PCR amplification; TP-a slowing of the amplification and accompanying fluorescence signal marks the transition phase; PP-the plateau phase is reached when there is little or no increase in fluorescent signal despite continued cycling.