

Transient luciferase expression in *Ostreococcus* (OTTH595, RCC809) and *Bathycoccus*

Francois-Yves Bouget, Valérie Vergé et Jean-Claude Lozano

Abstract

This protocol describes the quantification of luciferase activity in *Ostreococcus* and *Bathycoccus* cells that have been transformed with a firefly luciferase coding sequence upstream of the OTTH595 high affinity transporter promoter. Using this methods we were able to show that *O. tauri* promoters work in *Ostreococcus* RCC809 and *Bathycoccus*.

Citation: Francois-Yves Bouget, Valérie Vergé et Jean-Claude Lozano Transient luciferase expression in *Ostreococcus* (OTTH595, RCC809) and *Bathycoccus*. **protocols.io**

[dx.doi.org/10.17504/protocols.io.hcib2ue](https://doi.org/10.17504/protocols.io.hcib2ue)

Published: 14 Mar 2017

Protocol

Step 1.

1) Prepare CCLR Buffer

Solution stock	Final concentration
Potassium phosphate 1 M, pH 7.8	100 mM
EDTA 500 mM, pH 8	1 mM
Triton X 100	1%
Glycérol 50%	10%

Cool the buffer on ice.

2) Start rom the 40 ml transformation (see protocol on Transient transformation). Transfer cells in a 50 ml falcon tube

3) Centrifuge at 8000g for 10 min at 4°C

4) Discard the supernatant and resuspend cells in 300 µl of CCLR in 2ml eppendorf

5) Add 2 tungsten beads.

6) Break the cells in a tissu Lyser (2 x 45 sec. at 30 Hz).

7) Centrifuge for 5 minutes à 4°C, 6000g.

- 8) Transfer the supernatant in a new Eppendorf
- 9) Centrifuge for 20 minutes at 4°C, 13000g
- 10) Save the supernatant and Keep on ice.

Protein dosage using the bicinchoninic assay (BCA)

Step 2.

Principle : Proteins reduced Cu(II) to Cu(I) under alkaline conditions. Bicinchoninic acid is highly specific for Cu (I), which forms a purple complex that can be colorimetrically quantified at a 562 nm. Absorbance is proportional to the amount of proteins.

- 1) Mix 1 ml of bicinchoninic acid with 20 µl of copper sulfate at 4% (W/v) for each sample
- 2) For each sample, add 20 µl of protein extract to the 1m BCA reagent in 1.5 ml eppendorf.
- 3) incubate tubes at 37°C for 30 min (you should include a BSA standard between 0 and 2 mg/ml).
- 4) Transfer samples in cuvettes for measurement by spectrophotometry of absorbance at 562 nm.

Luciferase assay

Step 3.

- 1) Prepare Luciferase assay reagent (LAR)

Stock solution	Final concentration
Tricine 1M pH 7.8	100mM
MgCl ₂ 500mM	5 mM
EDTA 500mM pH8	100µM
DTT 1M	3.3mM
CoA sodium salt hydrate MM 767.53	270µM 2.07mg/10ml

D Luciférin 10mM	500μM
ATP dissodium salt hydrate MM 551 275.5 mg/ml	500μM

- 2) Keep LAR on ice in the dark until use
- 3) Transfer cell extracts in 96 wells white microplates (Nunc)
- 4) Place the microplate in a luminometer (such as Berthold Centro)
- 5) Read the luminescence 5 seconds after automated injection of LAR

Quantification of luciferase expression

Step 4.

Normalize the luminescence value (RLU) by the amount of protein. You should include a background control (no cell extract) .