

# The influence of agricultural tillage practices on soil biodiversity: Soil metagenomic methods, arthropod community

Jeff Strohm, Robert Hanner, Richard J Heck

## Abstract

This protocol provides the sampling and molecular biology lab methods used to prepare arthropod amplicons for MiSeq sequencing.

**Citation:** Jeff Strohm, Robert Hanner, Richard J Heck The influence of agricultural tillage practices on soil biodiversity: Soil metagenomic methods, arthropod community. **protocols.io**

[dx.doi.org/10.17504/protocols.io.efkbbkw](https://dx.doi.org/10.17504/protocols.io.efkbbkw)

**Published:** 28 Mar 2018

## Protocol

### Sample collection

#### Step 1.

Bulk soil samples were collected to a depth of 15cm using a trowel and placed in 2L plastic bags.

### Soil processing

#### Step 2.

Samples were placed in BioQuip Collapsible Berlese funnels with the only modification being an addition of a two-ply layer of cheese cloth to the metal grate.



### REAGENTS

Collapsible Berlese Funnel [2832](#) by [BioQuip](#)

### Soil processing

#### Step 3.

60W light bulbs were used to dry the soil samples until they felt completely desiccated (1 week).

### Soil processing

#### Step 4.

Emerging arthropods were collected in 70% ethanol.

### Soil processing

#### Step 5.

After fixing, 100% ethanol was added to increase the concentration to 95% and the samples were stored at -20 °C. Since only eight Berlese funnels were available, half of the bulk soil samples were stored at 4 °C for 1 week.

### Specimen sorting

#### Step 6.

Arthropods were visually sorted under a dissecting microscope and identified to the taxonomic levels of class or order

### Specimen sorting

### Step 7.

Larger arthropods such as Coleoptera or Diplopoda were sub-sampled by removing either a leg or a few body segments. For smaller Arthropods such as mites and Collembola, entire organisms were used.

### Specimen sorting

### Step 8.

All organisms from a sample were placed together in a 1.5mL Eppendorf tube in 95% ethanol and stored at -20 °C until DNA extraction.

### Arthropod DNA Extraction

### Step 9.

Ethanol was removed from all tubes using a heated SpeedVac centrifugal evaporator.

### Arthropod DNA Extraction

### Step 10.

The glass fiber insect DNA extraction protocol developed for the Canadian Center for DNA Barcoding was used (supplementary methods)

### PCR inhibitor removal

### Step 11.

The MO BIO PowerMax Soil DNA isolation kit was chosen as it contains two inhibitor precipitation reagents and was used for the microbial soil samples from the study



### REAGENTS

PowerMax® Soil DNA Isolation Kit [12988-10](#) by [Mobio](#)

### PCR inhibitor removal

### Step 12.

All reagent volumes were scaled down by a factor of 0.03 so the entire volume of DNA extract could be purified using Epoch Biolab spin columns and all centrifuge spin times and speeds were adopted from the bulk Arthropod DNA extraction protocols previously outlined

### PCR inhibitor removal

### Step 13.

All samples were eluted in 60µL of solution C6 (10mM Tris-HCL).

### PCR

### Step 14.

The high fidelity polymerase KAPA HiFi HotStart in the ReadyMix format was used.



### REAGENTS

HotStart ReadyMix (KAPA HiFi PCR kit) [KK2601](#) by [Kapa Biosystems](#)

### PCR

### Step 15.

Annealing temperature gradient experiments were used to select 52 °C, 55 °C and 58 °C as "low, medium and high" annealing temperatures

### PCR

### Step 16.

Reaction set-up and cycling conditions were identical to those used for microbial PCR with aside from the annealing temperatures and with 0.75µL primer cocktail (0.6µM) being used. As with microbial PCR, a negative control was set-up to monitor contamination.

### PCR Purification

### Step 17.

Aline PCRClean DX was used to purify all amplicons, a ratio of 0.55 (beads):1 (sample) was used.



## REAGENTS

PCRClean DX [C-1003-5](#) by [Aline Biosciences](#)

### PCR Purification

#### Step 18.

An Invitrogen E-gel was used to confirm purification success.

### PCR Purification

#### Step 19.

Amplicons were stored at 4 °C for index PCR the following day.

### Index PCR

#### Step 20.

A master mix tube was prepared for each index primer. Half the total volume of master mix was added to each reaction well from either master mix tube containing the forward or reverse index primer

### Index PCR

#### Step 21.

For PCR setup and cycling conditions see supplementary methods

### Purification

#### Step 22.

Aline PCRClean DX was used for purification, but this time, at the full concentration of 1.8X (beads): 1(sample).



## REAGENTS

PCRClean DX [C-1003-5](#) by [Aline Biosciences](#)

### Purification

#### Step 23.

Amplicons were stored at 4 °C to be quantified and normalized the following day.

### Quantification and Normalization

#### Step 24.

An Invitrogen Q-Bit Fluorometer was used to measure all amplicon concentrations

### Quantification and Normalization

#### Step 25.

All samples were normalized to 4nmol with the addition of 10mM Tris-HCL (pH 7.6).

### Quantification and Normalization

#### Step 26.

Amplicons were then pooled and submitted to the Advanced Analysis Center Genomics Facility at The University of Guelph for Illumina MiSeq sequencing using the 250bp paired end read chemistry.

### Sequencing

#### Step 27.

Amplicons were loaded into the MiSeq as a 6pM denatured library which was 15% PhiX.