



Jun 12,
2019

K and K+Si medium

Roscoff Culture Collection¹

¹CNRS-Sorbonne Université, Station Biologique, Place G. Tessier 29680 Roscoff FRANCE

Working

dx.doi.org/10.17504/protocols.io.s34egqw

Roscoff Culture Collection



Daniel Vaultot

Station Biologique, Roscoff, France



ABSTRACT

Medium to grow most phytoplankton species in small eukaryotes (Mamiellophyceae, Pelagophyceae etc...).

When using for diatoms and coccolithophorids add silica (Si).

Reference

Keller, M.D., Selvin, R.C., Claus, W. & Guillard, R.R.L. 1987. Media for the culture of oceanic ultraphytoplankton. J. Phycol. 23:633–8.

- 1 Filter 1L of old seawater of at least two months on prefilter and filter 0,2 microns
- 2 Heat seawater during 20min at 100°C
- 3 Under hood, to seawater, add these nutriments that have been autoclaved (excepted vitamin) :

Quantity	Compound	Stock Solution
1.0 mL	Sodium Nitrate (NaNO ₃)	75.0 g/L of H ₂ O
1.0 mL	Ammonium Chloride (NH ₄ Cl)	2.68 g/L of H ₂ O
1.0 mL	Na ₂ -Glycérophosphate (C ₃ H ₇ O ₆ PNa ₂ , 5 à 6 H ₂ O)	3.06 g/L of H ₂ O
1.0 mL	TRIS-Base (pH7.2)	121.1g/L of H ₂ O
1.0 mL	K Trace Metal Solution	(see recipe below)
0.1 mL	f/2 Vitamin Solution	(see recipe below)
For K+Si	Add	
1.0 mL	Sodium Metasilicate Nonahydrate (Na ₂ SiO ₃ .9H ₂ O)	30 g/L of H ₂ O

- 4 Filter the medium on 0.2 µm.

5 K Trace Metal Solution

- To 900 mL of H₂O add :

Quantity	Compound	Stock Solution
41,6g	EDTA Disodium Salt Dihydrate ($\text{Na}_2\text{EDTA} \cdot 2\text{H}_2\text{O}$)	
3,15g	Hexahydrated ferric chloride ($\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$)	1,5 cc (liquide)
1.0 ml	Sodium Molybdate Dihydrate ($\text{Na}_2\text{MoO}_4 \cdot 2\text{H}_2\text{O}$)	6,3 g/L of H_2O
1.0 ml	Zinc Sulfate Heptahydrate ($\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$)	22.0g/L of H_2O
1.0 ml	Cobalt Chloride Hexahydrate ($\text{CoCl}_2 \cdot 6\text{H}_2\text{O}$)	10.0g/L of H_2O
1.0 ml	Manganese (II) chloride, tetrahydrate ($\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$)	180.0g/L of H_2O
1.0 ml	Copper(II) sulfate pentahydrate ($\text{Cu SO}_4 \cdot 5\text{H}_2\text{O}$)	4.9g/L of H_2O
1.0 ml	Selenous acid (H_2SeO_3)	1,29g/L of H_2O

- Make final volume up to 1.0L using H_2O . Heat to dissolve

6 f/2 Vitamin solution

- In 100mL of distilled water, dissolve 0,05g of biotin (vit. H) and 0,05g of cyanocobalamin (vit. B_{12}) = solution A
- In 0,5mL of solution A, add 0,05g of thiamine HCl (vit. B_1) and complete final volume to 50mL of distilled water
- Filter sterilize on Millipore 0,22 μm
- Store in refrigerator or freezer



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited