

0.5M EDTA Version 2

Sarah Hessen-Schmidt

Abstract

Citation: Sarah Hessen-Schmidt 0.5M EDTA. protocols.io

dx.doi.org/10.17504/protocols.io.fjrbkm6

Published: 10 Aug 2016

Protocol

Step 1.

46.525g EDTA disodium 2H₂O in 200ml MilliQ



REAGENTS

Ethylenediaminetetraacetic acid disodium salt dihydrate E5134-1KG by [Sigma Aldrich](#)

Step 2.

Adjust to pH 8.0 with NaOH pellets (5g)

Step 3.

Autoclave