

# Chlamydomonas reinhardtii nuclear transformation by electroporation.

João Vitor Molino, Beth Rasala, Stephen Mayfield

## Abstract

This protocols describe the steps required for nuclear transformation of *Chlamydomonas reinhardtii* by electroporation.

**Citation:** João Vitor Molino, Beth Rasala, Stephen Mayfield Chlamydomonas reinhardtii nuclear transformation by electroporation.. **protocols.io**

dx.doi.org/10.17504/protocols.io.kfkctkw

**Published:** 24 Oct 2017

## Guidelines

Cell density for harvesting is important to overall transformant yields. It should be aimed to mid-log phase cells.

## Before start

- Separate cuvettes
- Prepare a ice bath
- Allow linearized vectors to melt

## Protocol

### DNA preparation

#### Step 1.

1. Digest 3-5 µg of DNA with appropriate enzymes, for 3-5h at 37°C.
2. Column purify digestion (Avoid gel purify, since vector backbone helps to prevent intracellular DNAses action)
3. Quantitate by absorbance measuring (Nanodrop)

📄 **AMOUNT**

5 µg Additional info: uncut vector

🌡 **TEMPERATURE**

37 °C Additional info:

### Cells preparation

## Step 2.

1. Aseptically inoculate 250 mL of TAP media with wild type cells. Either by scraping cells of a plate with a inoculating loop or from a previous cultured cells.
2. Incubate at 25°C, under constant shaking (150-180 RPM) and light (60  $\mu\text{mol}$ s de photons/m<sup>2</sup>s) until a cell density of  $3\text{-}6 \times 10^6$  cells/mL is reached
3. Pellet cells in centrifuge tubes. Separate culture in sufficient amount of 50mL tubes (Falcons), and centrifuge for 10 min at 2000xg
4. Gently resuspend cells at  **$3\text{-}6\text{-}10^8$  cells/mL** in **TAP/40mM sucrose**. (250 mL of culture usually yields 10 transformations)

### AMOUNT

250 ml Additional info: TAP media

### TEMPERATURE

25 °C Additional info: wild type cells culture

### NOTES

**João Vitor Molino** 24 Oct 2017

Usually cells will reach desired concentration in 3 days, depending on the amount of inoculum

## Transformation

### Step 3.

1. Add 250  $\mu\text{L}$  of resuspended cells (at  $3\text{-}6\text{-}10^8$  cells/mL) to each cuvette. Pipette up and down on DNA sample. Flick cuvette to mix DNA and cells. Shake cells to the bottom of the cuvette. Also add no DNA control.
2. Incubate cells with DNA **ON ICE** for 5-10 min
3. Wipe cuvette (to remove condensed water) and electroporate.
  - Voltage = 800 V
  - Capacitance = 25  $\mu\text{F}$
  - Resistance = infinity ( $\infty$ )
  - Gap = 4 mm
5. Let recover for 5-10 min on cuvette
6. Add 10 mL of Tap/40mM sucrose into one clean 50mL falcon tube per electroporation
7. Gently transfer cells from cuvette to TAP/40 mM sucrose. Rinse cuvette with TAP/40 mM sucrose and transfer remaining cells to Falcon tube
8. Incubate at room temperature on rocker overnight
9. Pellet cells by centrifuging for 10 min at 2000xg
10. Aseptically pour off supernatant. Add 600  $\mu\text{L}$  of TAP/40 mM sucrose to pellet. Gently resuspend cells and pipette onto 2 plates with appropriate antibiotics. ie. 400  $\mu\text{L}$  of cells per plate, and let it dry aseptically without plate cover.
11. Parafilm and place plates under constant light (60  $\mu\text{mol}$ s de photons/m<sup>2</sup>s), 25 °C. Colonies should be visible in 5-7 days.

### AMOUNT

600  $\mu\text{L}$  Additional info: TAP 40 mM sucrose to re-suspend cells

## Warnings

- High voltage is used in the electroporation, use EPIs and avoid contact with electrodes on the electroporator.