



Apr 02, 2019

Working

Transfection of construct containing kinetoplastid *Blastocrithidia* sp. p57 UTR's in *Diplonema papillatum* using AMAXA Nucleofactor apparatus.

Binnypreet Kaur¹

¹1 Institute of Parasitology, Biology Centre, Czech Academy of Sciences, 37005, České Budějovice, Czech Republic 2 Faculty of Sciences, University of South Bohemia, 37005, České Budějovice, Czech Republic.

[dx.doi.org/10.17504/protocols.io.zp9f5r6](https://doi.org/10.17504/protocols.io.zp9f5r6)

 Binnypreet Kaur   

ABSTRACT

Transfection of construct containing kinetoplastid *Blastocrithidia* sp. p57 UTR's in *Diplonema papillatum* using AMAXA Nucleofactor apparatus.

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

SAFETY WARNINGS

BEFORE STARTING

Construct Design: A construct containing kinetoplastid *Blastocrithidia* sp. p57 containing V5-tagged aminoglycoside 3'-phosphotransferase gene (conferring resistance to neomycin) surrounded by regulatory sequences derived from the hexokinase gene was linearized with *Sma*I and named as p57-V5+G418 plasmid.

- 1 Count the cells and plan to 1 nucleofection with 5×10^7 cells for each construct.
- 2 Harvest the cells by centrifugation at 1300xg for 5 min at room temperature in Swing Bucket Rotor.
- 3 Resuspend the cell pellet in 100ul of AMAXA Human T- cell solution at 4°C (from refrigerator combine 81.8ul of Human T-cell nucleofector solution + 18.2ul Supplement).
- 4 Add 3ug of linearized as p57-V5+G418 PLASMID into the cuvette(resuspend in 10ul of H₂O).Note: Do not add anything in case of Negative control.
- 5 Put everything into the cuvette, close the cap and place in the electroporator, cuvette should only fit in one direction, but metal sides should face towards you.
- 6 Press for the Program X-001/X-014 to electroporate.
- 7 Transfer the entire cell- DNA transfectant into 10ml of Diplonema media (chloramphenicol) with the supplied disposable micropipettes.Note: Make sure labeled flasks with media should be ready before electroporation.

- 8 Immediately observe the flasks containing transfectants under a microscope and place them in the incubator at 27C + shaker and let them grow for 6-8h.
- 9 After 6-8 hours took out the flasks and make controls which should be without antibiotics - put 1ml of the electroporated cells in the first well of 24 well plates.
- 10 In the rest - 9ml -transfectants add G418 FOR the selection.
- 11 Add 1.5ml of electroporated cells in the first row (6 wells) of 24 well plates. Put 900ul of the media (+ selectable drug) in the rest of the wells. Take 100ul from first well and transfer to the next to make 10x serial dilution and so on...up to the lowest concentration 5×10^2 cells/ml.
- 12 Do the same with the negative control plating.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited