

# Fluorescent Focus Unit Assay using LICOR Imaging System

## Version 2

**Bernardo Mainou**

### Abstract

Indirect immunofluorescence infectivity assay for reovirus using LICOR Imaging System

**Citation:** Bernardo Mainou Fluorescent Focus Unit Assay using LICOR Imaging System. **protocols.io**

[dx.doi.org/10.17504/protocols.io.k8fcztn](https://dx.doi.org/10.17504/protocols.io.k8fcztn)

**Published:** 11 Dec 2017

## Guidelines

This assay is standardized to for 96-well plates using reovirus polyclonal antiserum. New batches of polyclonal antiserum or monoclonal antibodies should be tested to determine optimal dilution.

## Before start

Plates and reagents:

Corning Black with clear bottom TC treated 96 well plates (Corning #3904)

DRAQ5 #4084 from Cell Signaling

Sapphire 700 #928-40022 from LICOR

Antibody 800 from LICOR (Life Technologies ones don't work for this)

## Protocol

### Step 1.

Remove media from wells.

### Step 2.

Wash 1x with Phosphate Buffered Saline (PBS).

### Step 3.

Add 100 uL ice-cold methanol per well.

### Step 4.

Store plate at -20°C for at least 30 min (can remain at -20°C for several weeks).

**Step 5.**

Remove methanol, allow plate to come to room temperature, and excess ethanol to evaporate.

**Step 6.**

Wash wells with 150 µl PBS

**Step 7.**

Add 150 µL Dulbecco's PBS (DPBS) with 0.5% Tween-20 (DPBS-T) to each well.

**Step 8.**

Remove immediately.

**Step 9.**

Add 50 µL rabbit anti-reovirus polyclonal antiserum (1:1000 dilution) in DPBS with 1% BSA (DPBS-BSA) to experimental wells.

Add 50 µL DPBS-BSA with no antibody to background control.

**Step 10.**

Incubate for 1 h at 37°C.

**Step 11.**

Remove DPBS-BSA.

**Step 12.**

Wash 3x with DPBS-T for 5 min each wash while shaking.

**Step 13.**

Add 50 µL of DPBS-BSA to all wells.

**Step 14.**

Incubate 1 h at 37°C.

**Step 15.**

Wash plates 3x with DPBS-T for 5 min each wash.

**Step 16.**

Prepare cell staining solution in DPBS-BSA:

Secondary antibody (e.g. Goat Anti Rabbit LICOR IRDye 800CW) - 1:1000

Draq5 - 1:10,000

Sapphire700 - 1:1000

Note: we have found that we get less background with the LICOR IRDye 800CW secondary antibody than with Goat Anti Rabbit Alexa 790)

**Step 17.**

Add 50 uL cell staining solution to all wells that were treated with primary.

Add 50 uL of DPBS-BSA with secondary antibody (1:1000) **only** to background control wells.

**Step 18.**

Remove solution.

**Step 19.**

Incubate 1 h at 37°C.

**Step 20.**

Wash 3x with 150 uL DPBS-T.

**Step 21.**

Add 50 uL of water.

**Step 22.**

Scan plates on LICOR Odyssey Imaging System.

a. Focus offset = 3.0 (depends on the plates)

b. 700nm Intensity = 6.5

c. 800nm Intensity = 7.5

**Step 23.****Plates and reagents:**

Corning Black with clear bottom TC treated 96 well plates (Corning #3904)

DRAQ5 #4084 from Cell Signaling

Sapphire 700 #928-40022 from LICORAntibody 800 from LICOR (Life Technologies secondary antibodies give extra background and should be avoided)