



Oct 01, 2019

Single worm lysis

[Ida Barlow](#)¹¹Imperial College London

1

Works for me

[dx.doi.org/10.17504/protocols.io.587g9zn](https://doi.org/10.17504/protocols.io.587g9zn)

Ida Barlow



ABSTRACT

Extraction of genomic DNA (gDNA) from single adult *C. elegans* suitable for downstream molecular biology applications such as PCR.

MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

[Proteinase K](#)

MATERIALS TEXT

Proteinase K from *Tritirachium album*

lyophilized powder, ≥ 30 units/mg protein

Made up to 20mg/ml in Tris-HCl

- 1 Prepare 1mL 1X PCR buffer with 0.4mg/ml Proteinase K:

▢ 855 µl water

▢ 95 µl 10X PCR buffer (pH8.3)

▢ 50 µl 20mg/ml Proteinase K solution

▢ 1000 µl Total volume

nb. volumes can be scaled as required



10X PCR buffer (pH8.3) preparation:

▢ 18.64 g Potassium chloride (KCl)

▢ 3.7 g Trizma base

▢ 3.07 g Tris-HCl

▢ 0.71 g Magnesium Chloride (MgCl₂)

▢ 500 ml water

Autoclave

Cool (temperature affects pH)

pH meter and add hydrochloric acid to achieve the correct pH

- 2 Prepare worm lysis tubes with PCR buffer and allow ▢ 15 µl 1X PCR buffer + PK per worm.

Pipette buffer into lid of PCR tube (up to 45ul can fit into tube lid)

- 3 Pick a single worm at a time into the liquid in the 200uL PCR tube. Make sure you can see the worm swimming in the buffer
- 4 Close the lids onto the tubes of the 200uL tube and spin in a microfuge so that the liquid and worms are at the bottom of the tube
- 5 Burst the cells open by freezing the tubes in the -80°C freezer for 10 minutes
- 6 Transfer tubes to thermocycler and incubate on worm lysis program:
60°C - 1hr 10 mins
95°C - 15 mins
- 7 Store gDNA in -20°C freezer until required. Use
2ul per 25uL PCR reaction.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited