

# Human T Cell Activation with anti-CD3 Antibodies (clone UCHT1, OKT3 or HIT3a) Version 2

BioLegend, Inc.

## Abstract

**Citation:** BioLegend, Inc. Human T Cell Activation with anti-CD3 Antibodies (clone UCHT1, OKT3 or HIT3a). [protocols.io](https://doi.org/10.17504/protocols.io.hv7b69n)  
[dx.doi.org/10.17504/protocols.io.hv7b69n](https://doi.org/10.17504/protocols.io.hv7b69n)

**Published:** 09 May 2017

## Guidelines

### Materials:

- Sterile PBS
- Anti-human CD3 Antibody
  - Clone UCHT1 (LEAF™ format, Cat. No. 300413/300414/300432; Ultra-LEAF™ format, Cat. No. 300437/300438)
  - Clone OKT3 (LEAF™ format, Cat. No. 317303/317304/317315; Ultra-LEAF™ format, Cat. No. 317325/317326)
  - Clone HIT3a (LEAF™ format, Cat. No. 300313/300314; Ultra-LEAF™ format, Cat. No. 300331/300332)
- Cell culture medium (e.g., RPMI-1640 or IMDM supplemented with 10% FBS and 2mM L-glutamine)
- Sterile single-cell suspension of Ficoll-Hypaque-purified peripheral blood mononuclear cells, isolated T cells, or T cell subsets
- 96-well flat-bottom tissue culture plates with lids (e.g., Costar® Cat. No. 3596)

\* Soluble forms of LEAF™ purified UCHT1 (1 µg/ml) or LEAF™ purified HIT3a (0.01 – 0.1 µg/ml) may be used to activate T cells from PBMC cell populations.

## Protocol

### Step 1.

Prepare a 10 µg/ml solution of anti-CD3 (clone UCHT1, OKT3, or HIT3a) in sterile PBS.

### Step 2.

Dispense 50 µl of the antibody solution to each microwell of the 96-well assay plate. For the unstimulated control wells, add 50 µl of sterile PBS.



### REAGENTS

96-well flat-bottom tissue culture plates with lids [3596](#) by [Corning](#)

### Step 3.

Seal plate. Incubate at 37°C for 2 hours or 4°C overnight.



### DURATION

02:00:00

**Step 4.**

Aseptically decant antibody solution from the microwell plate.

**Step 5.**

Wash plate microwells 3 times with sterile PBS (wash 1/3).

**Step 6.**

Wash plate microwells 3 times with sterile PBS (wash 2/3).

**Step 7.**

Wash plate microwells 3 times with sterile PBS (wash 3/3). Discard liquid.

**Step 8.**

Prepare single cell suspension of cells of interest in supplemented cell culture medium to  $1-2 \times 10^6$ /ml.

**Step 9.**

Aliquot 200  $\mu$ l cell suspension into plate microwells. Cover with lid. Incubate at 37°C in 5% CO<sub>2</sub> and 100% humidity for 3 days.