



Jun 12,
2019

CGAP Human Oesophagus Epithelium Dissociation - Tissue Stability

Forked from [CGAP Human Oesophagus Epithelium Dissociation](#)

[Anna Wilbrey-Clark](#)¹, [Adam Hunter](#)²

¹Wellcome Trust Sanger Institute, HCA, ²CGAP

Working [dx.doi.org/10.17504/protocols.io.34fgqtn](https://doi.org/10.17504/protocols.io.34fgqtn)

[Human Cell Atlas Method Development Community](#)

[Anna Wilbrey-Clark](#) ⚡

1	Material	Quantity	Supplier Info
	100mm Petri Dish	1	Corning (430591)
	Cold PBS	50ml	GIBCO (14190-144)
	Forceps	2	ThermoFisher UK Ltd (15232290)
	Scalpels	2	Swann-Morton Ltd (0507)
	50ml Falcon Tubes	2	Falcon (352098)
	15ml Falcon Tubes	1	Falcon (352097)
	Dissociation Agent (Trypsin-EDTA 0.25%)	40ml	GIBCO (25200-056)
	DNase I (Stock solution 10mg/ml)	400ul	Sigma (11284932001)
	RPMI + 20% FBS	25ml	Gibco (42401042) + Sigma (F7524-50ML)
	70um Cell Strainer	2	ThermoFisher UK Ltd (15370801)
	2.0ml Syringe	1	ThermoFisher UK Ltd. (10673555)
	1.5ml DNA LoBind Eppendorf Tubes	1	Eppendorf (0030 108.051)
	BSA	267µl 7.5% BSA in 50ml PBS-/-	Sigma (A8412)
	0.5ml DNA LoBind Eppendorf Tubes	2	Eppendorf (0030 108.035)
	Red Cell Lysis Buffer	5ml	Life Technologies Ltd. (00-4333-57)
	Trypan Blue	40ul	Fisher Scientific (11414815)
	C-Chips	2	Cambridge Bioscience (DHC-N01-50)

- 2 Receive oesophagus sample from mid-region in hypothermasol FRS solution (Sigma H4416).
- 3 Wash the samples with 10ml cold PBS to remove any residual contamination, stomach content and loose mucus.
- 4 Pour oesophagus onto 100mm glass petri dish and add another 10ml fresh cold PBS.
- 5 Open the samples longitudinally.
 - Epithelium/mucosa should be a relatively loose, yellowish layer on the lumen side.

- 6 Using two forceps (scissors if necessary) separate the mucosa/epithelium layer from the stroma. Place them onto separate 100mm petri dishes each with 10ml PBS (this is a relatively crude separation protocol and is not 'pure' but avoids long overnight incubations in dispase).
 - From this stage on only process the mucosa/epithelium part for dissociation.
- 7 Aspirate PBS.
- 8 In a few drops of PBS (~200µl), finely mince the epithelium using two scalpels simultaneously.
 - Too much PBS in the dish will make it more difficult to mince.
- 9 Add 19ml Trypsin-EDTA 0.25% containing 100µg/ml DNase I to the dish and transfer the tissue to a 50ml falcon tube. (DNase is critical for this protocol). Note later tests indicated that this volume can be reduced to 5ml total.
- 10 Wash the scalpel and the dish with a further 1ml of Trypsin-EDTA 0.25% and transfer it into the Falcon with the tissue.
- 11 Incubate the biopsies for 30 min at 37°C on a rocker.
- 12 Centrifuge sample 500g, 5 min at 4°C.
- 13 Add 20ml fresh 0.25% Trypsin-EDTA with 100µg/ml DNase and incubate at 37°C on rocker for a further 15min. (Note it was later determined that this second incubation step is not necessary for oesophagus, only for stomach, which causes a change in pH of the dissociation agent).
- 14 Pass through a 70µM cell strainer into a 50ml falcon tube.
- 15 Add 20ml of RPMI + 20% FBS to inactivate trypsin, mash the remaining tissue in the strainer with a syringe plunger and wash again with 5ml RPMI + 20% FBS. There is usually some undigested material remaining).
- 16 Centrifuge at 500g for 5 minutes, 4°C.
- 17 Wash cells with 5ml cold PBS-/-.
- 18 Centrifuge at 500g for 5 minutes, 4°C.
- 19 Add 5ml Red Cell Lysis buffer to the pellet and resuspend.
- 20 Place on ice 5 min with periodic agitation.
- 21 Add 10ml of Cold PBS-/-.

- 22 Filter through a 70µm cell strainer. (Often, small white filaments of undissociated material are still visible at this stage which must be removed before 10x loading).
- 23 Centrifuge at 500g for 5 minutes, 4°C. Remove supernatant.
- 24 Resuspend in 300-1000µl cold PBS with 0.04% BSA.
- 25 Count cells and viability manually using C-chip / trypan blue.
- 26 Proceed to 10X Preparation for scRNA sequencing.
- 27 Ensure all unused tissue, equipment and tubes that have been in contact with primary tissue are placed into Virkon in sweetie jar for a minimum of 1 hour. After this time aspirate and disposing in relevant sharps or waste routes.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited