

# Dispersing Labyrinthulomycete Cell Aggregates by Sonication

Mariana Rius, Kashyapa Bandaralage

## Abstract

Sonation procedure implemented to reduce labyrinthulomycete cell aggregation from culture for subsequent growth analyses.

**Citation:** Mariana Rius, Kashyapa Bandaralage Dispersing Labyrinthulomycete Cell Aggregates by Sonication. protocols.io

dx.doi.org/10.17504/protocols.io.hg2b3ye

**Published:** 31 Mar 2017

## Protocol

### Set up QSONICA Q800R Sonicator

#### Step 1.

Turn on QSONICA Q800R Sonicator. Add required distilled water and allow the water to reach a temperature of 4 C.

#### 📌 NOTES

**Mariana Rius** 30 Mar 2017

Output frequency: 20KHz

### Prepare cells

#### Step 2.

Pipette 1 - 1.2ml of cell culture into a 1.5ml eppendorf tube. Gently vortex. Invert the tube several times before placing in QSONICA adapter to maintain the cells in suspension.

### Sonicate cells

#### Step 3.

Program QSONICA to sonicate for 1 minute at 100% amplitude followed by a 30 second intermission and another 1 minute at 100% amplitude. Put on necessary protective ear wear. Place the 1.5ml tubes of suspended cells into QSonica adapter and then into QSONICA. Close the lid of the machine and immediately begin sonication procedure by hitting START.

#### 📌 NOTES

**Mariana Rius** 30 Mar 2017

*Aurantiochytrium limacinum* (ATCC MYA-1381)

Mass cell lysis was seen after 2 minutes of sonication at 100% amplitude

*Schizochytrium aggregatum* (ATCC 28209)

Can withstand up to 3 minutes sonication at 100% amplitude without significant cell lysis

*Aplanochytrium stocchinoi*

Can withstand 4 minutes at 100% amplitude with little to no visible cell lysis

#### Microscopy of post-sonicated cells

##### **Step 4.**

Identify cell mortality and remaining cell aggregation following sonication using microscopy and adjust sonication procedure as necessary.

This protocol has been successfully used to break up aggregates of *Aplanochytrium*, *Aurantiochytrium*, and *Schizochytrium* from cultures grown in liquid 790 medium.