

Saltwater BG-11 recipe Version 2

Jeric Harper

Abstract

Citation: Jeric Harper Saltwater BG-11 recipe. protocols.io

dx.doi.org/10.17504/protocols.io.fz6bp9e

Published: 05 Oct 2016

Protocol

Step 1.

Stock solutions:

Macronutrient Stock Solutions:

Prepare each stock solution separately in 32 ppt seawater.

150 g/L NaNO_3

4 g/L $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$

3.6 g/L $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$

0.6 g/L Citric Acid $\cdot 2\text{H}_2\text{O}$

0.6 g/L Ferric Ammonium Citrate

0.1 g/L $\text{Na}_2\text{EDTA} \cdot 2\text{H}_2\text{O}$

2 g/L NaCO_3

Prepare each stock solution separately in ddH₂O.

4 g/L K_2HPO_4

30 g/L NaSiO_3

Trace Metal Stock Solution:

Combine ingredients in ddH₂O.

2.86 g/L H_3BO_3

1.81 g/L $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$

0.22 g/L $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$

0.39 g/L $\text{Na}_2\text{MoO}_4 \cdot 2\text{H}_2\text{O}$

0.79 g/L $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$

0.50 g/L $\text{Co}(\text{NO}_3)_2 \cdot 6\text{H}_2\text{O}$

Vitamin Stock Solution:

Combine ingredients in ddH₂O.

0.135 g/L Vitamin B₁₂ (cyanocobalamin)

0.025 g/L Biotin

0.335 g/L Thiamine

12 g/L HEPES buffer pH 7.8

* Be sure to prepare K_2HPO_4 and NaSiO_3 stock solutions in ddH₂O to prevent precipitation.

Step 2.

Prepare media:

For each litre of BG-11:

- 750mL of 32ppt seawater
- 10 ml of K_2HPO_4 Macronutrient Stock Solution.
- 1mL of NaSiO_3 Macronutrient Stock Solution.
- 10 mL of each Macronutrient Stock Solution.
- 1ml of the Trace Metal Stock Solution.
- 1ml of the Vitamin Stock Solution.

Bring final volume to 1L.

Adjust pH of the final media to 8.2.

* Be sure to add K_2HPO_4 first, then $NaSiO_3$, then add the others stock solutions to prevent precipitation.

Step 3.

0.2 um filter sterilize into clean, autoclaved containers in a biosafety hood. Do not autoclave the media as the silica will form precipitate and be cloudy!

Step 4.

To Prepare Conjugation Plates (1/2 BG-11, 5% LB, 1% agar):

1. Dilute Saltwater BG-11 1:1 with ddH₂O.
2. Add 50mL of 1x LB for each litre of Conjugation Plate media being made.
3. Add 5g of Bacto Agar to a 1L Earlenmeyer flask for each 500mL of Conjugation Plate media being made.
4. Distribute 500 mL of Conjugation Plate media to each 1L Earlenmeyer flask.
5. Cover with flask opening with aluminum foil and autoclave at 121C for 30 minutes.
6. Allow media to cool while stirring on magnetic stir plate. Make sure the stirring is not too vigorous as to cause bubbles. It should cool enough to be able to touch, but not have the agar solidified.
7. Pour plates and allow to cool overnight.

Step 5.

To Prepare Selection Plates (1/2 BG-11, 1% agar, + antibiotics):

1. Dilute Saltwater BG-11 1:1 with ddH₂O.
2. Add 5g of Bacto Agarto to a 1L Earlenmeyer flask for each 500mL of Selection Plate media being made.
3. Distribute 500 mL of Selection Plate media to each 1L Earlenmeyer flask.
4. Cover flask opening with aluminum foil and autoclave at 121C for 30 minutes.
5. Allow media to cool while stirring on magnetic stir plate. Make sure the stirring is not too vigorous as to cause bubbles. It should cool enough to be able to touch, but not have the agar solidified.
6. Add 500ul of the required antibiotics to the Selection Plate media made. We make 500uL aliquots of 1000x antibiotic stocks. It is very important that the media is cool before adding the antibiotics or the heat will destroy them.
7. Pour plates and allow to cool overnight.