

# Cesium Chloride Protocol for Phage

Forest Rohwer

## Abstract

This protocol contains 4 sections:

1. Cesium Chloride Gradient to separate phage from bacteria. Marine viruses around 1.5g/ml
2. DNase I treatment of CsCl-purified phage
3. Formamide Extraction of DNA from phageE4.
4. CTAB

**Citation:** Forest Rohwer Cesium Chloride Protocol for Phage. **protocols.io**

[dx.doi.org/10.17504/protocols.io.c37yrm](https://dx.doi.org/10.17504/protocols.io.c37yrm)

**Published:** 25 Jan 2016

## Guidelines

Needed:

- Cesium solutions
- 0.02 µm filter
- CsCl
- Tubes
- Centrifuge @ 22,000 rpm, 4,000 rpm,
- 0.2 vol chloroform
- 10-15 µL DNase I
- RNase (if needed)
- Oak ridge tube
- 0.1 vol 2M Tris HCL (pH 8.5)
- 0.2M EDTA
- 100 µL 0.5M EDTA
- 10 Fl glycogen
- 2 vol 100% ethanol
- 70% ethanol
- Sorvall ST 8 Benchtop centrifuge
- 567 µL TE Buffer
- 30 µL 10% SDS
- 3 µL of 20 µg/mL proteinase K
- Incubator
- 100 µL 5M NaCl
- 80 µL CTAB/NaCl solution
- Microcentrifuge
- Phenol/Chloroform solution
- 0.7 vol Isopropanol
- 50 µL Sigma water

- Nanodrop

## Materials

SW41 Ti (with Beckman Ultra) Swinging Bucket [331362](#) by [Beckman Coulter](#)

Sorvall ST 8 Benchtop Centrifuge [75007200](#) by [Beckman Coulter](#)

## Protocol

Cesium Chloride Gradient to separate phage from bacteria. Marine viruses around 1.5g/ml

### Step 1.

Filter cesium solutions with 0.02µm filter

#### NOTES

**VERVE Team** 22 Jun 2015

Use filtered seawater to make solution. Mark top of first gradient layer with pen before adding next layers.

Cesium Chloride Gradient to separate phage from bacteria. Marine viruses around 1.5g/ml

### Step 2.

Mark top of first gradient layer with pen before adding next layers

Cesium Chloride Gradient to separate phage from bacteria. Marine viruses around 1.5g/ml

### Step 3.

Set up gradient as (from top to bottom):

- 1ml of 1.7 g/ml, mark top of layer
- 1ml of 1.5 g/ml, mark top of layer
- 1ml of 1.35 g/ml
- 8.5ml of 1.15 g/ml - this is phage concentrate.

Cesium Chloride Gradient to separate phage from bacteria. Marine viruses around 1.5g/ml

### Step 4.

Add CsCl to your sample.

Cesium Chloride Gradient to separate phage from bacteria. Marine viruses around 1.5g/ml

### Step 5.

Check balance of tubes.

Cesium Chloride Gradient to separate phage from bacteria. Marine viruses around 1.5g/ml

### Step 6.

Centrifuge at 22,000 rpm for 2 hours at 4°C, approximately 60-80,000xg.

#### DURATION

02:00:00

Cesium Chloride Gradient to separate phage from bacteria. Marine viruses around 1.5g/ml

### Step 7.

Pierce tube at 1.7/1.5 g/ml interface.

Cesium Chloride Gradient to separate phage from bacteria. Marine viruses around 1.5g/ml

### Step 8.

Bevel up, and collect 1.5 mls (should include 1.5 g/ml step and 1.5 to 1.35 interface)

Cesium Chloride Gradient to separate phage from bacteria. Marine viruses around 1.5g/ml

#### Step 9.

Check with slides for virus particles, etc.

Cesium Chloride Gradient to separate phage from bacteria. Marine viruses around 1.5g/ml

#### Step 10.

Use 20 µl from [viral] fraction and 500 µl from "upper" layers.

#### 📌 NOTES

**VERVE Team** 22 Jun 2015

Alternatively, can pierce at bottom of tube and take 0.5 ml fractions for first 3 mls, then make slides of each and extract from appropriate fraction.

Cesium Chloride Gradient to separate phage from bacteria. Marine viruses around 1.5g/ml

#### Step 11.

Store in fridge at 4°C until ready to extract.

DNase I treatment of CsCl-purified phage

#### Step 12.

Add 0.2 vol chloroform, mix and incubate 10 min. at room temperature

#### 📋 AMOUNT

0 µl Additional info:

#### 🕒 DURATION

00:10:00

DNase I treatment of CsCl-purified phage

#### Step 13.

Spin at 4000 rpm in table Beckman for 10 minutes to separate phases.

#### 🕒 DURATION

00:10:00

DNase I treatment of CsCl-purified phage

#### Step 14.

Save and transfer top phase to new 15 ml tubes.

DNase I treatment of CsCl-purified phage

#### Step 15.

Add 10-15 µl DNase I (1 mg/ml in H<sub>2</sub>O) to phage sample (1.2-1.5ml)

#### 📋 AMOUNT

12 µl Additional info:

DNase I treatment of CsCl-purified phage

#### Step 16.

Incubate 15 minutes at 37°C. Treat with RNase if RNA is to be extracted.

#### 🕒 DURATION

00:15:00

DNase I treatment of CsCl-purified phage

#### Step 17.

Inactivate for 15 minutes at 65°C.

#### 🕒 DURATION

00:15:00

DNase I treatment of CsCl-purified phage

#### Step 18.

Transfer all to new "oak ridge" tube for faster centrifugation later.

#### Formamide Extraction of DNA from phage

##### Step 19.

Add 0.1 volume 2 M TrisHCL (pH 8.5)/0.2 M EDTA

#### Formamide Extraction of DNA from phage

##### Step 20.

Add 100 µl 0.5 M EDTA per 10ml

 [AMOUNT](#)

100 µl Additional info:

#### Formamide Extraction of DNA from phage

##### Step 21.

Add 1 volume of formamide

#### Formamide Extraction of DNA from phage

##### Step 22.

Add 10 Fl glycogen

#### Formamide Extraction of DNA from phage

##### Step 23.

Incubate at room temperature for 30 minutes.

 [DURATION](#)

00:30:00

#### Formamide Extraction of DNA from phage

##### Step 24.

Add 2 volumes of room temperature 100% ethanol.

#### Formamide Extraction of DNA from phage

##### Step 25.

Pellet in Sorvall (12,000 rpm for 20 minutes)

 [DURATION](#)

00:20:00

#### Formamide Extraction of DNA from phage

##### Step 26.

Wash with 70% ethanol, two times.

#### Formamide Extraction of DNA from phage

##### Step 27.

Resuspend into 567 µl TE and continue with CTAB extraction.

 [AMOUNT](#)

567 µl Additional info:

 [PROTOCOL](#)

. [TE Buffer](#)

CONTACT: [VERVE Team](#)

##### Step 27.1.

Prepare 5ml of 1M Tris

 [AMOUNT](#)

61 g Additional info:

## [PROTOCOL](#)

### [. 1M Tris](#)

CONTACT: [VERVE Team](#)

#### **Step 1.1.**

60.57 g of Tris (hydroxymethyl) aminomethane in 0.5L Milli-Qwater

#### **Step 1.2.**

Bring pH to 8.0 using HCl

Making 0.5M EDTA

#### **Step 27.2.**

Prepare 1ml 0.5M EDTA

 [AMOUNT](#)

19 g Additional info:

 [PROTOCOL](#)

### [. 0.5M EDTA](#)

CONTACT: [VERVE Team](#)

#### **Step 2.1.**

Mix 18.6 g EDTA in 100ml Milli-Q water

 [AMOUNT](#)

19 g Additional info:

#### **Step 2.2.**

Bring pH to 8.0 using NaOH

 [NOTES](#)

**VERVE Team** 22 Jun 2015

EDTA will not be soluble until pH reaches 8.0–this will take time (hours)

Making 0.5M EDTA

#### **Step 27.3.**

Combine 5ml of 1M Tris and 1ml of 0.5M EDTA in 496ml dH<sub>2</sub>O

Making 0.5M EDTA

#### **Step 27.4.**

Autoclave to sterilize.

 [NOTES](#)

**VERVE Team** 22 Jun 2015

It is best to remove a working aliquot and not repeatedly access the stock

CTAB

#### **Step 28.**

Resuspend pellet in 567 µl TE.

 [AMOUNT](#)

567 µl Additional info:

CTAB

#### **Step 29.**

Add 30 µl of 10% SDS and 3 µl of 20 µg/ml proteinase K. Mix.

 [AMOUNT](#)

30 µl Additional info:

CTAB

### Step 30.

Incubate 1 hour at 37°C - 56°C

 DURATION

01:00:00

CTAB

### Step 31.

Add 100 µl of 5 M NaCl and mix thoroughly.

 AMOUNT

100 µl Additional info:

CTAB

### Step 32.

Add 80 µl CTAB/NaCl solution. Mix.

 AMOUNT

80 µl Additional info:

CTAB

### Step 33.

Incubate 10 minutes at 65°C.

 DURATION

00:10:00

CTAB

### Step 34.

Add equal volume of chloroform; mix.

CTAB

### Step 35.

Microcentrifuge for 2 minutes.

 DURATION

00:02:00

CTAB

### Step 36.

Transfer supernatant to separate tube.

 NOTES

**VERVE Team** 21 Jan 2016

Your DNA is in the supernatant (the top).

CTAB

### Step 37.

Add equal volume of phenol/chloroform to the supernatant fraction; mix.

 NOTES

**VERVE Team** 22 Jun 2015

Note: Your DNA is in the supernatant (the top)!

CTAB

### Step 38.

Microcentrifuge for 2 minutes.

 DURATION

00:02:00

CTAB

**Step 39.**

Transfer supernatant to separate tube.

CTAB

**Step 40.**

Add equal volume of chloroform to the supernatant fraction; mix.

CTAB

**Step 41.**

Microcentrifuge for 2 minutes.

 DURATION

00:02:00

CTAB

**Step 42.**

Transfer supernatant to separate tube.

CTAB

**Step 43.**

Add 0.7 volume isopropanol to the supernatant fraction.

 AMOUNT

1 µl Additional info:

CTAB

**Step 44.**

Mix gently until DNA precipitates.

CTAB

**Step 45.**

Centrifuge 15 minutes in cold.

 DURATION

00:15:00

CTAB

**Step 46.**

Wash with 70% ethanol.

CTAB

**Step 47.**

Resuspend in 50 µl Sigma water

 AMOUNT

50 µl Additional info:

CTAB

**Step 48.**

Check O.D. using the Nanodrop