



Sep 22, 2019

Glycerol Stock

[Alba Balletbó](#)¹¹Wageningen University

1

Works for medx.doi.org/10.17504/protocols.io.7j3hkqn

Alba Balletbó



MATERIALS TEXT

Overnight culture

80% Glycerol solution

- 1 Take 700 μ L of an overnight culture and add it to a sterile cryotube.
- 2 Add 300 μ L of a sterile 80% glycerol solution and invert 2-3 times.
- 3 Place the glycerol stock in -80°C until needed again.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited