

# TSS Transformation (after Chung et al., 1989, PNAS 86:2172-2175)

Lutz Berwanger

## Abstract

Transformation & Storage Solution (2X TSS) enables researchers to take advantage of the simple system described by Chung et al. <sup>1</sup> for the preparation, long-term storage and transformation of competent *E. coli*. Early log-phase cells are suspended in 1X TSS: a solution containing polyethylene glycol, dimethyl sulfoxide, and divalent cations in a bacterial growth medium.

**Citation:** Lutz Berwanger TSS Transformation (after Chung et al., 1989, PNAS 86:2172-2175). **protocols.io**

[dx.doi.org/10.17504/protocols.io.qtn dwme](https://dx.doi.org/10.17504/protocols.io.qtn dwme)

**Published:** 08 Jun 2018

## Protocol

### Preparing TSS

#### Step 1.

2 x TSS (Transformation and Storage Solution):

- 20 g PEG (3350 or 8000) 20% (w/v)
- 10 mL DMSO (after autoclave) 10% (v/v)
- 2,03 g MgCl<sub>2</sub> \* 6H<sub>2</sub>O 100 mM
- ad 100mL LB medium (→ pH 6.5)

### Transformation

#### Step 2.

- Inoculate starter culture of 1.5 mL LB medium
- After 1.5 - 2 hours culture should become turbid (→ exponential growth)

### Transformation

#### Step 3.

- Add 200 µL of culture to 1.5 mL tube and mix with 200 µL 2xTSS (work on ice)
- Add 0.5 - 1.0 µL plasmid (high copy: rather 0.3 µL)

### Incubation

#### **Step 4.**

- Incubate on ice for 20-30 min
- Incubate 45-60 min at 37°C on shaking device (> 200 rpm) – addition of LB is not requested
- Plate on LB-Agar (incl. appropriate antibiotic)