



Feb 25, 2019

In devel.

scTHS-seq

Kun Zhang¹, Brandon Sos¹, Dinh Diep¹, Thu E Duong¹¹University of California, San Diego

dx.doi.org/10.17504/protocols.io.xf8fjrw

 Dinh Diep
University of California, San Diego

ABSTRACT

Single cell transposome hypersensitive site sequencing (scTHS-seq) combined transposome hypersensitive site sequencing (THS-seq) with combinatorial cellular indexing using customized barcoded transposomes. Compared to a similar method, ATAC-seq, THS-seq have improved coverage of distal enhancer sites which is due to the utilization of an engineered super-mutant Tn5 transposase as well as in vitro transcription to perform linear amplification from a single insertion site of the Tn5 transposase instead of two insertion sites required in ATAC-seq. scTHS-seq was used to generate single nucleus accessibility profile of 36,869 cells with an average of 10,168 unique reads per cell that are also associated with a predicted accessible region. The doublet rate of scTHS-seq is estimated to be between 20-25% which is comparable to similar combinatorial indexing methods.

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

<https://www.nature.com/articles/nbt.4038>

PROTOCOL STATUS

In development

We are still developing and optimizing this protocol

MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
HiScribe T7 High Yield RNA Synthesis Kit - 50 rxns	E2040S	New England Biolabs
Taq 5X Master Mix - 500 rxns	M0285L	New England Biolabs
Pierce™ Dimethylformamide (DMF), Sequencing grade	20672	Thermo Fisher Scientific
2x Kapa HiFi Hotstart Readymix	KK2602	Kapa Biosystems
Tango Buffer	BY5	Thermo Fisher Scientific
Ultrapure SMART MMLV Reverse Transcriptase for RT-PCR	639524	Clontech
Advantage® UltraPure dNTP Combination Kit (100 mM each dNTP)	639125	Clontech
RNAse H	Y9220L	Enzymatics
KAPA Pure Beads (60 mL)	07983298001	Roche
twin.tec PCR Plate 96 LoBind semi-shirted clear 25 pcs.	30129504	Eppendorf

STEPS MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
KAPA Pure Beads (60 mL)	07983298001	Roche
Taq 5X Master Mix - 500 rxns	M0285L	New England Biolabs
KAPA Pure Beads (60 mL)	07983298001	Roche
KAPA Pure Beads (60 mL)	07983298001	Roche

NAME ▾	CATALOG # ▾	VENDOR ▾
KAPA Pure Beads (60 mL)	07983298001	Roche
KAPA Pure Beads (60 mL)	07983298001	Roche
KAPA Pure Beads (60 mL)	07983298001	Roche

MATERIALS TEXT

For purification of engineered super-mutant Tn5-059, see previous publications

(<https://genomebiology.biomedcentral.com/articles/10.1186/s13059-016-0882-7>) for details.

Tagmentation

- 1 Make 10 mL 75% glycerol with 0.3% tritonX-100

▢ 9.45 g Glycerol (100%)

▢ 300 µl TritonX-100, 10%

▢ 2200 µl nuclease free water

- 2 Make 3.045 mL transposons dilution buffer

▢ 2345 µl 75% glycerol/0.3% TritonX-100

▢ 700 µl Qiagen EB

- 3 Make 1000 uL transposase dilution buffer

▢ 667 µl 75% glycerol/0.3% TritonX-100

▢ 1 µl DTT, 1M

▢ 6.25 µl Tris-HCl, 4M, pH 7.5

▢ 325.75 µl nuclease free water

- 4 Make 10 mL of 1.5X D6 buffer

▢ 1500 µl Tango Buffer, 10X

▢ 2400 µl DMF, 100%

▢ 6100 µl nuclease free water

- 5 Make 10 mL of 1X lysis buffer

▢ 25 µl Tris-HCl, 4M, pH 7.5

▢ 20 µl Sodium chloride, 5M

▢ 15 µl Magnesium chloride, 2M

▢ 100 µl IGEPAL-630, 10%

▢ 2000 µl BSA, 10%

▢ 1000 µl PBS, 10X

 **6840 µl nuclease free water**

Finally, add 1 protease inhibitor tablet (EDTA-free).

6 Make 10 mL of 2X FACS buffer

 **1000 µl PBS, 10X**

 **2000 µl BSA, 10%**




 **80 µl EDTA, 0.5M**

 **6920 µl nuclease free water**

7 Generate 384 annealed transposons

 **7 µl Indexed Tn5 top strand (100uM)**

 **7 µl 5' phosphorylated mosaic end sequence bottom strand (100uM)**

Combine top and bottom strand oligos in a PCR tube. In a thermocycler, incubate for 2 minutes at  **95 °C** then ramp to  **14 °C** at a rate of  **-0.1 °C** per second.

8 Dilute annealed transposons

 **1 µl Annealed transposons**

 **6.7 µl transposon dilution buffer**

9 Make 420 uL of 4.8 uM transposase Tn5-059

 **252 µl transposase, Tn5-059, 8 uM**

 **168 µl transposase dilution buffer**

10 Load Tn5-059

 **1 µl 4.8 uM transposase Tn5-059**

 **1 µl diluted annealed transposons (6.5 uM)**








11 Incubate

 **22 °C Thermocycler**  **00:30:00**



12 Wash nuclei or cells suspension with ice-cold PBS



13 Count cells

14 Pellet 400,000 cells / nuclei by centrifuging at 330 RCF (cells) or 900 RCF (nuclei)









- 15 Resuspend 400,000 cells/nuclei in lysis buffer. Use 100 uL per 100,000 cells/nuclei
- 16 Incubate
-  **4 °C On ice**  **00:05:00**
- 17 Pellet by centrifuging at 900 RCF
- 18 Resuspend nuclei in 100 uL 1.5X D6 buffer initially. Count and add more 1.5X D6 buffer to create ~250 nuclei/uL
- 19 Aliquot 4.0 uL of cells in 1.5X D6 buffer into wells containing 2.0 uL of prepared transposomes
- 20 Incubation without shaking
-  **37 °C Thermocycler**  **00:30:00**
- 21 Terminate the tagmentation reaction with EDTA
-  **4 µl EDTA, 50 mM**
- 22 Incubate
-  **37 °C Thermocycler**  **00:15:00**
- 23 Store cells at -20C overnight or until ready to sort

FACS (BD Influx)

- 24 Coat 12 x 200 uL, 1 x 15 mL conical using 1X FACS buffer
-  **01:00:00 Let tubes sit in FACS buffer**
- 25 Add 10 uL PBS to each well of 10 to 15 x 96-wells plate for sorting
-  Cover and store all plates on ice.
- 26 Add 1 volume (10 uL) ice cold 2X FACS buffer to each sample well with cells in tagmentation mix
- 27 When 15 mL conical have been coated for 1 hr, pool all samples together using multichannel pipette and adding to 12 x 200 uL tubes

- 28 Centrifuge at 500 RCF  4 °C
- 29 Resuspend nuclei in 1.5 mL 1X FACS buffer
- 30 Add stain for sorting
-  75 µl **Propidium Iodide**
- 31 Sort 100 cells into each well using BD Influx, remove doublets using forward and side scatter plots and selecting for PI-staining events


In vitro transcription (10 plates of sorted nuclei)

- 32 Add 11 µL of 8M guanidine chloride to each well using multi-channel to lyse nuclei and release DNA
-  14 ml **Guanidine hydrochloride, 8M**
- 33 Add 30 µL of KAPA Pure Beads to each well to bind DNA
-  42.24 ml **KAPA Pure Beads**
-  KAPA Pure Beads (60 mL)
by Roche
Catalog #: 07983298001
- 34 Mix well and bind for at least 8 minutes
-  00:08:00 **Room temperature**
- 35 Place plate on magnet for at least 3 minutes then discard the supernatant
-  00:03:00 **Room temperature**
- 36 Prepare fresh 80% ethanol solution. Wash beads (while still on magnet) twice with 170 µL of 80% ethanol
-  384 ml **Freshly prepared 80% ethanol**
- 37 After last wash, completely remove all ethanol and transfer to a heated thermomixer to dry
-  37 °C  00:01:00
- 38 Add 8 µL of nuclease free water to each well using a multichannel pipette

 **8640 µl nuclease-free water**

39 Add 2 µL of 5X Taq Master Mix to perform gap filling

 **2160 µl 5X Taq Master Mix**



Taq 5X Master Mix - 500 rxns
by [New England Biolabs](#)
Catalog #: [M0285L](#)

40 Incubate plates for gap filling

 **72 °C**  **00:08:00**

41 Add 2 µL of HiScribe T7 IVT reaction components to each well

 **2112 µl T7 buffer, 10X**

 **2112 µl UTP**

 **2112 µl GTP**

 **2112 µl ATP**

 **2112 µl CTP**


 **2112 µl T7 RNA Polymerase**

42 Incubate

 **37 °C Thermocycler**  **16:00:00**

43 Add 40 µL of KAPA Pure Beads to each well to bind DNA

 **42.24 ml KAPA Pure Beads**



KAPA Pure Beads (60 mL)
by [Roche](#)
Catalog #: [07983298001](#)

44 Mix well and bind for at least 8 minutes

 **00:08:00 Room temperature**

45 Place plate on magnet for at least 3 minutes then discard the supernatant

 **00:03:00 Room temperature**

46 Prepare fresh 80% ethanol solution. Wash beads (while still on magnet) twice with 170 uL of 80% ethanol

 **384 ml** Freshly prepared 80% ethanol

47 After last wash, completely remove all ethanol and transfer to a heated thermomixer to dry

 **37 °C**  **00:01:00**

48 Add 9 uL of nuclease free water to each well using a multichannel pipette

 **9720 µl** nuclease-free water

cDNA synthesis (10 plates of sorted nuclei)

49 Add 2.5 uL of random hexamers

 **3000 µl** random hexamers, 20 uM

50 Incubate

 **70 °C Thermocycler**  **00:03:00**



Cool immediately on ice after.

51 Add 9.5 uL of MMLV RT mixture

 **4224 µl** First strand buffer, 5X

 **2112 µl** Advantage UltraPure dNTPs, 100 uM each

 **2112 µl** DTT, 100mM

 **528 µl** MMLV RT

52 Incubate in a thermocycler

 **22 °C**  **00:10:00**

 **42 °C**  **01:00:00**

 **70 °C**  **00:10:00**

53 Add 1 uL RNase H

 **1127.52 µl** nuclease-free water

 **125.28 µl** RNase H, 5 Units/uL

Incubate in thermocycler

54

 37 °C  00:20:00

55 Add 2.5 uL second strand synthesis primer

 3000 µl sss_scnXTv2, 20 uM

56 Incubate in thermocycler

 65 °C  00:02:00

Cool immediately on ice after incubation

57 Add 5.9 uL 5X Taq Master Mix

 7363 µl Taq Master Mix, 5X

58 Incubate in thermocycler

 72 °C  00:08:00

Tagmentation on beads (10 plates of sorted nuclei)

59 Add 40 uL of KAPA Pure Beads to each well to bind DNA

 42.24 ml KAPA Pure Beads

KAPA Pure Beads (60 mL)

by Roche

Catalog #: 07983298001

60 Mix well and bind for at least 8 minutes

 00:08:00 Room temperature

61 Place plate on magnet for at least 3 minutes then discard the supernatant

 00:03:00 Room temperature

62 Prepare fresh 80% ethanol solution. Wash beads (while still on magnet) twice with 170 uL of 80% ethanol

 384 ml Freshly prepared 80% ethanol

63 After last wash, completely remove all ethanol and transfer to a heated thermomixer to dry

 37 °C  00:01:00

64 Add 4 uL of 1.5X tagmentation buffer to each well using a multichannel pipette

 3172 µl nuclease-free water




 780 µl Tango Buffer, 10X

 1248 µl DMF, 100%

65 Generate annealed transposons

 40 µl i7 top strand (100uM)

 40 µl 5' phosphorylated mosaic end sequence bottom strand (100uM)

Combine top and bottom strand oligos in a PCR tube. In a thermocycler, incubate for 2 minutes at  95 °C then ramp to  14 °C at a rate of  -0.1 °C per second.

66 Dilute annealed transposons

 74 µl Annealed transposons

 1026 µl transposon dilution buffer

67 Make 420 uL of 2.4 uM transposase Tn5-059

 330 µl transposase, Tn5-059, 8 uM

 770 µl transposase dilution buffer

68 Load Tn5-059

 1100 µl 2.4 uM transposase Tn5-059

 1100 µl diluted annealed transposons (6.5 uM)

69 Incubate

 22 °C Thermocycler  00:30:00

70 Add 2 uL of loaded transposomes to each well

71 Incubate

 **55 °C Thermocycler**  **00:06:00**



Cool immediately on ice after.

72 Add 6 uL guanidine hydrochloride

 **7488 µl Guanidine hydrochloride, 8M**

73 Add 30 uL of KAPA Pure Beads to each well to bind DNA

 **31.68 ml KAPA Pure Beads**



KAPA Pure Beads (60 mL)
by Roche
Catalog #: 07983298001

74 Mix well and bind for at least 8 minutes

 **00:08:00 Room temperature**

75 Place plate on magnet for at least 3 minutes then discard the supernatant

 **00:03:00 Room temperature**

76 Prepare fresh 80% ethanol solution. Wash beads (while still on magnet) twice with 170 uL of 80% ethanol

 **384 ml Freshly prepared 80% ethanol**

77 After last wash, completely remove all ethanol and transfer to a heated thermomixer to dry

 **37 °C**  **00:01:00**

78 Add 16 uL nuclease free water

 **21.6 ml nuclease-free water**

Test amplification (10 plates of sorted nuclei)

79 Preheat hotstart polymerase mix

 **250 µl KAPA HiFi Hotstart ReadyMix, 2X**

 **98 °C Thermocycler**  **00:00:45**



Cool immediately on ice

80 Add primers to make a master mix to 417 nM final each

 **42 µl i5XX, 5 uM**

 **42 µl i7XX, 5 uM**

81 Add SYBR Green for qPCR to 0.8X final

 **4 µl SYBR Green, 100X**

82 Aliquot 8 uL of amplification mix to each well to qPCR plate

83 Sample 5 wells from each plate, adding 4 uL of each sample. Include an no nuclei control well.

84 Perform qPCR. Monitor the sample progress, terminate when all samples appears to be in exponential phase (usually 6-8 cycles). No nuclei control may appear to amplify but majority of those products are primer dimers.

1. 72C x 5 min
2. 98C x 30 s
3. 59C x 30 s
4. 72C x 1 min
- Plate read
5. Go to 2 for 20 cycles
6. 72C for 5 min

85 Assess PCR products on a 6% TBE PAGE gel. If lack of products, repeat test amplification on a different sampling of wells but adding 3-4 more cycles of amplification. If samples appear to be over-amplified, reduce the number of cycles by 1.

Library amplification (10 plates of sorted nuclei)

86 Preheat hotstart polymerase mix

 **12000 µl KAPA HiFi Hotstart ReadyMix, 2X**

 **98 °C Thermocycler**  **00:00:45**



Cool immediately on ice

87 Aliquot 12 uL of preheated KAPA HiFi ReadyMix to each well of a new PCR plate

88 Aliquot 2 uL of 16 different s5XX primers to each set of 2 plates (one per row)

 **150 µl each s5XX primer, 5 uM**

89 Aliquot 2 uL of 12 different s7XX primers to each column of every plate (one per column)

 **100 µl s7XX primer, 5 uM**

90 Transfer 8 uL of sample to each well of PCR plate

91 Perform PCR using the number of cycles determined previously with test PCR

1. 72C x 5 min
2. 98C x 30 s
3. 59C x 30 s
4. 72C x 1 min
5. Go to 2 for XX cycles
6. 72C for 5 min

92 Pool 4 uL of each well from 2 plates with non-overlapping s5XX primers. Final volume is ~800 uL per pool. 5 pools total

93 Add 800 uL of KAPA Pure Beads to each pool to bind DNA

 **3 ml KAPA Pure Beads**



KAPA Pure Beads (60 mL)
by Roche
Catalog #: 07983298001

94 Mix well and bind for at least 8 minutes

 **00:08:00 Room temperature**

95 Place plate on magnet for at least 3 minutes then discard the supernatant

 **00:03:00 Room temperature**

96 Prepare fresh 80% ethanol solution. Wash beads (while still on magnet) twice with 1000 uL of 80% ethanol

 **12 ml Freshly prepared 80% ethanol**

97 After last wash, completely remove all ethanol and transfer to a heated thermomixer to dry

 37 °C  00:01:00

98 Add 20 uL of nuclease free water to each pool

 100 µl **nuclease-free water**

99 Quantify each pool using dsDNA HS Qubit

100 Perform gel size selection for 220-1000 bp using 6% TBE PAGE. Use 5-wells gel, and loading ~1 ug of library per well.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited