

Recipe for 50x TAE buffer

Anna Behle, Alice Pawlowski

Abstract

Stock solution for 50x TAE. TAE buffer is a solution made up of Tris base, acetic acid and EDTA (Tris-acetate-EDTA). It is a common buffer for DNA separation using standard agarose gel electrophoresis.

Citation: Anna Behle, Alice Pawlowski Recipe for 50x TAE buffer. **protocols.io**

dx.doi.org/10.17504/protocols.io.gtvbwn6

Published: 06 Apr 2018

Materials

✓ EDTA by Contributed by users

Acetic acid, glacial [537020](#) by [Sigma Aldrich](#)

Tris (Tris Base) [T-400](#) by [Gold Biotechnology](#)

Protocol

Components required

Step 1.

Component	Molarity in 50x	Molecular weight
EDTA disodium salt	50 mM	372.24 g/mol
Tris	2 M	121.14 g/mol
glacial / acetic acid	1 M	60.05

Preparation of a 0.5 M EDTA stock solution

Step 2.

For 500 ml:

- weigh out 93.05 grams of EDTA disodium salt (MW=372.24 g/mol)
- Dissolve in 400 milliliter deionized water and adjust the pH with solid sodium hydroxide (NaOH) plates, EDTA will not go completely into solution until the pH is adjusted to about 8.0!
- Top up the solution to a final volume of 500 milliliters
- autoclave

Preparation of 50 x TAE buffer

Step 3.

- weigh out 242 grams of Tris-base (MW = 121.14 g/mol) and dissolve in approximately 700 milliliters of deionized water
- Carefully add 57.1 milliliters of 100 % glacial acid (or acetic acid) and 100 milliliters of 0.5 M EDTA (pH 8.0)
- adjust the solution to a final volume of 1 liter
- the pH of this buffer is not adjusted and should be about 8.5
- store stock solution at room temperature

Preparation of 1 x TAE working solution

Step 4.

20 ml 50 x TAE

ad 1000 ml a. dest.

final concentration in gel / running buffer: 40 mM Tris, 20 mM acetic acid, 1 mM EDTA