

# TCBS Agar Plates (500ml)

Sarah Hessen-Schmidt

## Abstract

Procedure for making TCBS Agar plates for selection of Vibrio bacteria

**Citation:** Sarah Hessen-Schmidt TCBS Agar Plates (500ml). **protocols.io**

dx.doi.org/10.17504/protocols.io.e7fbhjn

**Published:** 21 Jun 2016

## Materials

 TCBS Agar C7041 by Contributed by users

## Protocol

### Step 1.

In 500ml MilliQ, add 44g TCBS Agar. Mix thoroughly

 **AMOUNT**

44 g Additional info:

 **REAGENTS**

 TCBS Agar C7041 by Contributed by users

### Step 2.

Loosen cap. Bring to a boil in microwave to dissolve completely and boil for approximately 1 minute.  
DO NOT AUTOCLAVE

 **DURATION**

00:01:00

### Step 3.

Allow to cool to approximately 60°C by running container until MilliQ

### Step 4.

Under a flame, pour TCBS agar into sterile plates about half way. Quickly run flame over surface of agar to remove any air bubbles.