

Peroxidase-linked lectin assay to determine neuraminidase-inhibiting antibodies using reassortant influenza viruses

Version 2

Yulia A. Desheva, Tatiana A. Smolonogina, Ivan A. Sychev

Abstract

The NA inhibition test using large molecular fetuin substrate is believed to be the most specific method to detect antibodies against influenza virus neuraminidase. To avoid cross-reacting with hemagglutination-inhibition antibodies we have prepared several reassortant A/H7N1 viruses containing the NA of A/California/07/09(H1N1)pdm, A/New Caledonia/20/1999(H1N1) or A/Vietnam/1203/04(H5N1), and the hemagglutinin derived from A/equine/Prague/1/56(H7N7) influenza virus. Adjustment and validation of the NI test with the reassortant viruses included the determination of the viruses working concentration and selection of optimal conditions for the enzymatic reaction.

Citation: Yulia A. Desheva, Tatiana A. Smolonogina, Ivan A. Sychev Peroxidase-linked lectin assay to determine neuraminidase-inhibiting antibodies using reassortant influenza viruses. **protocols.io**

dx.doi.org/10.17504/protocols.io.n3adgie

Published: 27 Mar 2018

Before start

Preparation of reagents dilutions

1. Prepare a stock solution of fetuin at 10 mg/ml in coating buffer and store in 500 µl aliquots at -20 °C.
2. Prepare a working solution of fetuin (50 µg/ml) immediately before coating plates by diluting the stock solution 200-fold in coating buffer.
3. Prepare a stock solution of BSA at 100 mg/ml in sterile PBS (sPBS) and store in 10-25 ml aliquots at -20°C.
4. The diluent for virus and sample dilutions (PBS-BSA). Prepare the 1% work solution of BSA (10 mg/ml) right before assay using PBS as a diluent. After thawing aliquot with stock BSA, you may storage (1-2 months) it in a refrigerator at +4°C for next assays.
5. To prepare a stock solution of lectin, mix 1 mg of lectin and 5 ml of sPBS and store at +4°C.
6. Right before use, prepare a work solution of lectin at 2.5 µg/ml by diluting the stock solution 80-fold in PBS-BSA.
7. Prepare a work solution of RDE from original solution by 10-fold dilution at PBS-BSA.
8. Prepare a sufficient volume of wash buffer – 0.5% PBS-Tween 20 (T-PBS). Store at +4°C.
9. Prepare the peroxidase substrate TMB. Place at room temperature for 30 min before use to warm up.
10. Prepare the stop solution - 1 N H₂SO₄ - add 27.2 ml stock 98% H₂SO₄ to 973 ml distilled H₂O.

Serum samples

1. Heat all sera in a water bath at 56°C for 30 min. Store the sera at -20°C before and after the heat treatment.

Materials

- TMB Substrate Reagent Set [555214](#) by [BD Biosciences](#)
- Fetuin from fetal bovine serum [F3004](#) by [Sigma - Aldrich](#)
- High-binding 96-well microplates [655061](#) by [greiner bio-one](#)
- Bovine serum albumin (BSA) [A7030](#) by [Sigma - Aldrich](#)
- Peroxidase-conjugated lectin from Arachis hypogaea [L7759](#) by [Sigma - Aldrich](#)
- Receptor destroying enzyme (RDE) [View](#) by [Denka Seiken Co., Ltd](#)
- ✓ Sulfuric Acid (H₂SO₄) by Contributed by users
- ✓ 0,1M carbonate-bicarbonate buffer with pH=9.5-9.7 by Contributed by users
- ✓ 0,01M phosphate-buffered saline (PBS) with pH=7.3-7.5 by Contributed by users

Protocol

Step 1.

Adjustment of the virus working dose to use in enzyme-linked lectin assay.

Step 2.

Prepare fetuin-coated plates. Coat 96-well plate with 150 µl of the working solution of fetuin and incubate it at +4°C overnight (at least 18h).

AMOUNT

150 µl Additional info: per well

TEMPERATURE

4 °C Additional info:

NOTES

Yulia Desheva 07 Mar 2018

Still, It is possible to use plates after 2-3 days of incubation.

Step 3.

Prepare viruses dilutions. In 3 short rows of 96-well polymer plates for immunologic reactions with a U-shaped bottom serially two-fold dilutions of the interested virus at PBS-BSA. Titration starts from 130 µl of the virus at concentration of 1024 hemagglutination units (HAU) per 50 µl. Add 65 µl of PBS-BSA to have final volume 130 µl. Use 130 µl PBS-BSA in 4 wells as negative control. Use 130 µl of the work solution of RDE in 4 wells as positive control. Incubate for 30 min. Place the plate in a humidified incubator at 37°C. This step imitating contact of the virus with serum.

TEMPERATURE

37 °C Additional info: in a humidified incubator

Step 4.

Transfer samples with virus to the fetuin-coated plate. While contact is going, wash fetuin-coated plate 2x with sPBS 200 ml/well then invert and pat the plate on absorbent paper towels. After contact, transfer the virus dilutions into the corresponding well rows of a plate with sorbed fetuin using a multichannel pipette of 100 µl starting with the highest concentration. Place the plate in a humidified incubator at 37°C for 1h.

TEMPERATURE

37 °C Additional info: in a humidified incubator

NOTES

Yulia Desheva 07 Mar 2018

Primarily, It is better to transfer fluid from control wells, then change pipette tips and transfer the virus starting from highest dilution.

Step 5.

Add the lectin from *Arachis hypogaea*. When the incubation is complete discard liquid from wells. Wash 4x with T-PBS 200 ml/well. Add to each well 100 µl lectin. Place the plate at room temperature for 1h.

AMOUNT

100 µl Additional info: per well

TEMPERATURE

20 °C Additional info: at RT

Step 6.

Complete the assay. When the incubation is complete discard liquid from wells. Wash 4x with T-PBS 200 ml/well. Add to each well 100 µl TMB. Incubate for 2-3 min (maybe more or less) at room temperature. Stop the reaction by adding to each well 100 µl 1 N H₂SO₄. Measure the optical density (OD) at 450 nm using the universal plate reader.

Step 7.

Data analysis. Define the yield of the reaction product (RY) as follows: $RY = (OD-VD - OD-NC)/(OD-PC - OD-NC) \times 100\%$, where OD-VD – is a mean of OD from wells with virus dilution; OD-NC – is a mean of OD from wells with negative control; OD-PC – is a mean of OD from wells with positive control. Determination of the range of linear dependence of the yield of the sialidase reaction product on a concentration of a virus.

NOTES

Yulia Desheva 07 Mar 2018

Usually this range of linear dependence lasts from 8 HAU/50µl to 128 HAU/50µl. 128 HAU/50µl is the upper limit of a working concentration of a virus.

Step 8.

Enzyme-linked lectin assay.

Step 9.

Prepare fetuin-coated plates. See step 2.

Step 10.

Prepare sample dilutions. In short rows of 96-well polymer plates for immunologic reactions with a U-shaped bottom in 65 µl, serial two-fold dilutions (from 1/10 to 1/640) of the serum samples to be tested were made at PBS-BSA. In each plate, the same positive and negative reference sera are titrated. NB! Don't forget to vortex all sera samples before adding it in plates. Prepare a dilution of a virus at 128 HAU in 50 µl by PBS-BSA. Add 65 µl of the diluted virus into each well. The last row of 96-well plates (H-row) is used for positive and negative virus control (65 µl of the diluted virus + 65 µl PBS-BSA and 130 µl PBS-BSA respectively). Incubate for 30 min. Place the plate in a humidified incubator at 37 °C.

TEMPERATURE

37 °C Additional info: in a humidified incubator

Step 11.

Complete the assay by following steps 4 – 6.

Step 12.

Data analysis. The result of the enzymatic reaction estimated as OD data is determined for a series of dilutions of each blood serum. Use the OD set to plot the residual NA enzymatic activity in the presence of anti-NA antibodies, calculated according to the formula below: $\text{Activity} = (\text{OD}_{\text{sample}} - \text{OD}_{\text{NC}}) / (\text{OD}_{\text{PC}} - \text{OD}_{\text{NC}}) * 100\%$, where: OD-PC and OD-NC are a mean of OD from wells with positive control and negative control. Determine the titer of serum anti-NA antibodies as the reciprocal dilution of the sample causing 50% inhibition of NA activity. To calculate the 50% inhibitory concentration (IC₅₀) use regression analysis. Rough nonlinear regression curve fitting can be performed by Microsoft-Excel. Precise curve fitting can be performed by Wolfram-Alpha web-tool.

NOTES

Yulia Desheva 08 Feb 2018

If IC₅₀ was not achieved at any dilution, the titer will be less than the first dilution. If the first dilution is 1:10, then the titer will be <1:10.

Step 13.

Warnings

All live reassortant viruses must be handled using biosafety level 2 (BSL2)-enhanced practices in a laboratory.