

Sodium Dodecyl Sulfate PolyAcrylamide Gel Electrophoresis (SDS-PAGE) with Sypro straining

Audrey Dubourg, Dong Xia, John P Winpenny, Suha Al Naimi, Maha Bouzid, Darren W Sexton, Jonathan M Wastling, Paul R Hunter, Kevin M Tyler

Abstract

This protocol provides an efficient way to observe proteins in samples.

Citation: Audrey Dubourg, Dong Xia, John P Winpenny, Suha Al Naimi, Maha Bouzid, Darren W Sexton, Jonathan M Wastling, Paul R Hunter, Kevin M Tyler Sodium Dodecyl Sulfate PolyAcrylamide Gel Electrophoresis (SDS-PAGE) with Sypro straining. **protocols.io**

dx.doi.org/10.17504/protocols.io.mcdc2s6

Published: 04 Jan 2018

Before start

Prepare acrylamide gels a few hours before running the gels; acrylamide percentage depends on protein sizes. Prepare SDS-PAGE buffers beforehand as well.

Protocol

Protein sample preparation

Step 1.

Load proteins with loading buffer (containing loading dye) in a 3:1 ratio in Eppendorf tubes, then boil samples for 5 min at 100°C.

Protein migration by electrophoresis

Step 2.

Load 20 µl of protein samples and a protein standard molecular weight marker in gels placed in a vertical gel tank system filled with 1x running buffer. Run gels at 100 volts for 2 hours at room temperature.

📄 AMOUNT

20 µl Additional info: Protein samples

Protein fixation on gel

Step 3.

After electrophoresis, immerse gels in a fixing solution bath for 30 min on a rotary shaker. (fixation 1/2)

Protein fixation on gel

Step 4.

Repeat the fixation one more time; immerse gels in a fixing solution bath for 30 min on a rotary shaker. (fixation 2/2)

Protein staining with SYPRO solution

Step 5.

Immerse gels in 40 ml of pre-diluted SYPRO solution and incubate overnight, protected from the light on a rotary shaker at room temperature.

AMOUNT

40 ml Additional info: Pre-diluted SYPRO

Protein wash with washing buffer

Step 6.

Wash gels with a washing solution for 30 min on a rotary shaker at room temperature, in a container protecting gels from the light.

SAFETY INFORMATION

Collect SYPRO solution in a specific container waste and follow MSDS and health and safety regulation to discard.

Protein wash with ultrapure water

Step 7.

Immerse gels in ultrapure water and wash for 5 min at room temperature on a rotary shaker, in a container protecting gels from the light. (wash 1/2)

SAFETY INFORMATION

Collect washing buffer in a properly labeled container and follow health and safety regulation to discard.

Protein wash with ultrapure water

Step 8.

Immerse gels in ultrapure water and wash for 5 min at room temperature on a rotary shaker, in a container protecting gels from the light. (wash 2/2)

SAFETY INFORMATION

Collect washing buffer in a properly labeled container and follow health and safety regulation to discard.

Reading gels

Step 9.

Read gels under UV machine.

SAFETY INFORMATION

Wear required protecting gear and discard gels in a acrylamide-specific container to be discarded following health and safety regulation.

Warnings

Acrylamide is carcinogenic and should be handled and disposed of accordingly to MSDS sheet and Health and Safety regulation.