

# Overnight Bacterial Batch Culture

Nicholas Boudreau

## Abstract

Purpose of this protocol is to grow up single colonies of bacteria for use.

**Citation:** Nicholas Boudreau Overnight Bacterial Batch Culture. **protocols.io**

[dx.doi.org/10.17504/protocols.io.nv3de8n](https://doi.org/10.17504/protocols.io.nv3de8n)

**Published:** 17 Mar 2018

## Protocol

### Step 1.

Pipette 5 mL of LB Broth into bacterial culture tube.

### Step 2.

Add desired antibiotic at appropriate concentration (often 1:1000).

### Step 3.

Pick single colony/scrape glycerol stock and inoculate LB Broth.

### Step 4.

Incubate batch culture overnight at 37°C, shaking.