

# Transformation of supercompetent cells

Dave Lunt

## Abstract

Protocol to prepare *E. coli* supercompetent cells to transform with plasmid/ligation. This protocol is originally derived from Hanahan, D. (1983) J. Mol. Biol. 166:557-580 with minor changes. This version works exceptionally well for cloning PCR products where the number of colonies is not expected to be as high as from some ligations. These cells outperform commercially purchased supercompetent cells in my hands.

**Citation:** Dave Lunt Transformation of supercompetent cells. **protocols.io**

dx.doi.org/10.17504/protocols.io.ettbenn

**Published:** 04 Apr 2016

## Before start

### Prepare for day 2:

1. 42°C waterbath
2. 37°C waterbath
3. prechilled transformation tubes
4. 37°C SOC
5. LB amp<sup>+</sup> Xgal agar plates

## Buffers

### TFB Maniatis 1.78

10mM K-MES (pH 6.2) adjust with KOH  
100mM KCl  
45mM MnCl<sub>2</sub>  
10mM CaCl<sub>2</sub>  
3mM H<sub>3</sub>CoCl<sub>3</sub> (Haxamine cobalt III chloride)  
filter sterilize and store at 4°C (lasts a very long time)

### SOB

2% bacto-tryptone (20g l<sup>-1</sup>)  
0.5% yeast extract (5g l<sup>-1</sup>)  
10mM NaCl  
2.5mM KCl  
10mM MgCl<sub>2</sub>  
10mM MgSO<sub>4</sub>

SOC is SOB with 20mM glucose

## Protocol

### Step 1.

Inoculate two separate 1.5mls of LB medium with Invitrogen (InvαF') cells. Also include an uninoculated negative control tube. Shake at 250 rpm overnight at 37°C.

#### 🔗 NOTES

**Dave Lunt** 04 Apr 2016

These cells are similar to any *E. coli* sold for transformation. Inoculation is by scraping a pipette tip over the top of a microfuge tube of frozen glycerol stock and swirling in LB medium.

**Dave Lunt** 04 Apr 2016

Use two starter cultures in case one doesn't grow, which is rare but loses you a day. Include a negative control to check for contamination.

### Step 2.

Inoculate with 500ul of overnight *E. coli* culture two 50ml tubes each containing 25mls of sterile SOB media. Shake at 250 rpm for 4 hours

#### 🔗 NOTES

**Dave Lunt** 04 Apr 2016

Two replicate 25ml cultures are prepared, if one doesn't grow you haven't lost a day

**Dave Lunt** 04 Apr 2016

Use tape to keep tube lids loosely in place, but make sure you don't seal the tubes, they need air

### Step 3.

Ice the tubes of cells for 15 mins

### Step 4.

Pellet cells (approx 5 min spin, approx 3000 rpm) and discard supernatant

### Step 5.

Resuspend each pellet of cells in 8 mls of cold TFB by washing gently, then sucking up and down, with a wide-bore 10ml pipette tip

### Step 6.

Ice the tubes of cells for 15 mins

### Step 7.

Pellet cells (approx 5 min spin, approx 3000 rpm) and discard supernatant

### Step 8.

Resuspend one pellet of cells in 2 mls of cold TFB and use this liquid to resuspend the second pellet. You should now have all cells in a single aliquot of 2mls cold TFB.

### Step 9.

Add 105ul of DMSO swirl and ice for 5 mins.

### Step 10.

Add 80ul of 100mM DTT, swirl, and ice for 10 mins.

### Step 11.

Add 105ul of DMSO swirl and ice for 5 mins.

### Step 12.

Aliquot 180ul of cells into pre-chilled, labelled 5ml polypropylene growth tubes with loose-cap lids. Include positive and negative controls.

### Step 13.

Add ligation DNA to each sample (in 10ul volume) and ice for 30 mins

📌 **NOTES**

**Dave Lunt** 04 Apr 2016

Now would be a great time to check the waterbath is at exactly 42°C and that you have some warm SOC

**Step 14.**

Heat pulse tubes at exactly 42°C for exactly 90 secs.

📌 **NOTES**

**Dave Lunt** 04 Apr 2016

Use a floating tube rack that can be moved rapidly between waterbath and ice bucket.

**Step 15.**

Ice for 2 mins

**Step 16.**

Add 400ul of 37°C SOC to each tube and incubate at 37°C for 10 mins

**Step 17.**

Incubate a further 45 mins at 37°C shaking at 225 rpm

📌 **NOTES**

**Dave Lunt** 04 Apr 2016

Some think this step is unnecessary

**Step 18.**

Place tubes in ice to halt growth

**Step 19.**

Transfer the cells from each tube to a labelled individual 1.5ml microfuge tube

**Step 20.**

Give the microfuge tubes a very short (5 sec) spin to very loosely pellet the cells and then remove and discard 370ul of the supernatant

📌 **NOTES**

**Dave Lunt** 04 Apr 2016

This step increases the cell density, omit this step if you expect a lot of colonies anyway

**Step 21.**

Plate 40ul of the cell solution onto LB/ Amp+/ Xgal plates and grow overnight at 37°C

📌 **NOTES**

**Dave Lunt** 04 Apr 2016

InvαF' cells do not need IPTG

**Dave Lunt** 04 Apr 2016

Spread 3 replicate plates with 40ul as redundancy is good