

TFB buffer

James Thornton

Abstract

Citation: James Thornton TFB buffer. **protocols.io**

dx.doi.org/10.17504/protocols.io.eq9bdz6

Published: 13 Jun 2016

Guidelines

Components:

10 mM K-MES, pH 6.2, 100 mM KCl, 45 mM $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$, 10 mM $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$, 3 mM Na_2CO_3 .

Protocol

Step 1.

Weigh out the components.

Step 2.

Add the MES to d- H_2O and adjust the pH with KOH.

Step 3.

Add the remaining components (in order), waiting until one component is in solution before adding the next.

Step 4.

Adjust the volume to the final volume.

Step 5.

Filter sterilize.

Step 6.

Store frozen at -20°C in 15 ml aliquots.