

# Preparation of Genomic DNA of Microsporidia

ShiNan Dong

## Abstract

To prepare the genomic DNA of microsporidia, 400  $\mu$ l suspension of purified spores (1010 spores/ml) was mixed with 40  $\mu$ l KOH (2 mol/l) in a 1.5 ml Eppendorf tube and incubated at 27°C for 1 h, TEK buffer (1 mmol/l Tris-HCl, 10 mmol/l EDTA, 0.17 mol/l KCl, pH 8.0) was added and continued incubating at 27°C for 1 h. Adjusted pH to 8.0 with 1 mol/l HCl, added 10% SDS in order to attain 0.5% in the mixture and kept it in ice-bath for 15 min. 20 mg/ml proteinase K (TaKaRa Biotechnology Co. LTD) was added in order to attain 200  $\mu$ g/ml in the mixture and incubated at 50°C for 4 h. Subsequently, isovolumetric tris-phenol extraction was performed twice, and washed with chloroform: isoamyl alcohol (24:1 v/v). 10% NaOAC was added to the recovered aqueous phase after centrifuging (10,000 r/min, 5 min), and then the DNA was precipitated by 2.5 times volume of cold ethanol at -20°C for 40 min before centrifuging (12,000 r/min, 10 min). In addition, the precipitated DNA was rinsed twice with 500  $\mu$ l cold 70% ethanol before drying at 37°C for 5–10 min. Then, the extracted DNA was stored at -20°C after dissolution in 50  $\mu$ l TE buffer at 65°C for 10 min.

**Citation:** ShiNan Dong Preparation of Genomic DNA of Microsporidia. **protocols.io**

[dx.doi.org/10.17504/protocols.io.igacbse](https://dx.doi.org/10.17504/protocols.io.igacbse)

**Published:** 14 Jun 2017

## Protocol