

Oct 04. 2019

a protein precipitation extraction method 👄

PLOS One

Idha Arfianti¹

¹Universiti Malaya



dx.doi.org/10.17504/protocols.io.3xrgpm6



EXTERNAL LINK

https://doi.org/10.1371/journal.pone.0221774

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Wiraagni IA, Mohd MA, Rashid RbA, Haron DEbM (2019) Validation of a simple extraction procedure for bisphenol A identification from human plasma. PLoS ONE 14(10): e0221774. doi: 10.1371/journal.pone.0221774

- 1 The frozen plasma was thawed at room temperature (25 ± 1 °C).
- 2 The thawed plasma was vortexed to ensure the sample was homogenous.
- To each 100 μL plasma sample, 50 μL of internal standard (IS) (containing 2 μg/mL of IS) was added, followed by the addition of 250 μL of acetonitrile (ACN).
- 4 The mixture was vortexed, shaken for 5 s, and centrifuged for 2 min at 14800 rpm.
- The supernatant was filtered with a 0.2 μm syringe filter, then transferred to new a vial.
- 6 Two microliters were then injected into the LC-MS/MS system.

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited