

# Analysis and characterization of the carbohydrate fraction of Aiptasia and coral tissue using targeted GC-MS

Hagen Gegner, Michael Ochsenkuehn, Christian Voolstra

## Abstract

The description and identification of metabolites using Gas Chromatography - Mass Spectrometry (GC-MS) is a powerful tool to study the cnidarian-dinoflagellate symbiosis. The applications range from first descriptions of metabolite profiles to the identification of specific metabolic responses (e.g. biomarkers) of the cnidarian metaorganism ('holobiont'), including its associated symbionts, under stress.

The following step-by-step protocol is optimized to characterize the carbohydrate fraction from Aiptasia anemones or coral fragments, as well as their algal symbionts, using a targeted GC-MS approach. The protocol was previously used in Ochsenkühn et al. (2017) where it identified the following carbohydrates (floridoside, inositol, mannitol, glucose, glycerol, galactose, ribose and fructose) as well as some amino acids (glycine, alanine, valine and proline).

In addition to the step-by-step protocol for sample preparation and derivatization, we provide detailed settings for the Agilent GC-MS system (GC (Agilent 7890A) and MS (Agilent 5975C)) in a separate document.

**Citation:** Hagen Gegner, Michael Ochsenkuehn, Christian Voolstra Analysis and characterization of the carbohydrate fraction of Aiptasia and coral tissue using targeted GC-MS. **protocols.io**

dx.doi.org/10.17504/protocols.io.s7fehjn

**Published:** 03 Sep 2018

## Guidelines

It is crucial that all steps, until the derivatization, are performed on ice to minimize degradation of metabolites. Normalization is essential for GC-MS data, as such, two options (dry-weight and total protein content) are included in this protocol.

## Materials

- 🧴 double distilled water (ddH<sub>2</sub>O) by Contributed by users
- 🧴 Ethanol, absolute 99.8% 10342652 by Contributed by users
- 🧴 Methoxamine (MOX) Reagent [TS-45950](#) by [Thermo Fisher Scientific](#)

 N-Methyl-N-(trimethylsilyl)trifluoroacetamide (MSTFA) [69479](#) by [Sigma Aldrich](#)

 4-Hydroxybenzoic acid [240141 Aldrich](#) by [Sigma Aldrich](#)

## Protocol

### Disclaimer

#### Step 1.

This protocol is a step-by-step guide to extract and analyze the carbohydrate fraction of Aiptasia and coral tissue. We included further notes to separate the algal symbiont fraction from the host fraction.

Normalization is essential for GC-MS data, as such two options (dry-weight and total protein content) are included in this protocol. Depending on the method chosen follow the respective 'STEP-CASE' normalization to dry-weight and/or to protein content.

To minimize degradation of metabolites, it is advised to work quickly and reduce waiting times. Further, it is crucial that all steps, until the derivatization, are performed on ice/in the cold.

If you are interested in other protocols related to the model organism Aiptasia:

 [PROTOCOL](#)

. [Getting started with the Aiptasia-Symbiodinium Model System](#)

CONTACT: [Aiptasia Model](#)

### Materials

#### Step 2.

### Chemicals

- MOX reagent 2% methoxamine HCL in pyridine (e.g. Thermo scientific)
- MSTFA, 1% TMCS (10x1 ml. (e.g. Sigma 69479-10X1ML))
- 4-Hydroxybenzoic acid (ISTD = Internal standard)
- Absolute ethanol  $\geq 99.8\%$  (GC)

- Standards for quantification  $\geq 99.8\%$  (GC) (This depends on the compounds targeted)
  - We use: Glucose, Mannose, Sucrose, Glycine

## Consumables

- Falcon tubes (50ml)
- GC sample vials with Teflon cap (2ml)
- Glass syringe e.g. Hamilton® syringe, 700 series
- Syringe filters - 0.2 micron (13 mm diameter)
- Beckman bottles for Beckmann Coulter centrifuge
- GC vial glass inlets (0.2ml)
- Plastic syringes (1ml)
- Parafilm
- Kim wipes (or any other tissue)

## Equipment

- GC/MS system
  - GC (Agilent 7890A)
  - MS (Agilent 5975C)
  - Autosampler (Agilent 7693)
- Concentrator system (Labconco Centrивap Complete)
- Ultracentrifuge (Beckmann Coulter Avanti J-26 XP)
- FreezeDryer (Ultradry)
- Ultrasonicator (Branson digital sonifier)
- Homogenizer (MicroDisTec MDT 125) (if needed)
- Centrifuge (Eppendorf 5415 R)
- Thermoblock fitting for GC vials
- Fume hood

## Sample collection

### Step 3.

At the end of your experiment:

- Rinse anemones or coral fragments with ddH<sub>2</sub>O to reduce the salt load
- Transfer anemones to cryotubes or in the case of coral fragments to a falcon tube or wrap in aluminum foil and snap freeze in liquid nitrogen

## Sample collection

### Step 4.

Store at -80 °C or continue with the protocol

Decide on the normalization method by following the step-case below:

- to protein content, directly follow the rest of the protocol
- to dry-weight, follow the steps and then go back to step-case 'normalization to protein content' to continue with the rest of the protocol

#### Sample preparation (dry-weight)

##### **Step 5 - Normalization to dry-weight.**

Thaw samples on ice and add ddH<sub>2</sub>O

- For Aiptasia add 1ml ddH<sub>2</sub>O
- For corals use an airbrush to remove the tissue or crush/pulverize fragments in liquid nitrogen and add up to 15ml ddH<sub>2</sub>O (the more you use the longer it takes in the next steps)

#### Sample preparation (dry-weight)

##### **Step 6 - Normalization to dry-weight.**

Transfer samples to a new, pre-weighed, falcon tube

#### Sample preparation (dry-weight)

##### **Step 7 - Normalization to dry-weight.**

Lyophilize samples over night using an type of lyophilizer (e.g. Ultradry)

#### Sample preparation (dry-weight)

##### **Step 8 - Normalization to dry-weight.**

After lyophilizing weigh the falcon tube + dried sample for later normalization to dry-weight

#### Sample preparation (dry-weight)

##### **Step 9 - Normalization to dry-weight.**

Continue with STEP-CASE: Normalization to total protein content and follow the rest of the protocol.

#### Sample preparation

##### **Step 5 - Normalization to protein content.**

Thaw samples on ice and add ddH<sub>2</sub>O

- For Aiptasia add 1ml ddH<sub>2</sub>O
- For corals use an airbrush to remove the tissue or crush/pulverize fragments in liquid nitrogen (if not done already)
  - add up to 15ml ddH<sub>2</sub>O (the more you use the longer it takes in the next steps)

#### Sample preparation

##### **Step 6 - Normalization to protein content.**

Sonicate using a tip ultrasonicator (e.g. Branson digital sonifier) on ice at 7 watts for 2min until

everything is dissolved.

#### Sample preparation

##### **Step 7 - Normalization to protein content.**

Remove cell debris by centrifuging at 4000g for 20min at 4 °C

Keep all samples on ice.

#### Sample preparation

##### **Step 8 - Normalization to protein content.**

Move supernatant to Beckman tube (50ml) and add 9 parts of absolute ethanol (cold) to 1 part of sample

#### Sample preparation

##### **Step 9 - Normalization to protein content.**

Leave them for 20min in a 4 °C fridge or on ice to precipitate DNA/RNA/Protein

#### SAFETY INFORMATION

**Balance the tubes for the ultracentrifuge, (+-0.05mg is ok) **

#### Sample preparation

##### **Step 10 - Normalization to protein content.**

Centrifuge precipitate at 20.000g for 20min at 4 °C in Beckman centrifuge

#### PROTOCOL

##### **. Protein extraction from Aiptasia**

CONTACT: [Jason Presnell](#)

##### **Step 10.1.**

Prepare anemones by placing in the incubator in artificial seawater (ASW) 3-4 days prior to extraction. Be sure to change out the (ASW) each day to remove any residual debris.

##### **Step 10.2.**

Transfer 4-5 large anemones to a small tissue grinder on ice with 0.5-1 mL of chilled extraction buffer. Upon removal from their original container, anemones should be blotted on Kimwipes/weigh paper to remove as much water as possible.

##### **Step 10.3.**

Transfer homogenate to a 1.5 mL tube and centrifuge at 14,000xg for 15 minutes at 4°C. This

step will pellet the dinoflagellates and cell debris, while the protein will be in the supernatant.

#### **Step 10.4.**

Remove the supernatant (be careful not to disturb the white lipid layer!) and place in a new tube.

#### **Step 10.5.**

Determine the concentration of your protein (we used the Bradford assay).

#### **Step 10.6.**

Aliquot your protein into 100 µl volumes and store in the -80°C freezer until further use. Freeze-thaw cycles should be avoided and protein can be kept for a short period in the fridge on ice.

### Sample preparation

#### **Step 11 - Normalization to protein content.**

Transfer the supernatant to a new falcon tube

### Sample preparation

#### **Step 12 - Normalization to protein content.**

Put a kim wipe on the tubes lid, make holes in the lid and cool them down to -80C

(optinal: use Kim wipes and parafilm)

After preparation start the lyophilization (e.g. Ultradry lyophilizer)

#### DURATION

24:00:00 : minimum lyophilization

#### SAFETY INFORMATION

**Cool samples down (-80 C) - otherwise the falcon tubes will overflow in the lyophilizer**

### Sample preparation

#### **Step 13 - Normalization to protein content.**

Possibility to pause the protocol and store dried samples

### Derivatization (preparation)

#### **Step 14 - Normalization to protein content.**

Preparation of:

1. GC vials with smaller labels
2. Heated Thermo Block for GC vials
3. Concentrator system (max. 40 °C) Coldtrap (on) time 99=continuous (Labconco Centrivap)

Complete)

4. Internal standard hydroxyl benzoic acid (ISTD) (1 ug/ul HBA in ddH2O)
5. Standards for quantification (that depends on the targeted compounds for quantification)

#### Derivatization (preparation)

##### **Step 15 - Normalization to protein content.**

Suspend dried samples in 500µl of ddH2O (its up to you, vials can hold 2ml)

- vortex vigorously

#### Derivatization (preparation)

##### **Step 16 - Normalization to protein content.**

Transfer suspended sample to a pre-labelled GC Vial

#### Derivatization (preparation)

##### **Step 17 - Normalization to protein content.**

Add 10µl of internal ISTD

- Add the internal Standard to all samples (including your calibration curves of standards for quantification)

#### Derivatization (preparation)

##### **Step 18 - Normalization to protein content.**

Dry samples using a Concentrator System (e.g. Labconco Centrivap Complete)

#### Derivatization (preparation)

##### **Step 19 - Normalization to protein content.**

Store samples

- Keep it dry, make sure lid is on and store in a vacuum

#### Derivatization (reaction)

##### **Step 20 - Normalization to protein content.**

Preparation of:

1. Thermoblock for GC vials
2. MOX (2% methoxamine HCL in pyridine)
3. MSTFA (MSTFA, 1% TMCS) (you can also use BSTFA))
4. Hamilton syringes (100µl)
5. Plastic syringes
6. Syringe Filters
7. Glass inserts for GC vials
8. Teflon lids for vials

## SAFETY INFORMATION

**MOX and MSTFA are volatile and toxic (corrosive) wear goggles, two layers of gloves and work under a fume hood. If a spill occurs immediately switch gloves.**

### Derivatization (reaction)

#### **Step 21 - Normalization to protein content.**

Add 50µl of MOX using a Hamilton syringe to the dried samples and close vials tightly

### Derivatization (reaction)

#### **Step 22 - Normalization to protein content.**

put samples in the thermo block for 1h at 75 °C

### Derivatization (reaction)

#### **Step 23 - Normalization to protein content.**

Remove samples from the thermo block and add 100µl of MSTFA using a Hamilton syringe

- Have a second syringe to break open clumps so everything is accessible for the MSTFA
- vortex

### Derivatization (reaction)

#### **Step 24 - Normalization to protein content.**

Close lids tightly and leave again for 1h in the 75 °C thermo block

### Derivatization (reaction)

#### **Step 25 - Normalization to protein content.**

After 1 h, particles can be removed by syringe filtration through 0.2-micron (13 mm diameter) filters and transfer 120µL to glass inserts placed inside GC vials

### Derivatization (reaction)

#### **Step 26 - Normalization to protein content.**

The vials are closed with PTFE septa screw caps or crimp caps and stored at -80 °C until measurement

### GC-MS measurement

#### **Step 27 - Normalization to protein content.**

- On the day of the measurement remove samples from the freezer and let them reach room temperature before measuring.
- Measure in batches of ca. 10 samples and start with quality control and calibration curves first.
- **Find the settings for the GC-MS attached in an extra document.**



Good luck.

## Sample preparation (Protein content normalization)

### **Step 5.**

## **Warnings**

The derivatization reaction steps have to be done in a fume hood as both reagents (MOX and MSTFA) are volatile and toxic (corrosive).

Wear goggles and a double layer of gloves and remove directly in case of any spill.