

Blasticidin S HCl in Mammalian Culture Determine Working Concentration

Sean Seaver

Abstract

The working concentration of [blasticidin S HCl](#) for use in mammalian cell gene selection can range anywhere from 2-100µg/mL depending on the cell lines used.

Optimal doses typically range from 2-20µg/mL

However, we recommend performing a kill curve to determine the optimal dose for your cell lines.

Citation: Sean Seaver Blasticidin S HCl in Mammalian Culture Determine Working Concentration. **protocols.io**
dx.doi.org/10.17504/protocols.io.cikucv

Published: 01 Nov 2014

Protocol

Step 1.

Plate 5×10^4 blasticidin sensitive host cells in 6-12 wells of a 24-well plate using appropriate media (without blasticidin) and incubate overnight.



REAGENTS



Blasticidin S HCl [TE-B001](#) by [P212121](#)



DURATION

24:00:00

Step 2.

Replace media with fresh media containing 6-12 different concentrations (one for each well) of blasticidin S HCl.



REAGENTS



Blasticidin S HCl [TE-B001](#) by [P212121](#)



DURATION

12:00:00

Step 3.

Continue replacing media with fresh blasticidin S HCl containing media every 3-4 days and monitor cells for normal activity.



REAGENTS



Blasticidin S HCl [TE-B001](#) by [P212121](#)



DURATION

12:00:00

Step 4.

After 5-7 days, determine the lowest concentration of blasticidin S HCl containing media that killed

100% of blasticidin sensitive host cells.

This is the optimal dose to use for your selection procedure.



REAGENTS



Blasticidin S HCl [TE-B001](#) by [P212121](#)



DURATION

24:00:00