

Ascorbic-EDTA Buffer

B.R. Cunningham, J.R. Brum, S.M. Schwenck, M.B. Sullivan, S.G. John

Abstract

For use in [Wet-mount Method for Enumeration of Aquatic Viruses](#)

Citation: B.R. Cunningham, J.R. Brum, S.M. Schwenck, M.B. Sullivan, S.G. John Ascorbic-EDTA Buffer. **protocols.io**
dx.doi.org/10.17504/protocols.io.c8xzxm

Published: 07 Jan 2016

Guidelines

Note: An alternative Ascorbate-EDTA Buffer can be made with MgCl₂ and Na₂EDTA if Mg₂EDTA is unavailable:

- Combine equal parts of 0.125 M Tris-base, 0.1 M Na₂-EDTA dihydrate, 0.2 M MgCl₂ hexahydrate, and 0.2 M ascorbic acid, adjust with 10 N NaOH to reach a pH of 6-7.

- 0.125 M Tris-base prepared by dissolving 0.151 g Tris-base in 10 mL ultrapure water

- 0.1 M Na₂-EDTA dehydrate prepared by dissolving 0.372 g Na₂-EDTA dehydrate in 10 mL ultrapure water

- 0.2 M MgCl₂ hexahydrate prepared by dissolving 0.407 g MgCl₂ hexahydrate in 10 mL ultrapure water

- 0.2 M ascorbic acid prepared by dissolving 0.352 g ascorbic acid in 10 mL ultrapure water

Protocol

Step 1.

Combine equal parts of 0.4 M Mg₂EDTA and 0.8 M ascorbic acid, adjust with 10 N NaOH to reach a pH of 6-7.

Step 2.

Prepare fresh within 48 hours of use