

Fixation of marine samples for flow cytometry analysis

Daniel Vaultot

Abstract

Protocol to fix marine samples for flow cytometry analysis of phytoplankton and bacteria with SYBR-Green.

Fix at least 2 samples per depth sampled and at least 6 to 10 depths per vertical profile

Reference

Marie, D., Rigaut-Jalabert, F. & Vaultot, D. (2014). An improved protocol for flow cytometry analysis of phytoplankton cultures and natural samples. *Cytometry*. 85. p.pp. 962-968.

Citation: Daniel Vaultot Fixation of marine samples for flow cytometry analysis. **protocols.io**

dx.doi.org/10.17504/protocols.io.dkk4uv

Published: 05 Aug 2015

Before start

Necessary equipment

- Vortex mixer
- Cryotubes 2 mL
- Liquid nitrogen tank

Materials

Glutaraldehyde EM Grade 25% G5882-50ML by [Sigma Aldrich](#)

Pluronic 10% [P5556-100mL](#) by [Sigma Aldrich](#)

Protocol

Step 1.

Prefilter seawater sample onto 200 µm mesh

 [AMOUNT](#)

2 ml Additional info:

Step 2.

Add 1.5 mL of prefiltered seawater sample to a 2 mL cryotube

Step 3.

Add Glutaraldehyde

 [AMOUNT](#)

15 µl Additional info:

 [REAGENTS](#)

Glutaraldehyde EM Grade 25% G5882-50ML by [Sigma Aldrich](#)

Step 4.

Add Pluronic (facultative)

 **AMOUNT**

2 µl Additional info:

 **REAGENTS**

Pluronic 10% [P5556-100mL](#) by [Sigma Aldrich](#)

Step 5.

Vortex

Step 6.

Wait 10 min

 **DURATION**

00:10:00

Step 7.

Flash freeze in liquid nitrogen tank

Warnings

Samples must be stored either in liquid nitrogen or at -80°C, not at -20°C because degradation will take place at the latter temperature