



Jun 05,  
2019

## Mammalian Cell Culture: Freezing

Kenneth Schackart<sup>1</sup>

<sup>1</sup>University of Arizona

Working

[dx.doi.org/10.17504/protocols.io.3szgnf6](https://doi.org/10.17504/protocols.io.3szgnf6)

 Kenneth Schackart 

### ABSTRACT

This protocol describes how to freeze cells that are being cultured in a tissue culture flask (T-75 or T-25).



### GUIDELINES

- Gloves must be worn at all times.
- Perform all tasks within biosafety cabinet.
- Anything entering biosafety cabinet must be generously sprayed with 70% ethanol (even you).
- When finished, wipe biosafety cabinet with 70% ethanol, and UV for at least 15 minutes.



### MATERIALS TEXT



- Cell culture media (at least 4 mL for T-75 or 1 mL for T-25)
- Fetal Bovine Serum (FBS)
- DPBS
- Trypsin-EDTA
- Dimethylsulfoxide (DMSO)
- 15 mL centrifuge tube
- Cyrovial(s)

#### Wash Cells



- 1 Remove media from flask.
- 2 Using serological pipet, add  4 ml DPBS to flask [ 1 ml for T-25].
- 3 Using serological pipet, remove DPBS and dispose into waste beaker.
- 4 Repeat the above 2 steps, for a total of 2 washes.

#### Trypsinize



- 5 Add  4 ml warmed Trypsin-EDTA to flask [ 1 ml for T-25].

- 6 Wait for cells to detach.
- 7 Add  **4 ml** warmed cell culture media to flask to neutralize Trypsin-EDTA [ **1 ml** for T-25]

#### Pellet Cells


- 8 Using a serological pipet, transfer cell suspension into 15 mL centrifuge tube.
- 9 Centrifuge cell suspension at  **1500 rpm** for  **00:03:00**.

#### Resuspend

- 10 Remove bulk of supernatant with serological pipet.
- 11 Remove remaining supernatant with 1000  $\mu$ L pipette. For small cell pellets, it is best to leave a little supernatant to avoid disturbing the pellet.
- 12 Add  **950  $\mu$ L** of either FBS or cell culture media and allow to sit for  **00:01:00** to make resuspension easier.




Some researchers prefer to use FBS as a freezing medium while others prefer whole media.

- 13 Gently pipette mix cell pellet until resuspended.
- 14 Add  **50  $\mu$ L** DMSO to cell suspension.



5% DMSO solution prevents ice crystals from forming in the liquid, minimizing cell rupture during freezing.

#### Freeze

- 15 Place cryovial in deep freezer cell vial container.
- 16 Allow at least  **04:00:00** in deep freezer before transferring cryovial to liquid nitrogen tank.



This prevents too rapidly freezing, which may cause cell rupture.

## 17 Update Lab Frozen Storage Inventory to reflect the new cell vial.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited