

Aug 17,  
2019

## Total RNA extraction from transgenic flies misexpressing foreign genes to perform Next generation RNA sequencing

PLOS One

Abijeet Mehta<sup>1</sup>, Agustin Luz-Madrigal<sup>2</sup>, Jian-Liang Li<sup>3</sup>, Panagiotis A. Tsonis<sup>4</sup>, Amit Singh<sup>5</sup><sup>1</sup>(Department of Biology, University of Dayton, Dayton, OH, USA ), <sup>2</sup>(Department of Biology, University of Dayton, Dayton, O, USA), <sup>3</sup>(Sanford Burnham Prebys Medical Discovery Institute at Lake Nona, Orlando, Florida, USA), <sup>4</sup>(Deceased), <sup>5</sup>(Department of Biology, University of Dayton, Dayton, OH, USA)

1 Works for me dx.doi.org/10.17504/protocols.io.5bng2me



Abijeet Mehta



### ABSTRACT

Due to absence of transgenic approaches in *Notophthalmus Viridescens* (*newt*), and conservation of genetic machinery across species, we generated transgenic *Drosophila melanogaster* to misexpress unique genes from newt. Novel newt genes cloned, and inserted at attP site in *Drosophila* were misexpressed ubiquitously using tubulin Gal-4. Sample (total RNA) for RNA sequencing was collected at 3<sup>rd</sup> instar larval stage during which major developmental events takes place in *Drosophila*. Total RNA was extracted, and purified using RNA clean and Concentrator<sup>TM</sup>. RNA quality was quantitated by calculating absorbance at 260 nm (A260) and 280 nm (A280) wavelengths using Nanodrop 2000 spectrophotometer. Good quality samples had A260/ A280 ratio greater than 2 and a peak at 260 nm. Our results show that following this protocol high quality of RNA was obtained. These high quality RNA samples were used for downstream processes e.g. Next generation RNA sequencing. Of the total 36,099 transcripts in *Drosophila*, 34,967 transcripts were detected, and 2775 transcripts were significantly regulated by misexpressing foreign gene (Unique gene from newt) in *Drosophila*. Genes involved in the developmental process, cell cycle, apoptosis, and immune response are among those that are highly enriched. Wingless/Wnt was one of the important evolutionarily conserved pathway that was differentially regulated.

### EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0220416>

### THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Mehta AS, Luz-Madrigal A, Li J, Tsonis PA, Singh A (2019) Comparative transcriptomic analysis and structure prediction of novel Newt proteins. PLoS ONE 14(8): e0220416. doi: [10.1371/journal.pone.0220416](https://doi.org/10.1371/journal.pone.0220416)

### GUIDELINES

Obtaining high quality RNA from the tissue samples is the first and most important step in performing molecular techniques like RNA sequencing. To achieve such quality it is important to include critical steps, and take precautions during sample collection, sample storage, RNA isolation, RNA purification, and RNA storage. Many of which are discussed in this protocol.

### MATERIALS

NAME	CATALOG #	VENDOR
<a href="#">RNA Clean &amp; Concentrator<sup>TM</sup>-5</a>	<a href="#">R1015</a>	<a href="#">Zymo Research</a>
<a href="#">RNAlater</a>	<a href="#">AM7020</a>	<a href="#">Thermo Fisher Scientific</a>
<a href="#">RNaseZap® RNase Decontamination Solution</a>	<a href="#">AM9780</a>	<a href="#">Life Technologies</a>
<a href="#">Chloroform</a>	<a href="#">319988</a>	<a href="#">Sigma</a>
<a href="#">TRIZOL reagent</a>	<a href="#">15596-026</a>	<a href="#">Invitrogen - Thermo Fisher</a>

### MATERIALS TEXT

Materials used were RNaseZap® (Thermofisher scientific Cat No. AM9780) as RNase decontamination solution. RNAlater (Thermo

Fisher, Cat. No. # AM7024) to store tissue sample before lysis, and RNA isolation. TRIzol Reagent (Thermo Fisher, Cat. No. # 15596926) was used to homogenate tissue sample and isolate good quality RNA. Chloroform (Sigma Catalog #: 319988/W205702) was used to separate homogenate into a clear upper aqueous layer containing RNA. RNA clean and Concentrator™ (Zymo research, Cat. No. R1080) kit was used to isolate high quality RNA recommended to perform RNA sequencing

#### SAFETY WARNINGS

When working with Trizol and/or similar chemical reagents while purifying nucleic acids, always work in a chemical fume hood and wear lab coat, disposable gloves and eye protection. Also strongly recommended to avoid mixing bleach with RNA extraction kit waste, and/or trizol. These kits, and trizol contain guanidine salts (e.g. guanidine thiocyanate and guanidine hydrochloride) that may produce hazardous gases, and/or strong acids when combined with bleach (sodium hypochlorite).

#### BEFORE STARTING

Before starting any step it is recommended to clean working space (RNA isolation bench), and all the equipments (e.g. scalpels, forceps, pipettes) with RNase decontamination solution e.g. use RNaseZap® (ThermoFisher scientific Cat No. AM9780). Simply spray RNaseZap® Solution onto the surface to be decontaminated and rinse it off with RNase-free water.

## 1 Sample collection and Storage

Transgenic *Drosophila melanogaster* were generated by genetivision (<https://www.genetivision.com/>). Briefly, transgenic flies were generated using microinjection-based  $\phi$ C31 integrase mRNA-mediated method. A cloned candidate gene using pUAST-attB plasmid containing both a transgene and donor sequence (attB) is coinjected along with  $\phi$ C31 integrase mRNA into attP-containing recipient embryos, resulting in the site-specific insertion of the transgene (Fish et al., 2007; Groth et al., 2004). Following this procedure 5 independent transgenic flies were generated with insertion of newt gene at Chromosome III. Foreign gene of interest was misexpressed ubiquitously using tubulin Gal-4. The third instar larvae were selected and were stored in RNeasy lysis buffer (Thermo Fisher, Cat. No. # AM7024) solution. The samples were then stored at 4°C. For long term storage we isolated larvae into empty RNeasy tubes and store them at -70°C.

## 2 RNA isolation and Purification

The tubes were briefly centrifuged and RNeasy lysis buffer was completely removed. We added 500  $\mu$ l of TRIzol Reagent (Thermo Fisher, Cat. No. # 15596926), which is used to isolate good quality RNA from tissue samples. TRIzol is a monophasic solution, primarily consisting of phenol and guanidine isothiocyanate along with other proprietary components, which disrupt samples, stabilize nucleic acids, and are compatible with downstream analysis. Tissue was homogenized in TRIzol using a 150 handheld homogenizer Motor (Fisher scientific, Cat. No. # 15-340-167). The solution was incubated for 30 minutes in ice. Centrifuged for 10 minutes at 4°C and supernatant was transferred to clean tube. Chloroform (Sigma Catalog #: 319988/W205702) was added and homogenate was allowed to separate into a clear upper aqueous layer. 200  $\mu$ l of the aqueous phase was transferred to RNeasy spin columns (Zymo research, Cat. No. R1080) columns, and the recommended protocol was followed, i.e., solution was passed through the RNA binding buffer, which binds RNA to the desired columns. Then RNA wash buffer was added to remove all the impurities from the column. Finally RNA was eluted in 20  $\mu$ l of molecular grade water (DNase/RNase free) and collected in a separate tube. The molecular grade water serves as the elution buffer. It releases the RNA from the column.

## 3 RNA quantitation, and Storage

RNA quantitation is an important and necessary step prior to most RNA analysis methods. Quality of RNA as well as concentration of RNA was determined by calculating absorbance at 260 nm (A260) and 280 nm (A280) wavelengths using Nanodrop 2000 spectrophotometer (Thermo Scientific). Good quality samples had A260/ A280 ratio greater than 2 and a peak at 260 nm. High quality total RNA was shipped in Dry ice to our collaborator for Next generation RNA sequencing.

#### 4 RNA Sequencing and Analysis

Illumina reads were mapped to the *Drosophila* genome dm6 using TopHat splice-aware aligner ([Kim et al., 2013](#)). Expectation-Maximization (EM) approach was used to estimate transcript abundance ([Jiang and Wong, 2009](#)). Reads per kilobase per million mapped reads (RPKM) approach was applied for within sample normalization ([Mortazavi et al., 2008](#)). Between sample normalization and differentially expressed test were performed by BioConductor DESeq package (v 1.20) which allows analysis of non-replicate experiments ([Anders and Huber, 2010](#)). The significant criteria were the detected transcript in at least one sample (RPKM>1), fold change over 2 and adjusted p-value less than 0.05.

Enriched gene ontology (GO) terms were identified using gene ontology enrichment analysis and visualization tool (GORilla) ([Eden et al., 2009](#)). Single ranked list of gene was chosen as a running mode criteria, and for enriched GO terms searched P-value threshold was equal to  $10^{-3}$ . Using reduce + visualize gene ontology (REVIGO) tool long lists of gene ontology terms obtained by running Gorilla were summarized by removing redundant GO terms ([Supek et al., 2011](#)). The remaining terms were visualized as tree map/ bar graph. We also used protein analysis through evolutionary relationships (PANTHER) 14.1 version as a tool to select set of enriched gene ontology terms for classifications by molecular function, cellular component, and Protein class ([Mi et al., 2019](#)).



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited