Human T Cell Activation with anti-CD3 Anithodies (clone UCHT1, OKT3 or HIT3a)

BioLegend, Inc.

Abstract

Citation: BioLegend, Inc. Human T Cell Activation with anti-CD3 Anitbodies (clone UCHT1, OKT3 or HIT3a). protocols.io

dx.doi.org/10.17504/protocols.io.extbfnn

Published: 26 May 2016

Guidelines

Materials:

- Sterile PBS
- Anti-human CD3 Antibody
- Clone UCHT1 (LEAF™ format, Cat. No. 300413/300414/300432; Ultra-LEAF™ format, Cat. No. 300437/300438)
- Clone OKT3 (LEAF™ format, Cat. No. 317303/317304/317315; Ultra-LEAF™ format, Cat. No. 317325/317326)
- Clone HIT3a (LEAF™ format, Cat. No. 300313/300314; Ultra-LEAF™ format, Cat. No. 300331/300332)
- Cell culture medium (e.g., RPMI-1640 or IMDM supplemented with 10% FBS and 2mM L-glutamine)
- Sterile single-cell suspension of Ficoll-Hypaque-purified peripheral blood mononuclear cells, isolated T cells, or T cell subsets
- 96-well flat-bottom tissue culture plates with lids (e.g., Costar® Cat. No. 3596)
- * Soluble forms of LEAFTM purified UCHT1 (1 μ g/ml) or LEAFTM purified HIT3a (0.01 0.1 μ g/ml) may be used to activate T cells from PBMC cell populations.

Protocol

Step 1.

Prepare a 10 μg/ml solution of anti-CD3 (clone UCHT1, OKT3, or HIT3a) in sterile PBS.

Step 2.

Dispense 50 μ l of the antibody solution to each microwell of the 96-well assay plate. For the unstimulated control wells, add 50 μ l of sterile PBS.



REAGENTS

96-well flat-bottom tissue culture plates with lids 3596 by Corning

Step 3.

Seal plate. Incubate at 37°C for 2 hours or 4°C overnight.

© DURATION

02:00:00

Step 4.

Aseptically decant antibody solution from the microwell plate.

Step 5.

Wash plate microwells 3 times with sterile PBS (was 1/3).

Step 6.

Wash plate microwells 3 times with sterile PBS (was 2/3).

Step 7.

Wash plate microwells 3 times with sterile PBS (was 3/3). Discard liquid.

Step 8.

Prepare single cell suspension of cells of interest in supplemented cell culture medium to $1-2 \times 10^{67}$ ml.

Step 9.

Aliquot 200 μ l cell suspension into plate microwells. Cover with lid. Incubate at 37°C in 5% CO2 and 100% humidity for 3 days.

© DURATION

12:00:00