

Qiagen - RNeasy mini kit for cells

Maysa Silva, Maryana Branquinho, Maria Cármén Sales

Abstract

Citation: Maysa Silva, Maryana Branquinho, Maria Cármén Sales Qiagen - RNeasy mini kit for cells. **protocols.io**
dx.doi.org/10.17504/protocols.io.ssbeean

Published: 21 Aug 2018

Before start

Clean the benches and all the material that will be used with alcohol 70.




Use tips with filter.

Add 4 volumes of 100% ethanol to the RPE buffer.

To lyse the cells, you can use β -mercaptoethanol or 2 M dithiothreitol. For each 1ml of the RLT buffer, add 10 μ l of β -mercaptoethanol. For each 1ml of the RLT buffer, add 20 μ l of 2 M dithiothreitol.

You must use a maximum of 1×10^7 cells.

Materials

- ✓ Buffer RPE by Contributed by users
-  RLT Buffer by [Qiagen](#)
-  Ethyl Alcohol [E7023](#) by [Sigma](#)
- ✓ RNase-free water by Contributed by users
-  RW1 buffer 74106 by [Qiagen](#)
- ✓ Ethyl alcohol 70% by Contributed by users

Protocol

RNA extraction

Step 1.

If you are using less than 5×10^6 add 350 μ l of the RLT buffer prepared above. If the number of cells is greater than this, use 700 μ l.

RNA extraction

Step 2.

Homogenize the samples using vortex or use QIAshredder. If necessary use an insulin syringe to break the cells.

RNA extraction

Step 3.

Add 1 volume of 70% ethanol to the lysate and homogenize with the pipette.

RNA extraction

Step 4.

Transfer up to 700 µl of the sample, including any precipitate, to an RNeasy Mini spin column placed in a 2 ml collection tube.

RNA extraction

Step 5.

Centrifuge for 15 seconds at 8000 g. Discard the flow-through.

RNA extraction

Step 6.

Add 700 µl of the RW1 buffer to the column and centrifuge for 15 seconds at 8000 g. Then discard the flow-through.

RNA extraction

Step 7.

Add 500 µl of the RPE buffer to the column and centrifuge for 15 seconds at 8000 g. Then discard the flow-through.

RNA extraction

Step 8.

Add 500 µl of the RPE buffer to the column and centrifuge for 2 minutes at 8000 g.

RNA extraction

Step 9.

This step is optional. Place the column in a new 2 mL collection tube and centrifuge at full speed for one minute to dry the membrane.

RNA extraction

Step 10.

Place the column in a new 1.5 mL collector tube and add 30 to 50 µl RNase Free water and centrifuge for 1 minute to 8000 g to elute the RNA.

RNA extraction

Step 11.

If you expect to have more than 30 µg of RNA, repeat step10 again using 30 to 50 µl of RNase-free water. Or, use the elution you acquired in step 10. Reuse the pickup tube from step 10.

RNA extraction

Step 12.

Stock the sample at -80 ° C.