

Diplonema Genomic DNA isolation

Binnypreet Kaur^{1,2}, Drahomíra Faktorová^{1,2}, , Priscila Peña-Díaz¹ and Julius Lukeš^{1,2}

Abstract

E-mail: binny@paru.cas.cz

Citation: Binnypreet Kaur^{1,2}, Drahomíra Faktorová^{1,2}, , Priscila Peña-Díaz¹ and Julius Lukeš^{1,2} Diplonema Genomic DNA isolation. **protocols.io**

dx.doi.org/10.17504/protocols.io.hfyb3pw

Published: 11 Jul 2018

Protocol

Step 1.

Spin down cells (10^6) for 5 mins at 1800g to get clear pellet and try to remove all the traces of sea salt

Step 2.

Resuspend cells in 0.5 ml (per 10^8 cells) lysis buffer (10mM Tris-HCl pH 8.0, 5mM EDTA, 200mM NaCl, 0.2% w/v SDS), add 5µl RNase A (20 mg/ml) (Fermentas) and mix by inverting the tube several times

Step 3.

Incubate for 10 min at 60°C and add 2.5µl Proteinase K (10 mg/ml) (Fermentas). Incubate 1 hour at 60°C.

Step 4.

Add 250 µl ice cold 5 M NaCl, and incubate on ice for 10 mins (protein precipitation step—salting out of protein).

Step 5.

Spin down at 16000-20000 RPM for 15 min and transfer supernatant to new tube

Step 6.

Add an equal volume (about 700-750 µl) of isopropanol and mix by inversion (genomic DNA precipitation). Pending on amount of gDNA, incubate 15 mins at room temperature (RT) or 1 hour at -20°C to improve the recovery

Step 7.

Spin down 10 min at RT to pellet gDNA precipitate

Step 8.

Wash pellet 2X times with 500 µl 70% (v/v) ethanol by centrifuge at 16000-20000 RPM for 15 min

Step 9.

Remove the supernatant completely by aspiration and air dry gDNA pellet for 10-25 min at room temperature

Step 10.

Resuspend gDNA in 10mM Tris-HCl (pH 8.0) leave overnight at RT for gDNA to resolve completely
