

Lysis Buffer (10mM Tris-HCl, 2mM EDTA, 1% SDS)

Sarah Hessen-Schmidt

Abstract

Citation: Sarah Hessen-Schmidt Lysis Buffer (10mM Tris-HCl, 2mM EDTA, 1% SDS). **protocols.io**

dx.doi.org/10.17504/protocols.io.fjqbkwmw

Published: 10 Aug 2016

Protocol

Step 1.

Add 2.5ml 1M Tris-HCl pH 7.5 to an autoclaved bottle

PROTOCOL


. 1M Tris-HCl

CONTACT: [Sarah Hessen-Schmidt](#)

Step 1.1.

Dissolve 30.275g Trizma base in 200ml MilliQ

REAGENTS

 Trizma Base 93362 by Contributed by users

Step 1.2.

Adjust to pH 7.5 with HCl (35ml)

Step 1.3.

Adjust total volume to 250ml with MilliQ

Step 1.4.

Autoclave

Step 2.

Add 1ml 0.5M EDTA

PROTOCOL

. **0.5M EDTA**

CONTACT: [Sarah Hessen-Schmidt](#)

Step 2.1.

46.525g EDTA disodium 2H₂O in 200ml MilliQ



REAGENTS

Ethylenediaminetetraacetic acid disodium salt dihydrate E5134-1KG by [Sigma Aldrich](#)

Step 2.2.

Adjust to pH 8.0 with NaOH pellets (5g)

Step 2.3.

Autoclave

Step 3.

Add 2.5g SDS



REAGENTS

Sodium dodecyl sulfate 436143-25G by [Sigma Aldrich](#)

Step 4.

Adjust total volume to 250ml with MilliQ