

Checking DNA Concentration with Agarose Gel

Matthew Sullivan

Abstract

Citation: Matthew Sullivan Checking DNA Concentration with Agarose Gel. **protocols.io**

dx.doi.org/10.17504/protocols.io.c35yq5

Published: 25 Jan 2016

Guidelines

This protocol comes from a group of other protocols. This protocol is (4) of (4):

1. ['Large Volume Marine Cyanophage Phage Protocols'](#)
2. ['DNA Extraction Protocol'](#)
3. ['DNA Precipitation Protocol'](#)
4. ['Checking DNA Concentration with Agarose Gel'](#)

Needed:

- 0.7% Agarose
- (Agarose)
- (TAE Buffer)
- (HINDIII (lambda) DNA)
- Nanodrop

Materials

- ✓ SYBR Gold Nucleic Acid Gel Stain [S-11494](#) by Contributed by users
- ✓ Lambda DNA/HindIII Marker [SM0101](#) by Contributed by users

Protocol

Step 1.

Prepare 0.7% agarose

✓ [PROTOCOL](#)

• [0.7% Agarose](#)

CONTACT: [VERVE Team](#)

Preparation of 0.7% Agarose

Step 1.1.

Mix 0.35 g agarose in 50 ml TAE running buffer

 [AMOUNT](#)

0 g Additional info:

Step 2.

Post-stain with SYBR Gold

Step 3.

While the gel is cooling, prepare samples:

1 µl of 10x loading buffer per 10 µl of sample

Step 4.

Prepare the appropriate dilutions of the appropriate marker (commonly 10ng of HINDIII [lambda] DNA)

 [AMOUNT](#)

10 ng Additional info:

 [REAGENTS](#)

✓ Lambda DNA/HindIII Marker [SM0101](#) by Contributed by users

Step 5.

Add running buffer to cover gel

Step 6.

Gently rock the comb before pulling it out

Step 7.

Run at 50V about 2-4 hours to get clean band separations

 [DURATION](#)

03:00:00

 [NOTES](#)

VERVE Team 17 Jun 2015

Note: (if in hurry, watch for bromophenol blue marker to run about 3-4 cm)

Step 8.

Post-stain the gel in SYBR Gold 'bath' (20ul SYBR: 200ml TAE for 30 minutes)

 [AMOUNT](#)

20 µl Additional info:

 [DURATION](#)

00:30:00

Step 9.

Compare gel quantification to that obtained using Nanodrop to get a feeling for how 'clean' the DNA prep is

 [NOTES](#)

VERVE Team 21 Jan 2016

2 µl volume needed; only good down to the ~4ng/µl concentrations; A260/A280 estimates are very approximate in that range.