

In vitro Transcription assay for cyanobacterial DNA-dependent RNA Polymerase Version 2

Lutz Berwanger

Abstract

In vitro transcriptions assay for cyanobacterial RNA polymerase. This protocol can be used either to do an in vitro transcription with small DNA fragments with desired promoter region or to use it with a promoter region integrated in a plasmid.

Citation: Lutz Berwanger In vitro Transcription assay for cyanobacterial DNA-dependent RNA Polymerase. **protocols.io** dx.doi.org/10.17504/protocols.io.mpvc5n6

Published: 17 Jan 2018

Guidelines

RNA experiments! Always wear gloves and prepare probes under RNase free and cooled conditions!

Before start

Prepare all Buffer, measure RNAP and DNA concentrations

Materials

✓ in vitro transcription buffer by Contributed by users

Protocol

Preparations

Step 1.

In vitro transcription Buffer:

20 mM Tris-HCl pH 7,9

40 mM KCl,

15 mM MgCl₂

Prepartions

Step 2.

Polyacrylamid Gel:

10-20% 19:1 acrylamide:bisacrylamide

8M urea

1x TBE

📌 NOTES

Lutz Berwanger 17 Jan 2018

For smaller fragments use higher concentrations of 19:1 acrylamide:bisacrylamide.

Prepartions

Step 3.

Loading dye should contain at least 50% formamide or 8M urea. You can use commercial RNA loading dye, also.

Example recipie:

8M urea saturated with formamide

10mM EDTA

1xTBE,

bromophenol blue and xylene cyanol (amounts are empirical, just to make it visible on the gel)

In vitro Transcription

Step 4.

Mix:

100 nM RNAP

10 nM promoter-containing DNA

100µM NTPs

NOTES

Lutz Berwanger 17 Jan 2018

If desired one could use P32 labelled NTPs. Then it is recommendable to decrease the concentration of the same cold NTP to 10uM, to increase labelling efficiency.

In vitro Transcription

Step 5.

Incubate the reaction at 34°C for at least 3 hours.

Stop with equal volume of the loading dye. Before loading onto denaturing PAA gel boil the Samples for 5 min at 95°C and let them cool on ice.

NOTES

Lutz Berwanger 17 Jan 2018

If desired, one could expand reaction duration.

Visualization

Step 6.

Visualization of the products on the gel, can be done with nucleic acid gel stain of your choice. For small products it is recommended to stain the gel over night.

REAGENTS

GelRed™ Nucleic Acid Gel Stain, 10,000X in Water [G-725](#) by [Gold Biotechnology](#)

✓ SYBR Gold Nucleic Acid Gel Stain [S-11494](#) by Contributed by users

NOTES

Lutz Berwanger 17 Jan 2018

Two common used stains are indicated.

Warnings

Harmless