



Purification of total RNA (microRNA and mRNA) from liquid samples [↗](#)

PLOS One

[atsuko mizoguchi](#)¹

¹toray industries Inc.,

[dx.doi.org/10.17504/protocols.io.vu3e6yn](https://doi.org/10.17504/protocols.io.vu3e6yn)

[atsuko mizoguchi](#)
Toray Industries, Inc

ABSTRACT

Total RNA was extracted from cells and the cultured supernatant and exosome which was extracted from the same amount of supernatant using 3D-Gene® RNA extraction reagent. Comprehensive miRNA and mRNA expression analysis were performed using 3D-Gene® Human miRNA Oligo Chip (miRBase ver.21) and mRNA Oligo Chip (Toray Industries, Inc.) which are featured with the columnar structure and bead-mixing for high sensitivity.

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0209750>

PROTOCOL STATUS

Working

1

For 300ul of cultured supernatant and exosome in 1.7ml centrifuge tube,
add 900 µl 3D-Gene RNA extraction (Toray Industries, Inc., Kanagawa, Japan) reagent from a liquid sample kit.

2

Homogenize by vortexing for 1min to mix.
Place the tube containing the homogenate on the benchtop

00:03:00 incubation 25 °C

3

00:03:00 incubation 0 °C

4

Centrifuge for 10 min at 12,000 x g at 4°C without brakes.
After centrifugation, heat the centrifuge up to room temperature.

00:03:00 incubation 25 °C

5

Transfer the 300ul of upper aqueous phase to a new collection tube.

6

Add 450ul of 100% ethanol and mix thoroughly by pipetting up and down several times. Do not centrifuge.

- 7 Transfer the 750ul of sample into the silica-based membrane column and purification according to general nucleic acid purification protocol.
- 8 Elute the RNA with RNase free water, and concentrate to 4ul.
- 9 Comprehensive miRNA and mRNA expression analysis were performed using 3D-Gene® Human miRNA Oligo Chip (miRBase ver.21) and mRNA Oligo Chip(Toray Industries, Inc.) which are featured with the columnar structure and bead-mixing for high sensitivity.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited